FAQ

MGIEasy FS DNA Library Prep Set V2.0

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Version: A0

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1.1 Which sample can be used in MGIEasy FS DNA Prep Set?

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A: MGIEasy FS DNA Prep Set is suitable for human (saliva, blood, tissue), plant, animal, fungus, Meta gDNA.

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1.2 Do I have to prepare any special equipment for library construction?

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A: Thermocycler, Pipets, Magnetic rack DynaMag is needed. Rack with high magnetism is recommended like ALPAQUA,

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1.3 What's the requirement of the sample quality?

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A: Integrity and Purity are needed for MGIEasy FS DNA Library Prep Set. It is strongly recommended to use high quality genomic DNA (gDNA) samples ($A_{260}/A_{280}=1.8-2.0$, $A_{260}/A_{230}>2.0$, no degradation or degraded slightly) for fragmentation. That fragmentation test using part of gDNA is recommended for gDNA with low quality to adjust the fragment time, but risk of failure remains.

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1.4 Is there any recommendation for different input in MGIEasy FS DNA Library Prep Set?

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A: The recommendation of different input shows as below. The input could be increased if the quality is pool.

Sample type	Input Range	Recommended	Recommended
		Input	Concentration
Complex genome	50-400 ng	200 ng	≥15 ng/μL
Simple genome	5-400 ng	100 ng	≥7.2 ng/μL
Microbiome	5-400 ng	100 ng	≥7.2 ng/μL
Meta	5-400 ng	100 ng	≥7.2 ng/μL
Pathogenic samples	5-100 ng	100 ng	≥7.2 ng/μL

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1.5 How to identify Complex/Simple genome?

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A: Simple genome means genome size less than 20 Mb. For example: Microbiome/Meta genome/Virus. Complex genome means genome size more than 20 Mb. For example: human, polyploidy plant.

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1.6 Is there any stopping point during library Preparation? How long can the product be saved?

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A: The stopping points are after each cleanup step. The purified DNA product can be stored at -20°C for up to 6 months, the circularized ssDNA libraries can be stored for up to 3 months.

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1.7 What's the exact temperature when it refers to "heated lid: On"?

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A: The heated lid should be set 5-10°C above the reaction when using the Thermocycler with adjustable lid. For example, if the reaction condition is 37°C, the heated lid should be set 42°C-47°C.

We recommend the heated lid remain unclosed when the condition below 25°C when using the Thermocycler without adjustable lid. 105°C could be set under the other conditions.

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1.8 Can I use fragment with insert size exceed 250-350bp range?

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A: Fragment with insert size ± 50 bp can be used for normal preparation. 400bp can be accepted for PE100 sequencing, but size less than 250bp is not recommended for PE150 because of the overlap.

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A: there is no PE200 sequencing data yet. Library for PE200 can be done using less fragmentation time and less size-selection beads. For example, the fragmentation condition could be set to "32°C for 3minutes", and the size selection condition could be set to "0.55+0.15"

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1.10 Can I use the WGS library to do the SE sequencing?

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A: PE sequencing is recommended for complex genome species sequencing. SE sequencing can also be done with good sequencing quality but MGI doesn't commit the alignment performance.

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A: The size selection step can be changed to one-turn purification in the FS library prep set. Ranged size may influence the sequencing quality.

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1.12 Can I use the amplicon product to make libraries?

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A: Make the library by fragmenting the size as we recommended when dealing with the amplicon product.

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1.13 Can I use the adapter from the other vendor?

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A: Adapter with a "T" base extrude is all available. MGI doesn't commit any performance of the adaptor from the other vendor.

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1.14 What is the difference between MGIEasy "FS" and MGIEasy "Universal"?

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A: the "FS" one includes Fra Enzyme to do the fragmentation for libraries without fragmentation machine and automatic platform.

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1.15 Can reagents that shared the same name but from different set (such as Universal/FS) be used as substitutions?

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A: Reagents with the same item number could be used freely.

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2.1 How many frag-reaction can MGIEasy FS DNA Library Prep Set supply?

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A: we recommend costumers to do several test before formal experiment. MGIEasy FS DNA Prep Set(Lot: 1000005254) have 18 reaction for real.(2 for extra experiment), MGIEasy FS DNA Prep Set (Lot: 1000005256) have 120 reaction for real.(24 for extra experiment), If the set is used in MGI SP-960, the test should be no more than 2 samples.

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2.2 What will influence the fragmentation result?

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A: ① Under 32°C, the more fragmentation time being cost, the smaller the size is. ② more EDTA decreases the fragmentase activity. ③ High concentration of divalent cation improves the fragmentase activity. ④ the size increased slightly when more gDNA was input into fragmentation. ⑤ Wrong storage condition will significantly decrease the activity.⑥ Vortex is necessary for fragmentation Mix preparing, or the fragmentation will be effected seriously.

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A: If the DNA yield is considered to be more than 100ng, the quantification step should be done and choose no more than 100ng to do next ER step. Otherwise the quantification can be skipped.

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2.4 What is the concentration of the MGIEasy DNA Adapters?

MGIEasy DNA Adapters 的浓度是多少?

A: 10 μ M

答: 10 μ M

2.5 How to make sure the buffer mix is homogeneous as the ligation buffer is sticky?

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A: Thaw the Ligation Buffer at room temperature and vortex 6 times (6 s each), which ensuring homogeneous. If the mix remain sticky after sample added, 6 times (6 s each) vortexing should be done again, or the yield will be influenced.

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2.6 Why is 20 ul TE Buffer added after ligation step? What will happen when forgot?

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A: 20 ul TE Buffer can mitigate the adaptor residue while the efficiency remains good. If this step is missed: ① during step 3.5.2, 20 ul TE could be added into the beads. ② If the purification step is finished, 1X beads purification once more is recommended.

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2.7 What is the concentration of the PCR Primer Mix

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A: 20 uM

答: 20 μ M

2.8 What happened if I got no yield of single circle DNA when the PCR Product is normal?

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A: Most likely there is something wrong in the circularization step. It is recommended to add a positive control combine

with the product to do the circularization again.

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2.9 Can I use the DNB to do sequencing if the circularization efficiency is very low?

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A: We don't recommend doing that.

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2.10 Can I pool single circle with various insert size?

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A: if the size range >100nt, the sequencing quality may be influenced. We don't recommend to pool libraries like that.

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2.11 Is MGIEasy FS DNA library Prep set compatible for automatic machine?

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A: MGIEasy FS DNA Library Prep Set is compatible for MGISP-100, MSP-960 and most kind of automatic machine (need additional test to confirm). Set 1000017572 is compatible for preparing 96 single circle libraries each time on MGISP-960. Set 100006988 is compatible for preparing 96 PCR product using MGISO-960 and make 16 single circle using MGISP-100.

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