



Part No.: SOP-013-B01-080

DNBSEQ-T7RS

High-throughput Sequencing Set

User Manual

Version: 6.0

Address 1: Building B13, No.818, Gaoxin Avenue, East Lake High-Tech Development Zone,
430075, Wuhan, P.R.China

Address 2: Building 24, Stage 3.1, BioLake Accelerator, No.388, 2nd Gaoxin Road, East Lake
High-Tech Development Zone, 430075, Wuhan, P.R.China

Address 3: Main Building and Second floor of No. 11 Building, Beishan Industrial Zone, Yantian
District, 518083 Shenzhen, Guangdong, P.R.China

E-mail: MGI-service@genomics.cn

Website: en.mgi-tech.com

Research Use
Only

MGI Tech Co., Ltd.
Wuhan MGI Tech Co., Ltd.

About the user manual

This user manual is applicable to DNBSEQ-T7RS High-throughput Sequencing Set. The manual version is 6.0.

This user manual and the information contained within are proprietary to MGI Tech Co., Ltd./Wuhan MGI Tech Co., Ltd. (hereinafter called MGI), and are intended solely for the contractual use of its customer in connection with the use of the product described herein and for no other purpose. Any person or organization cannot entirely or partially reprint, copy, revise, distribute or disclose to others the user manual without the prior written consent of MGI. Any unauthorized person should not use this user manual.

MGI does not make any promise of this user manual, including (but not limited to) any commercial of special purpose and any reasonable implied guarantee. MGI has taken measures to guarantee the correctness of this user manual. However, MGI is not responsible for any missing parts in the manual, and reserves the right to revise the manual and the device, so as to improve the reliability, performance or design.

Figures in this user manual are all illustrations. The contents might be slightly different from the kit, please refer to the kit purchased.

Unless otherwise informed, all sequencers and sequencing reagents are not available in Germany, USA, Spain, UK, Hong Kong, Sweden, Belgium and Italy.

NextSeq™ and TruSeq™ are trademarks of Illumina, Inc., or its subsidiaries. Qubit™ is the trademark of Thermo Fisher Scientific, or its subsidiaries. Other company, product names, and other trademarks are the property of their respective owners.

©2019-2022 MGI Tech Co., Ltd./Wuhan MGI Tech Co., Ltd. All rights reserved.

Revision history

Version	Date	Description
6.0	February 2022	<ul style="list-style-type: none">● Updated logo of MGI.● Updated transport temperature.● Deleted the description of the small RNA adapter in App-D.
5.0	July 2021	<ul style="list-style-type: none">● Updated disclaimer.● Added the validity of reagents.
4.0	July 2021	<ul style="list-style-type: none">● Added App-D PE150.● Updated the DNB preparation and load DNB for PE150.● Updated the sequencing cartridge well No.9 and No.10 reagent adding.
A2	December 2020	Updated the logo, website address and mailbox, and template of the manual .
A1	November 2020	<ul style="list-style-type: none">● Added the SE35, SE50, SE100 and PE150 read length.● Added Dual barcode PE sequencing.● Added the stLFR PE100.● Added the App-A PE100 and App-A PE150.● Updated part of PUI figures.● Revised DNB pooling.● Added an attachment for quantify DNB.
A0	December 2019	Initial release

Sequencing set

Catalog number	Name	Version
1000020243	DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE35)	1.0
1000016102	DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE50)	1.0
1000016103	DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE100)	1.0
1000028455	DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE100)	2.0
1000028454	DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE150)	2.0
1000019251	DNBSEQ-T7RS High-throughput Sequencing Set (stLFR FCL PE100)	1.0
940-000005-00	DNBSEQ-T7RS High-throughput Sequencing Set (App-A FCL PE100)	2.0
940-000003-00	DNBSEQ-T7RS High-throughput Sequencing Set (App-A FCL PE150)	2.0
1000020834	CPAS Barcode Primer 3 Reagent Kit	2.0
1000014047	High-Throughput Barcode Primer 3 Reagent Kit (App-A)	1.0
1000028550	High-Throughput Sequencing Primer Kit (App-D)	1.0

---This page is intentionally left blank.---

Contents

Chapter 1 Introduction	1
1.1 Applications	1
1.2 Sequencing principle	1
1.3 Data analysis	1
1.4 Sequencing read length	1
1.5 Sequencing time	2
1.6 Precautions and warnings	3
 Chapter 2 Sequencing sets and self-prepared consumables	 4
2.1 List of sequencing set components	4
2.2 Self-prepared equipment and consumables	14
 Chapter 3 Sequencing workflow	 16
 Chapter 4 Making DNB	 16
4.1 Insert size recommendation	16
4.2 Library concentration and amount requirement	18
4.3 Library pooling	18
4.3.1 Number of samples that can be pooled together	18
4.3.2 Verifying the base balance for barcode	20
4.4 Making DNB	20
4.4.1 Preparing DNB for the FCL SE35, FCL SE50, FCL SE100 and FCL PE100	20
4.4.2 Preparing DNB for the FCL PE150 kit	24
4.4.3 Preparing DNB for the App-A FCL PE100 kit	27
4.4.4 Preparing DNB for the App-A FCL PE150 kit	31
4.4.5 Preparing DNB for the App-D FCL PE150 kit	34
4.4.6 Preparing DNB for the stLFR FCL PE100 kit	37

4.5 Quantifying DNB and pooling	40
4.5.1 Quantifying DNB	40
4.5.2 DNB pooling	41

Chapter 5 Loading DNB 42

5.1 Preparing the post load plate and buffers	42
5.1.1 Preparing the post load plate and buffer for FCL SE35/ SE50/SE100/PE100, App-A FCL PE100, stLFR FCL PE100 sequencing	43
5.1.2 Post load plate and buffers preparation for FCL PE150, App-A FCL PE150 and App-D FCL PE150 sequencing	44
5.2 Preparing the sequencing flow cell	46
5.3 DNB loading	46

Chapter 6 Preparation before sequencing 54

6.1 Preparing the sequencing reagent cartridge	54
6.2 Preparing the washing cartridge	58
6.3 Filling the pure water container	60

Chapter 7 Sequencing 61

7.1 Loading the reagent cartridge	61
7.2 Entering sequencing interface	62
7.3 Loading the flow cell	62
7.4 Sequencing parameters	63
7.5 Reviewing parameters	66
7.6 Starting sequencing	67

Chapter 8 Device maintenance 69

8.1 Terminology and definition	69
8.2 Preparing wash reagents	70
8.3 Washing cartridge	70

8.4 Wash procedures	71
8.4.1 MGIDL-T7RS manual wash	71
8.4.2 DNBSEQ-T7RS manual wash	71
Chapter 9 Troubleshooting	72
9.1 Low DNB concentration	72
9.2 Abnormal negative pressure	72
9.3 Bubbles	72
9.3.1 Bubbles in MGIDL-T7RS	72
9.3.2 Bubbles in DNBSEQ-T7RS	73
9.4 Impurities	73
9.5 Pump fails	73
9.6 Reagent kit storage	74
Chapter 10 Warnings and precautions	74
Appendix 1 Qubit ssDNA assay kit	75
Appendix 2 Manufacturer	77

---This page is intentionally left blank.---

Chapter 1 Introduction

This manual explains how to perform sequencing by using the DNBSEQ-T7RS High-throughput Sequencing Set and includes instructions on sample preparation, Flow Cell preparation, sequencing kit storage, the sequencing protocol and device maintenance.

1.1 Applications

DNBSEQ-T7RS High-throughput Sequencing Set is specifically designed for DNA or RNA sequencing on DNBSEQ-T7RS. This sequencing set is intended to be used for scientific research only and cannot be used for clinical diagnosis.

1.2 Sequencing principle

This sequencing set utilizes DNBSEQ technology. A sequencing run starts with the hybridization of a DNA anchor, then a fluorescent probe is attached to the DNA Nanoball (DNB) using combinatorial probe anchor sequencing (cPAS) chemistry. Finally, the high-resolution imaging system captures the fluorescent signal. After digital processing of the optical signal, the sequencer generates high quality and highly accurate sequencing information.

1.3 Data analysis

During the sequencing run, the control software automatically operates base calling analysis software and delivers raw sequencing data outputs for secondary analysis.

1.4 Sequencing read length

In the sequencing run, the number of sequencing cycles depends on the sequencing read length. For example, a PE100 cycle run performs reads of 100 cycles from each end, for a total of 200 (2×100) cycles. At the end of the insert sequencing run, an extra 10 cycles of index read can be carried out, if required.


Table 1 Sequencing cycle

Sequencing read length	Read 1 read length	Read 2 read length	Barcode read length	Total read length	Maximum cycles
SE35	35	/	10	35+10	55
SE50	50	/	10	50+10	70
SE100	100	/	10	100+10	120
PE100	100	100	10	200+10	220
PE150	150	150	10	300+10	320
App-A PE100	100	100	10	200+10	220
App-A PE150	150	150	10	300+10	320
App-D PE150	150	150	10	300+10	320
stLFR PE100	100	100	42+10	200+52	252

1.5 Sequencing time

Table 2 Theoretical sequencing time (hours) for various sequencing kits on DNBSEQ-T7RS

Read length	Single flow cell (hours)	Four flow cells (hours)	DNB preparation (hours)	DNB loading (hours)
SE35	4.5	5.0	1	2
SE50	5.5	6.0	1	2
SE100	9.0	10.5	1	2
PE100	17.5	20.5	1	2
PE150	24.5	29.0	1	2
App-A PE100	17.5	20.5	1	2
App-A PE150	24.5	29.0	1	2
App-D PE150	24.5	29.0	1	2
stLFR PE100	21.0	24.5	1	2

-  **NOTE**
- Sequencing run time for both single flow cell and four flow cells only refer to the time elapsing from the “start” to the “finish” of the sequencing run. The time used for DNB preparation, DNB loading and Write FQ is not included. Write FQ for a single flow cell will take about 1.5 hours.
 - Two flow cells can be loaded with DNB concurrently using one MGIDL-T7RS instrument. Total time is about 2 hours.
 - Sequencing run time is based on the DNBSEQ-T7RS instrument with standard model, actual sequencing run time could vary among different instruments.
 - Sequencing run time includes the time for the single barcode (10 cycles) sequencing, except for the stLFR PE100 where the time for 42+10 barcode cycles run is included.

1.6 Precautions and warnings

- These products are used for scientific research only. Please read this manual carefully before use.
- Please also read the user manual of the DNBSEQ-T7RS and MGIDL-T7RS instrument before the experiment.
- Personal injury can occur through ingestion, inhalation, skin contact and eye contact with these reagents. If it happens, please rinse with large amount of water, then see doctor immediately.
- Please follow local governmental safety guidance to dispose of the sequencing run waste and unused contents.
- All components in the sequencing kits are designed for single use. Do not use them repeatedly.
- Do not use expired sequencing reagents.

Chapter 2 Sequencing sets and self-prepared consumables

2.1 List of sequencing set components

It is worth reminding that:

- For stLFR FCL PE100 and single barcode FCL SE35 / SE50 / SE100 / PE100 / PE150 sequencing run, only its corresponding sequencing set is required.
- A dual barcode FCL PE100 / PE150 sequencing run requires cPAS Barcode Primer 3 Reagent Kit (Cat. 1000020834), in addition to its corresponding sequencing set.
- A dual barcode App-A FCL PE100 / PE150 sequencing run requires High-Throughput Barcode Primer 3 Reagent Kit (App-A), in addition to its corresponding sequencing set.
- An App-D FCL PE150 sequencing run requires High-Throughput Sequencing Primer Kit (App-D), in addition to High-throughput Sequencing Set (FCL PE150).

Table 3 DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE35)
Catalog number: 1000020243

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit Catalog number 1000016115				
Low TE Buffer	960 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 µL×1 tube			
Make DNB Enzyme Mix I	800 µL×1 tube			
Make DNB Enzyme Mix II (LC)	80 µL×1 tube			
Stop DNB Reaction Buffer	400 µL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit Catalog number: 1000016114				
DNB Load Buffer I	300 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	150 µL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL SE35) Catalog number: 1000019813				
dNTPs Mix II	4.50 mL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix IV	1.70 mL×1 tube			
Sequencing Enzyme Mix	3.20 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL SE35) Catalog number: 1000019812				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 4 DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE50)
Catalog number: 1000016102

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit Catalog number: 1000016115				
Low TE Buffer	960 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 μL×1 tube			
Make DNB Enzyme Mix I	800 μL×1 tube			
Make DNB Enzyme Mix II (LC)	80 μL×1 tube			
Stop DNB Reaction Buffer	400 μL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit Catalog number: 1000016114				
DNB Load Buffer I	300 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	150 μL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL SE50) Catalog number: 1000016108				
dNTPs Mix II	2.70 mL×2 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix IV	2.00 mL×1 tube			
Sequencing Enzyme Mix	3.80 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL SE50) Catalog number: 1000016117				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 5 DNBSEQ-T7RS High-throughput Sequencing Set (FCL SE100)
Catalog number: 1000016103

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit Catalog number: 1000016115				
Low TE Buffer	960 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 μL×1 tube			
Make DNB Enzyme Mix I	800 μL×1 tube			
Make DNB Enzyme Mix II (LC)	80 μL×1 tube			
Stop DNB Reaction Buffer	400 μL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit Catalog number: 1000016114				
DNB Load Buffer I	300 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	150 μL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL SE100) Catalog number: 1000016109				
dNTPs Mix II	4.05 mL×2 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix IV	3.00 mL×1 tube			
Sequencing Enzyme Mix	2.85 mL×2 tubes			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL SE100) Catalog number: 1000016118				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 6 DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE100)
Catalog number: 1000028455

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit Catalog number: 1000016115				
Low TE Buffer	960 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 µL×1 tube			
Make DNB Enzyme Mix I	800 µL×1 tube			
Make DNB Enzyme Mix II (LC)	80 µL×1 tube			
Stop DNB Reaction Buffer	400 µL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit Catalog number: 1000028452				
DNB Load Buffer I	300 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	150 µL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate (T7 FCL)	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL PE100) Catalog number: 1000028450				
dNTPs Mix II	3.40 mL×3 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix V	3.90 mL×1 tube			
Sequencing Enzyme Mix	3.65 mL×2 tubes			
MDA Reagent	4.20 mL×1 tube			
MDA Enzyme Mix	0.60 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL PE100) Catalog number: 940-000004-00				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 7 DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE150)
Catalog number: 1000028454

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Rapid Make Reagent Kit Catalog number: 1000028453				
Low TE Buffer	960 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 μL×1 tube			
Make DNB Rapid Enzyme Mix II	800 μL×1 tube			
Make DNB Enzyme Mix II (LC)	80 μL×1 tube			
Stop DNB Reaction Buffer	400 μL×1 tube			
DNBSEQ-T7RS DNB Rapid Load Reagent Kit Catalog number: 1000028451				
DNB Load Buffer IV	200 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Micro Tube 0.5 mL (Empty)	1 tube			
Rapid Post Load Plate (T7 FCL PE150) V2.0	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL PE150) Catalog number: 1000028449				
dNTPs Mix II	6.90 mL×2 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix V	5.30 mL×1 tube			
Sequencing Enzyme Mix	4.95 mL×2 tubes			
MDA Reagent	4.20 mL×1 tube			
MDA Enzyme Mix	0.60 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL PE150) Catalog number: 940-000006-00				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 8 DNBSEQ-T7RS High-throughput Sequencing Set (stLFR FCL PE100)
Catalog number: 1000019251

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit (stLFR) Catalog number: 1000019257				
Low TE Buffer	480 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
stLFR Make DNB Buffer	160 µL×1 tube			
Make DNB Enzyme Mix III	320 µL×1 tube			
Make DNB Enzyme Mix IV	42 µL×1 tube			
Stop DNB Reaction Buffer	200 µL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit (stLFR) Catalog number: 1000019256				
DNB Load Buffer I	500 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	500 µL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate (stLFR)	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (stLFR FCL PE100) Catalog number: 1000019252				
dNTPs Mix II	4.90 mL×3 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix IV	5.40 mL×1 tube			
Sequencing Enzyme Mix	5.15 mL×2 tubes			
MDA Reagent	4.20 mL×1 tube			
MDA Enzyme Mix	0.60 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (stLFR FCL PE100) Catalog number: 1000019254				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Table 9 DNBSEQ-T7RS High-throughput Sequencing Set (App-A FCL PE100)
Catalog number: 940-000005-00

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Make Reagent Kit Catalog number: 1000016115				
Low TE Buffer	960 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 μL×1 tube			
Make DNB Enzyme Mix I	800 μL×1 tube			
Make DNB Enzyme Mix II (LC)	80 μL×1 tube			
Stop DNB Reaction Buffer	400 μL×1 tube			
DNBSEQ-T7RS DNB Load Reagent Kit Catalog number: 1000028452				
DNB Load Buffer I	300 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
DNB Load Buffer II	150 μL×1 tube			
Micro Tube 0.5 mL (Empty)	1 tube			
Post Load Plate (T7 FCL)	1 EA			
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL PE100) Catalog number: 1000028450				
dNTPs Mix II	3.40 mL×3 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix V	3.90 mL×1 tube			
Sequencing Enzyme Mix	3.65 mL×2 tubes			
MDA Reagent	4.20 mL×1 tube			
MDA Enzyme Mix	0.60 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL PE100) Catalog number: 940-000004-00				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
High-Throughput Pair-End Sequencing Primer Kit (App-A) Catalog number: 1000020832				
App-A Make DNB Buffer	400 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
1 µM App-A Insert Primer 1	2.20 mL×1 tube			
1 µM App-A Insert Primer 2	4.20 mL×1 tube			
1 µM App-A MDA Primer	4.20 mL×1 tube			
1 µM App-A Barcode Primer 2	3.50 mL×1 tube			

Table 10 DNBSEQ-T7RS High-throughput Sequencing Set (App-A FCL PE150) Catalog number: 940-000003-00

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS Sequencing Flow Cell Catalog number: 1000016269				
Sequencing Flow Cell (T7 FCL)	1 EA	0 °C to 30 °C	0 °C to 30 °C	10 months
DNBSEQ DNB Rapid Make Reagent Kit Catalog number: 1000028453				
Low TE Buffer	960 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Make DNB Buffer	400 µL×1 tube			
Make DNB Rapid Enzyme Mix II	800 µL×1 tube			
Make DNB Enzyme Mix II (LC)	80 µL×1 tube			
Stop DNB Reaction Buffer	400 µL×1 tube			
DNBSEQ-T7RS DNB Rapid Load Reagent Kit Catalog number: 1000028451				
DNB Load Buffer IV	200 µL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
Micro Tube 0.5 mL (Empty)	1 tube			
Rapid Post Load Plate (T7 FCL PE150)	1 EA			

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
DNBSEQ-T7RS High-throughput Sequencing Kit (FCL PE150) Catalog number: 1000028449				
dNTPs Mix II	6.90 mL×2 tubes	-80 °C to -15 °C	-25 °C to -15 °C	12 months
dNTPs Mix V	5.30 mL×1 tube			
Sequencing Enzyme Mix	4.95 mL×2 tubes			
MDA Reagent	4.20 mL×1 tube			
MDA Enzyme Mix	0.60 mL×1 tube			
Sequencing Reagent Cartridge	1 EA			
Transparent Sealing film	2 sheets			
DNBSEQ-T7RS Cleaning Reagent Kit (FCL PE150) Catalog number: 940-000006-00				
Washing Cartridge	1 EA	below 40 °C	0 °C to 30 °C	12 months
High-Throughput Pair-End Sequencing Primer Kit (App-A) Catalog number: 1000020832				
App-A Make DNB Buffer	400 μL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
1 μM App-A Insert Primer 1	2.20 mL×1 tube			
1 μM App-A Insert Primer 2	4.20 mL×1 tube			
1 μM App-A MDA Primer	4.20 mL×1 tube			
1 μM App-A Barcode Primer 2	3.50 mL×1 tube			

Table 11 CPAS Barcode Primer 3 Reagent Kit
Catalog number: 1000020834


Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
Primer for dual barcode sequencing (Pair End Sequencing use only)				
1 µM AD153 Barcode Primer 3	3.50 mL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months

Table 12 High-Throughput Barcode Primer 3 Reagent Kit (App-A)
Catalog number: 1000014047

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
Primer for dual barcode sequencing (Pair End Sequencing use only)				
1 μ M App-A Barcode Primer 3	3.50 mL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months

Table 13 High-Throughput Sequencing Primer Kit (App-D)
Catalog number: 1000028550

Component	Spec & quantity	Transportation temperature	Storage temperature	Expiry date
Primer for dual barcode sequencing (FCL PE150 Sequencing use)				
1 μ M App-D Insert Primer 1	2.20 mL×1 tube	-80 °C to -15 °C	-25 °C to -15 °C	12 months
1 μ M App-D MDA Primer	4.20 mL×1 tube			
1 μ M App-D Insert Primer 2	4.20 mL×1 tube			
1 μ M App-D Barcode Primer 2	3.50 mL×1 tube			
1 μ M App-D Barcode Primer 3	3.50 mL×1 tube			
App-D Make DNB Buffer	400 μ L×1 tube			

 **NOTE** The kit supports the sequencing of mixed libraries with Truseq, Nextera and AD153 adapters and should be used together with DNBSEQ-T7RS High-throughput Sequencing Set (FCL PE150).

2.2 Self-prepared equipment and consumables

Table 14 Self-prepared equipment and consumables

Equipment and consumables	Recommended brand	Catalog number
Qubit 3.0 fluorometer	Thermo Fisher	Q33216
Thermal cycler	Bio-Rad	/
MPC2000 96-well plate centrifuge	/	/
Pipette	Eppendorf	/

Equipment and consumables	Recommended brand	Catalog number
Electronic pipette	Labnet	FASTPETTEV-2
Mini centrifuge	Major Laboratory Supplier (MLS)	/
Vortex mixer	MLS	/
2 °C to 8 °C Refrigerator	MLS	/
-25 °C to -15 °C Freezer	MLS	/
Qubit ssDNA assay kit	Thermo Fisher	Q10212
2 M NaOH solution	Aladdin	S128511-1L
100% Tween-20	BBI	A600560-0500
5 M NaCl solution	SIGMA	S5150-4L
75% Ethanol	/	/
Power dust remover	MATIN	M-6318
Sterile pipette tip (box)	AXYGEN	/
5 mL Sterile pipette tip (box)	AXYGEN	/
200 µL Wide-bore pipette tips	AXYGEN	T-205-WB-C
Qubit assay tubes	Thermo Fisher	Q32856
0.2 mL PCR 8-tube strip	AXYGEN	/
1.5 mL microcentrifuge tube	AXYGEN	MCT-150-C
Ice box	MLS	/
Ice machine	/	/
100 mL Serological pipet	CORNING	4491
25 mL Serological pipet	CORNING	4489
10 mL Serological pipet	CORNING	4488
15 mL Sterile tube	SARSTEDT	60.732.001
Microfiber clean wiper	DUSTFREE TECHNOLOGY CO.,LTD	LJ618180B1
5 mL Transport tubes	AXYGEN	/
Lint-free paper	MLS	/
Ziplock bag	/	/

Chapter 3 Sequencing workflow

1	Making DNB: use reagents from DNB Making Kit and DNA library from user to make DNB
2	Loading DNB: load DNB into the flow cell using reagents from DNB Load Reagent Kit at MGIDL-T7RS loader
3	Preparing sequencing reagent kit: inspect, thaw the reagent kit and then add the required reagents, as well as check the pure water container and waste liquid container
4	Sequencing
5	Data analysis


Chapter 4 Making DNB

4.1 Insert size recommendation

These sequencing sets are compatible with the libraries prepared by using MGI Library Prep Kits. For general purpose, library refers to single stranded circular DNA (ssDNA). For the best sequencing quality, it is recommended that the insert size of the library should be between 50-500 nucleotides (nt), and the main band is centered within ± 100 nt. For the stLFR library prepared with MGIEasy stLFR Library Prep Kit, the library is circularized dsDNA. It is recommended that the insert size of the library be between 200-1500 nt. If there are special requirements or specifications from the library preparation kit, then the requirements of the kit should be followed.

Table 15 Recommended insert size and theoretical throughput for each flow cell

Sequencing kit	Suggested insert distribution (bp)	Applications	Mappable reads (M)	Data output (Gb)
FCL SE35	50 - 230	NIPT	5000	175
FCL SE50	50 - 230	NIPT, PMSEQ	5000	250
FCL SE100	200 - 400	PMSEQ	5000	500
FCL PE100	200 - 400	WGS, WES, RNAseq, Single Cell	5000	1000
FCL PE150	300 - 500	WGS, WES, RNAseq	5000	1500
App-A FCL PE100	200 - 400	WGS, WES, RNAseq	5000	1000
App-A FCL PE150	300 - 500	WGS, WES, RNAseq	5000	1500
App-D FCL PE150	300 - 500	WGS, WES, RNAseq	4000	1200
stLFR FCL PE100	200 - 1500	stLFR	4000	800

-  **NOTE**
- Consider the insert size and required data output when selecting sequencing kits.
 - Average data output will vary with different library type and applications.

4.2 Library concentration and amount requirement

Table 16 Recommended insert size and theoretical throughput for each flow cell

Libraries	Library concentration
general libraries	ssDNA library concentration ≥ 2 fmol/ μ L
PCR-free	ssDNA library concentration ≥ 3.75 fmol/ μ L
stLFR	dsDNA library concentration ≥ 1.9 ng/ μ L

NOTE • If the library concentration is unknown, it is recommended to perform ssDNA library quantitation (ng/ μ L) using the Qubit ssDNA Assay Kit and the Qubit Fluorometer. Use the following equation to convert the concentration of the ssDNA library from ng/ μ L to fmol/ μ L:

$$\text{Concentration (fmol}/\mu\text{L)} = 3030 \times \text{Concentration (ng}/\mu\text{L)} / N$$

- If there are special requirements or specifications of the library preparation kit, then the requirements of the kit should be followed.

4.3 Library pooling

4.3.1 Number of samples that can be pooled together

The DNBSEQ-T7RS sequencer can simultaneously perform sequencing of 4 flow cells and each flow cell can theoretically produce 5000M reads. For PE100 sequencing, one flow cell can produce 1 Tb of data in theory. The number of samples that can be pooled together for each flow cell depends on the required data output, read length, and specific application.

As a guide, do not pool more samples with their total data output larger than 90% of the theoretical data output as described in *Table 15 on Page 17*, due to variation in pooling and the fact that not all barcodes will generate the same amount of the data output from the same amount of DNB.

$$\text{Maximum number of samples pooled} = \frac{\text{Total data output of one flow cell} \times (1 - \text{poolingCV})}{\text{required data per sample} \times 2}$$

Example 1: Human Whole-genome Sequencing (WGS).

When using the PE100 sequencing kit, 9 samples on each flow cell are recommended.

Example 2: stLFR sample.

When using the PE100 sequencing kit, if the required sequencing depth is 40X, then 6 samples are recommended to be pooled for each flow cell.

Example 3: 50G are required for each sample.

When using the PE100 sequencing kit, if 50G are required for each sample, then 18 samples are recommended to be pooled for each flow cell.

Example 4: Pooling samples with various applications.

When using the PE150 sequencing kit, if samples to be sequenced include WGS (100G/sample) and RNASeq (50G/sample), it is recommended to pool 4 WGS samples and 18 RNASeq samples for each flow cell.

Table 17 Examples of various sample pooling

Index	Read length	Minimum data for each sample	Pooling sample number	Theoretical data output range for each sample
1	PE100	100 Gb	9	100 - 122 Gb
2	stLFR PE100	120 Gb	6	120 - 146 Gb
3	PE100	50 Gb	18	50 - 61 Gb
4	PE150	50 Gb	18 RNAseq	52 - 64 Gb
		100 Gb	4 WGS	102 - 122 Gb

 **NOTE** Assuming pooling variation are within $\pm 10\%$.

4.3.2 Verifying the base balance for barcode

- A balanced base composition in each sequencing cycle is very important for high sequencing quality. It is strongly recommended that the minimum base composition of A, C, G, T for each position in the barcode is not lower than 12.5%. For a given pooling of samples, if the minimum base composition of A, C, G, T within the barcode is between 5% and 12.5%, the barcode split rate may be compromised. If the minimum base composition of A, C, G, T in any position of the barcode is less than 5%, it is strongly suggested to re-design the pooling strategy to have a more balanced base composition in the barcode.
- It is also important to note that two or more samples with an identical barcode should not be pooled together, otherwise, it is impossible to assign the read correctly.

4.4 Making DNB

Six DNB making protocols are listed below, please select the appropriate one according to the sequencing kit used.

- Section 4.4.1 *Preparing DNB for the FCL SE35, FCL SE50, FCL SE100 and FCL PE100 on Page 20.*
- Section 4.4.2 *Preparing DNB for the FCL PE150 kit on Page 24.*
- Section 4.4.3 *Preparing DNB for the App-A FCL PE100 kit on Page 27.*
- Section 4.4.4 *Preparing DNB for the App-A FCL PE150 kit on Page 31.*
- Section 4.4.5 *Preparing DNB for the App-D FCL PE150 kit on Page 34.*
- Section 4.4.6 *Preparing DNB for the stLFR FCL PE100 kit on Page 37.*

4.4.1 Preparing DNB for the FCL SE35, FCL SE50, FCL SE100 and FCL PE100

4.4.1.1 Calculating the required amount of ssDNA library

- 270 μ L of DNB is required to load one flow cell for the FCL SE35, FCL SE50, FCL SE100 and FCL PE100.


- One DNB making reaction can make either 100 μL or 50 μL of DNB. The volume of the DNB making reaction system depends on the amount of data required for sequencing per sample and the types of DNA libraries.
- The required ssDNA library volume to make either 100 μL or 50 μL of DNB are shown in the table below.

Table 18 FCL SE35/SE50/SE100/PE100 required ssDNA volume

Sample type	Required ssDNA volume: V (μL)	
	100 μL DNB reaction	50 μL DNB reaction
Regular library	$V=40 \text{ fmol} / C$	$V=20 \text{ fmol} / C$
PCR free	$V=75 \text{ fmol} / C$	$V=37.5 \text{ fmol} / C$

 **NOTE** C refers to the concentration of the ssDNA library ($\text{fmol} / \mu\text{L}$).

- For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V=a/b \times 270$ (μL).
- If the total sample number pooled is <6 , it is suggested to select the volume of 100 μL DNB reaction, and the number of 100 μL DNB making reactions is equal to round $(V/100)+1$. (for example: If $V=80$, it requires one 100 μL DNB making reaction; If $V=120$, it requires two 100 μL DNB making reactions)
- If the total sample number pooled is ≥ 6 , it is suggested to select the volume of 50 μL DNB reaction, and the number of 50 μL DNB making reactions is equal to round $(V/50)+1$.

 **NOTE** All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.1.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.
2. Take out the Make DNB Buffer, Low TE Buffer and Stop DNB Reaction Buffer from the DNB Make Reagent Kit and thaw reagents at room temperature.
3. Thaw Make DNB Enzyme Mix I on ice for approximately 0.5 hours.

- After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.1.3 Making DNB

Perform the steps below:

- Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.1.1 *Calculating the required amount of ssDNA library on Page 20*.

Table 19 FCL SE35/SE50/SE100/PE100 make DNB reaction mixture 1

Component	Volume / 100 μ L DNB reaction (μ L)	Volume / 50 μ L DNB reaction (μ L)
Low TE Buffer	20 - V	10 - V
Make DNB Buffer	20	10
ssDNA libraries	V	V
Total volume	40	20

 **NOTE** V is the volume of ssDNA library determined from *Table 18 on Page 21* for the 100 μ L and 50 μ L DNB reaction, respectively.

- Mix the make DNB reaction mixture 1 thoroughly by vortex and centrifuge for 5 seconds by using a mini centrifuge.
- Place the mixture into a thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 20 FCL SE35/SE50/SE100/PE100 primer hybridization reaction condition

Temperature	Time
Heated lid (105 $^{\circ}$ C)	On
95 $^{\circ}$ C	1 min
65 $^{\circ}$ C	1 min
40 $^{\circ}$ C	1 min
4 $^{\circ}$ C	Hold

NOTE As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

- Take out the Make DNB Enzyme Mix II (LC) from storage and place it on ice. Centrifuge briefly for 5 seconds and place it on ice until use. Mix the Make DNB Enzyme Mix II (LC) with Make DNB Enzyme Mix I by gently pipetting 6 to 8 times before use as described in the table below.

NOTE

- Do not place Make DNB Enzyme Mix II (LC) at room temperature.
- Avoid holding the tube for a prolonged time.

Table 21 FCL SE35/SE50/SE100/PE100 make DNB reaction mixture 2

Component	Volume / 100 μ L DNB reaction (μ L)	Volume / 50 μ L DNB reaction (μ L)
Make DNB Enzyme Mix I	40	20
Make DNB Enzyme Mix II (LC)	4	2

- Take the PCR tube out of the thermal cycler when the temperature reaches 4 °C . Centrifuge briefly for 5 seconds, then place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.
- Mix the mixture thoroughly by vortex, centrifuge for 5 seconds by using mini centrifuge.
- Place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 22 FCL SE35/SE50/SE100/PE100 rolling circle amplification condition

Temperature	Time
Heated lid (35 °C)	On
30 °C	25 min
4 °C	Hold

NOTE It is recommended to set the temperature of the heated lid to 35 °C or as close as possible to 35 °C .

8. Add 20 μL (for 100 μL DNB reaction) or 10 μL (for 50 μL DNB reaction) Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 $^{\circ}\text{C}$. Mix gently by pipetting 5 to 8 times using a wide bore tip.



NOTE

- It is very important to mix DNBs gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
- Store the DNBs at 4 $^{\circ}\text{C}$ and perform sequencing within 48 hours.

4.4.2 Preparing DNB for the FCL PE150 kit

4.4.2.1 Calculating the required amount of ssDNA library

300 μL of DNB is required to load one flow cell for the FCL PE150 kit. One DNB making reaction can make 90 μL of DNB. The required ssDNA library volume needed to make 90 μL of DNB (one DNB reaction) is shown in the table below.

Table 23 FCL PE150 required ssDNA volume

Sample type	Required ssDNA volume: V (μL)
Regular Library	$V = 40 \text{ fmol} / C$
PCR free	$V = 75 \text{ fmol} / C$



NOTE C refers to the concentration of the ssDNA library ($\text{fmol} / \mu\text{L}$).

For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V = a / b \times 300$ (μL).

The number of the 90 μL DNB making reactions is equal to round $(V / 90) + 1$.



NOTE All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.2.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.

2. Take out the Make DNB Buffer, Low TE Buffer and Stop DNB Reaction Buffer from the DNB Rapid Make Reagent Kit and thaw reagents at room temperature.
3. Thaw Make DNB Rapid Enzyme Mix II on ice for approximately 0.5 hours.
4. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.2.3 Making DNB

Perform the steps below:

1. Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.2.1 *Calculating the required amount of ssDNA library on Page 24*.

Table 24 FCL PE150 make DNB reaction mixture 1

Component	Volume / 90 μ L DNB reaction (μ L)
Low TE Buffer	20 - V
Make DNB Buffer	20
ssDNA libraries	V
Total volume	40

 **NOTE**


- V is the volume of ssDNA library determined from *Table 23 on Page 24* for the 90 μ L DNB reaction.
- Do not throw away the Low TE Buffer after this step. The remaining Low TE buffer will be used in DNB loading.

2. Mix the Make DNB reaction mixture 1 thoroughly by vortexing and centrifuge for 5 seconds by using a mini centrifuge.


- Place the mix into a thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 25 FCL PE150 primer hybridization reaction condition

Temperature	Time
Heated lid (105 °C)	On
95 °C	1 min
65 °C	1 min
40 °C	1 min
4 °C	Hold

 **NOTE** As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

- Take out the Make DNB Enzyme Mix II (LC) from storage and place on ice. Centrifuge briefly for 5 seconds and place on ice until use.
- Mix the Make DNB Enzyme Mix II (LC) with Make DNB Rapid Enzyme Mix II by gently pipetting 6 to 8 times before use as described in the table below.

 **NOTE**

- Do not place Make DNB Enzyme Mix II (LC) at room temperature.
- Avoid holding the tube for a prolonged time.


Table 26 FCL PE150 make DNB reaction mixture 2

Component	Volume / 90 µL DNB reaction (µL)
Make DNB Rapid Enzyme Mix II	40
Make DNB Enzyme Mix II (LC)	1.6


- Take the PCR tube out of the thermal cycler when the temperature reaches 4 °C . Centrifuge briefly for 5 seconds.
- Place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.
- Mix thoroughly by vortexing, centrifuge for 5 seconds by using a mini centrifuge and place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 27 FCL PE150 rolling circle amplification condition

Temperature	Time
Heated lid (35 °C)	On
30 °C	10 min
4 °C	Hold

 **NOTE** It is recommended to set the temperature of the heated lid to 35 °C or as close as possible to 35 °C .

- Add 10 µL Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 °C. Mix gently by pipetting 5 to 8 times using a wide bore tip.

 **NOTE**

- It is very important to mix DNB gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
- Store the DNB at 4 °C and perform sequencing within 8 hours.

4.4.3 Preparing DNB for the App-A FCL PE100 kit

4.4.3.1 Calculating the required amount of ssDNA library

- 270 µL of DNB is required to load one flow cell for the App-A FCL PE100 kit. One DNB making reaction can make either 100 µL or 50 µL of DNB. The required ssDNA library volume to make either 100 µL or 50 µL of DNB is shown in the table below.


Table 28 App-A FCL PE100 required ssDNA volume

Sample type	Required ssDNA volume: V (µL)	
	100 µL DNB reaction	50 µL DNB reaction
Regular library	$V = 60 \text{ fmol} / C$	$V = 30 \text{ fmol} / C$
PCR free	$V = 75 \text{ fmol} / C$	$V = 37.5 \text{ fmol} / C$

 **NOTE** C refers to the concentration of the ssDNA library (fmol/µL).

- For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V = a / b \times 270$ (µL).


- If the total sample number pooled is <6 , it is suggested to select the volume of 100 μL DNB reaction, and the number of 100 μL DNB making reactions is equal to $\text{round}(V/100)+1$. (For example: If $V=80$, it requires one 100 μL DNB making reaction; If $V=120$, it requires two 100 μL DNB making reactions)
- If the total sample number pooled is ≥ 6 , it is suggested to select the volume of 50 μL DNB reaction, and the number of 50 μL DNB making reactions is equal to $\text{round}(V/50)+1$.

 **NOTE** All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.3.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.
2. Take out the App-A Make DNB Buffer from High-Throughput Pair-End Sequencing Primer Kit (App-A) and thaw reagents at room temperature.
3. Take out the Low TE Buffer and Stop DNB Reaction Buffer from the DNB Make Reagent Kit and thaw reagents at room temperature.
4. Thaw Make DNB Enzyme Mix I on ice for approximately 0.5 hours.
5. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.3.3 Making DNB

Perform the steps below:

1. Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.3.1 *Calculating the required amount of ssDNA library on Page 27*.

Table 29 App-A FCL PE100 make DNB reaction mixture 1


Component	Volume / 100 μ L DNB reaction (μ L)	Volume / 50 μ L DNB reaction (μ L)
Low TE Buffer	20 - V	10 - V
App-A Make DNB Buffer	20	10
ssDNA libraries	V	V
Total volume	40	20

 **NOTE** V is the volume of ssDNA library determined from *Table 28* on *Page 27* for the 100 μ L and 50 μ L DNB reaction, respectively.

- Mix the Make DNB reaction mixture 1 thoroughly by vortex and centrifuge for 5 seconds by using a mini centrifuge.
- Place the mixture into a thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 30 App-A FCL PE100 primer hybridization reaction condition

Temperature	Time
Heated lid (105 $^{\circ}$ C)	On
95 $^{\circ}$ C	1 min
65 $^{\circ}$ C	1 min
40 $^{\circ}$ C	1 min
4 $^{\circ}$ C	Hold

 **NOTE** As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

- Take out the Make DNB Enzyme Mix II (LC) from storage and place on ice. Centrifuge briefly for 5 seconds and place on ice until use.

- Mix the Make DNB Enzyme Mix II (LC) with Make DNB Enzyme Mix I by gently pipetting 6 to 8 times before use as described in the table below.


-  **NOTE**
- Do not place Make DNB Enzyme Mix II (LC) at room temperature
 - Avoid holding the tube for a prolonged time.


Table 31 App-A FCL PE100 make DNB reaction mixture 2

Component	Volume / 100 μ L DNB reaction (μ L)	Volume / 50 μ L DNB reaction (μ L)
Make DNB Enzyme Mix I	40	20
Make DNB Enzyme Mix II (LC)	4	2


- Take the PCR tube out of the thermal cycler when the temperature reaches 4 °C . Centrifuge briefly for 5 seconds.
- Place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.
- Mix thoroughly by vortexing, centrifuge for 5 seconds by using mini centrifuge.
- Place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 32 App-A FCL PE100 rolling circle amplification condition

Temperature	Time
Heated lid (35 °C)	On
30 °C	25 min
4 °C	Hold

-  **NOTE** It is recommended to set the temperature of the heated lid to 35 °C or as close as possible to 35 °C .

- Add 20 μ L (for a 100 μ L DNB reaction) or 10 μ L (for a 50 μ L DNB reaction) of Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 °C. Mix gently by pipetting 5 to 8 times using a wide bore tip.

-  **NOTE**
- It is very important to mix DNBs gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
 - Store the DNBs at 4 °C and perform sequencing within 48 hours.

4.4.4 Preparing DNB for the App-A FCL PE150 kit

4.4.4.1 Calculating the required amount of ssDNA library

- 300 μL of DNB is required to load one flow cell for the FCL PE150 kit. One DNB making reaction can make 90 μL of DNBs. The required ssDNA library volume needed to make 90 μL of DNBs (one DNB reaction) is shown in the table below.


Table 33 App-A FCL PE150 required ssDNA volume

Sample type	Required ssDNA volume: V (μL)
Regular Library	$V = 60 \text{ fmol} / C$
PCR free	$V = 75 \text{ fmol} / C$

 **NOTE** C refers to the concentration of the ssDNA library ($\text{fmol} / \mu\text{L}$).

- For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V = a / b \times 300$ (μL).

The number of the 90 μL DNB making reactions is equal to round $(V / 90) + 1$.

 **NOTE** All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.4.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.
2. Take out the App-A Make DNB Buffer from the High Throughput Pair-End Sequencing Primer Kit (App-A).
3. Take out the Low TE Buffer and Stop DNB Reaction Buffer from the DNB Rapid Make Reagent Kit and thaw reagents at room temperature.
4. Thaw Make DNB Rapid Enzyme Mix II on ice for approximately 0.5 hours.

- After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place it on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.4.3 Making DNB

Perform the steps below:

- Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.4.1 *Calculating the required amount of ssDNA library on Page 31*.

Table 34 App-A FCL PE150 make DNB reaction mixture 1

Component	Volume / 90 μ L DNB reaction (μ L)
App-A Make DNB Buffer	20
Low TE Buffer	20 - V
ssDNA libraries	V
Total volume	40

-  **NOTE**
- V is the volume of ssDNA library determined from *Table 33 on Page 31* for the 90 μ L DNB reaction.
 - Do not throw away the Low TE Buffer after this step. The remaining Low TE buffer will be used in DNB loading.

- Mix the Make DNB reaction mixture 1 thoroughly by vortex and centrifuge for 5 seconds by using a mini centrifuge. Place the mix into a thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 35 App-A FCL PE150 primer hybridization reaction condition

Temperature	Time
Heated lid (105 $^{\circ}$ C)	On
95 $^{\circ}$ C	1 min
65 $^{\circ}$ C	1 min
40 $^{\circ}$ C	1 min
4 $^{\circ}$ C	Hold

NOTE As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

3. Take out the Make DNB Enzyme Mix II (LC) from storage and place on ice. Centrifuge briefly for 5 seconds and place on ice until use.
4. Mix the Make DNB Enzyme Mix II (LC) with Make DNB Rapid Enzyme Mix II by gently pipetting 6 to 8 times before use as described in the table below.

NOTE

- Do not place Make DNB Enzyme Mix II (LC) at room temperature.
- Avoid holding the tube for a prolonged time.

Table 36 App-A FCL PE150 make DNB reaction mixture 2

Component	Volume / 90 μ L DNB reaction (μ L)
Make DNB Rapid Enzyme Mix II	40
Make DNB Enzyme Mix II (LC)	1.6


5. Take the PCR tube out of the thermal cycler when the temperature reaches 4 $^{\circ}$ C . Centrifuge briefly for 5 seconds.
6. Place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.
7. Mix thoroughly by vortex, centrifuge for 5 seconds by using mini centrifuge and place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 37 App-A FCL PE150 rolling circle amplification condition

Temperature	Time
Heated lid (35 $^{\circ}$ C)	On
30 $^{\circ}$ C	10 min
4 $^{\circ}$ C	Hold

NOTE It is recommended to set the temperature of the heated lid to 35 $^{\circ}$ C or as close as possible to 35 $^{\circ}$ C .

8. Add 10 μL Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 $^{\circ}\text{C}$. Mix gently by pipetting 5 to 8 times using a wide bore tip.

-  **NOTE**
- It is very important to mix DNB gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
 - Store the DNB at 4 $^{\circ}\text{C}$ and perform sequencing within 8 hours.

4.4.5 Preparing DNB for the App-D FCL PE150 kit

4.4.5.1 Calculating the required amount of ssDNA library

- 300 μL of DNB is required to load one flow cell for the FCL PE150 kit. One DNB making reaction can make 90 μL of DNB. The required ssDNA library volume needed to make 90 μL of DNB (one DNB reaction) is shown in the table below.


Table 38 App-D FCL PE150 required ssDNA volume

Sample type	Required ssDNA volume: V (μL)
Regular Library	$V=60 \text{ fmol} / C$
PCR free	$V=75 \text{ fmol} / C$

-  **NOTE** C refers to the concentration of the ssDNA library ($\text{fmol} / \mu\text{L}$).

- For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V=a/b \times 300$ (μL).

The number of the 90 μL DNB making reactions is equal to round $(V/90)+1$.

-  **NOTE** All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.5.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.
2. Take out the App-D Make DNB Buffer from the High Throughput Pair-End Sequencing Primer Kit (App-D).
3. Take out the Low TE Buffer and Stop DNB Reaction Buffer from the DNB Rapid Make Reagent Kit and thaw reagents at room temperature.
4. Thaw Make DNB Rapid Enzyme Mix II on ice for approximately 0.5 hours.
5. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.5.3 Making DNB

Perform the steps below:

1. Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.5.1 *Calculating the required amount of ssDNA library on Page 34*.

Table 39 App-D FCL PE150 make DNB reaction mixture 1

Component	Volume (μL) for 90 μL DNB reaction
App-D Make DNB Buffer	20
Low TE Buffer	20-V
ssDNA libraries	V
Total volume	40


 **NOTE**

- V is the volume of ssDNA library determined from *Table 38 on Page 34* for the 90 μL DNB reaction.
- Do not throw away the Low TE Buffer after this step. The remaining Low TE buffer will be used in DNB loading.


- Mix the Make DNB reaction mixture 1 thoroughly by vortex and centrifuge for 5 seconds by using a mini centrifuge. Place the mix into a Thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 40 App-D FCL PE150 primer hybridization reaction condition

Temperature	Time
Heated lid (105 °C)	On
95 °C	1 min
65 °C	1 min
40 °C	1 min
4 °C	Hold

 **NOTE** As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

- Take out the Make DNB Enzyme Mix II (LC) from storage and place on ice. Centrifuge briefly for 5 seconds and place on ice until use.
- Mix the Make DNB Enzyme Mix II (LC) with Make DNB Rapid Enzyme Mix II by gently pipetting 6 to 8 times before use as described in the table below.

 **NOTE**

- Do not place Make DNB Enzyme Mix II (LC) at room temperature.
- Avoid holding the tube for a prolonged time.

Table 41 App-D FCL PE150 make DNB reaction mixture 2


Component	Volume (μL) for 90 μL DNB reaction
Make DNB Rapid Enzyme Mix II	40
Make DNB Enzyme Mix II (LC)	1.6

- Take the PCR tube out of the thermal cycler when the temperature reaches 4 °C . Centrifuge briefly for 5 seconds.
- Place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.


- Mix thoroughly by vortex, centrifuge for 5 seconds by using mini centrifuge and place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 42 App-D FCL PE150 rolling circle amplification condition

Temperature	Time
Heated lid (35 °C)	On
30 °C	10 min
4 °C	Hold

 **NOTE** It is recommended to set the temperature of the heated lid to 35 °C or as close as possible to 35 °C .

- Add 10 µL Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 °C. Mix gently by pipetting 5 to 8 times using a wide bore tip.

 **NOTE**

- It is very important to mix DNB gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
- Store the DNB at 4 °C and perform sequencing within 8 hours.

4.4.6 Preparing DNB for the stLFR FCL PE100 kit


4.4.6.1 Calculating the required amount of dsDNA library

- 270 µL of DNB is required to load one flow cell for the stLFR FCL PE100. One DNB making reaction can make 80 µL of DNB. 30 ng dsDNA libraries are needed to make 80 µL of DNB; Therefore, the volume of stLFR library needed for each 80 µL DNB preparation reaction is defined as follows: $V (\mu\text{L}) = 30 \text{ ng} / C$

 **NOTE** C refers to the concentration of stLFR dsDNA library (ng/µL).

- For a given sample A, if it requires “a” million base data output and the total theoretical expected data output for this flow cell is “b” million bases, then the required DNB volume (V) in the pooling for sample A is as follows: $V = a / b \times 270$ (μL).

The number of the 80 μL DNB making reactions is equal to round $(V / 80) + 1$.

 **NOTE** All samples should be considered potentially infectious and should be handled in accordance with relevant national regulations.

4.4.6.2 Preparing reagents for DNB making

Perform the steps below:

1. Place the library on ice until use.
2. Take out the stLFR Make DNB Buffer, Low TE Buffer and Stop DNB Reaction Buffer from the DNB Make Reagent Kit (stLFR) and thaw reagents at room temperature.
3. Thaw Make DNB Enzyme Mix III on ice for approximately 0.5 hours.
4. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

 **NOTE** Mixed use of reagent components from different batches is not recommended.


4.4.6.3 Making DNB

Perform the steps below:

1. Take a 0.2 mL 8-tube strip or PCR tubes and prepare the reaction mix on ice following the table below. The following table only illustrates one make DNB reaction. The required number of make DNB reactions is determined by the actual application as described in 4.4.6.1 *Calculating the required amount of dsDNA library on Page 37*.

Table 43 stLFR FCL PE100 make DNB reaction mixture 1


Component	Volume / 80 μL DNB reaction (μL)
stLFR Make DNB Buffer	16
Low TE Buffer	16 - V
dsDNA libraries	V
Total volume	32

 **NOTE** V is the volume of dsDNA library determined from 30/C. C refers to the concentration of stLFR dsDNA library in ng/μL.

2. Mix the Make DNB reaction mixture 1 thoroughly by vortex and centrifuge for 5 seconds by using a mini centrifuge.
3. Place the mix into a thermal cycler and start the primer hybridization reaction. Thermal cycler settings are described in the table below:

Table 44 stLFR FCL PE100 primer hybridization reaction condition

Temperature	Time
Heated lid (105 °C)	On
95 °C	3 min
40 °C	3 min
4 °C	Hold

 **NOTE** As some thermal cyclers are slow in temperature adjustment, when the lid is being heated or cooled, the sample block may remain at room temperature and the procedure is not performed. For these types of thermal cyclers, pre-heating of the heated lid is required to ensure that the heated lid is at working temperature during the DNB reaction.

4. Take out the Make DNB Enzyme Mix IV from storage and place on ice. Centrifuge briefly for 5 seconds and place on ice until use.
5. Mix the Make DNB Enzyme Mix IV with Make DNB Enzyme Mix III by gently pipetting 6 to 8 times before use as described in the table below.


 **NOTE** Do not place Make DNB Enzyme Mix IV at room temperature and avoid holding the tube for a prolonged time.

Table 45 stLFR FCL PE100 make DNB reaction mixture 2


Component	Volume / 80 μL DNB reaction (μL)
Make DNB Enzyme Mix III	32
Make DNB Enzyme Mix IV	3.2

6. Take the PCR tube out of the thermal cycler when the temperature reaches 4 °C . Centrifuge briefly for 5 seconds.
7. Place the tube on ice and add all the Make DNB Reaction Mixture 2 to the tube.


8. Mix thoroughly by vortex, centrifuge for 5 seconds by using mini centrifuge and place the tubes into a thermal cycler for the next reaction. The condition is shown in the table below:

Table 46 stLFR FCL PE100 rolling circle amplification condition

Temperature	Time
Heated lid (35 °C)	On
30 °C	30 min
4 °C	Hold

 **NOTE** It is recommended to set the temperature of the heated lid to 35 °C or as close as possible to 35 °C .

9. Add 16 µL Stop DNB Reaction Buffer to the tube immediately once the temperature reaches to 4 °C. Mix gently by pipetting 5 to 8 times using a wide bore tip.

 **NOTE**

- It is very important to mix DNB gently by using a wide bore pipette tip. Do not centrifuge, vortex, pipette vigorously or shake the tube.
- Store the DNB at 4 °C and perform sequencing within 48 hours.

4.5 Quantifying DNB and pooling

4.5.1 Quantifying DNB

- When the DNB making is completed, use the Qubit ssDNA Assay Kit and Qubit Fluorometer to quantify the DNBs.

Table 47 Requires the DNB concentration


Sequencing kit	Regular library	Other library
FCL PE100 and below	≥15 ng/µL	≥8 ng/µL
FCL PE150	≥10 ng/µL	≥5 ng/µL

 **NOTE** If the concentration is lower than 8 ng/µL or 5 ng/µL, make a new DNB.

- If the concentration exceeds 40 ng/µL, the DNB need to be diluted to 20 ng/µL before loading, please refer to the following table for DNB dilution.

Table 48 DNB dilution buffer

Sequencing kit	Dilution reagent	DNB storage temperature	Maximum DNB storage time (hour)
FCL SE35	DNB Load Buffer I	4 °C	≤48
FCL SE50	DNB Load Buffer I	4 °C	≤48
FCL SE100	DNB Load Buffer I	4 °C	≤48
FCL PE100	DNB Load Buffer I	4 °C	≤48
App-A FCL PE100	DNB Load Buffer I	4 °C	≤48
stLFR FCL PE100	DNB Load Buffer I	4 °C	≤48
FCL PE150	Low TE Buffer	4 °C	≤8
App-A FCL PE150	Low TE Buffer	4 °C	≤8
App-D FCL PE150	Low TE Buffer	4 °C	≤8

 **NOTE** Only Low TE Buffer can be used as DNB dilution reagent for FCL PE150, App-A FCL PE150 and App-D FCL PE150 kits.

4.5.2 DNB pooling

Amount of DNB (μL) needed for each sample in the pool depends on the relative amount for this sample and the total amount of DNB needed for loading one flow cell which is defined by the specific type of sequencing kit.

4.5.2.1 Calculating the relative amount for each sample

Assuming there are 8 samples (A to H) in the pool, the relative amount for each sample is defined as:

The relative amount of A sample (A1)=data output required for sample A / the concentration of DNB for sample A.

The relative amount of B sample (B1)=data output required for sample B / the concentration of DNB for sample B.

...

The relative amount of H sample (H1)=data output required for sample H / the concentration of DNB for sample H.

4.5.2.2 Calculating the total relative amount (V) for all sample

$$V = A1 + B1 + \dots + H1$$

4.5.2.3 Calculating the DNB volume needed for each sample

For each FCL flow cell used for SE35/SE50/SE100/PE100, App-A PE100 and stLFR PE100 requiring 270 μ L DNB, the DNB volume for pooling is calculated as follows:

DNB volume for sample A: $A2 = 270 \times A1 / V$

DNB volume for sample B: $B2 = 270 \times B1 / V$

...

DNB volume for sample H: $H2 = 270 \times H1 / V$


For each FCL flow cell used for PE150, App-A PE150 and App-D PE150 requiring 300 μ L DNB, the DNB volume for pooling is calculated as follows:

DNB volume for sample A: $A2 = 300 \times A1 / V$

DNB volume for sample B: $B2 = 300 \times B1 / V$

...

DNB volume for sample H: $H2 = 300 \times H1 / V$

 **NOTE** Use normal pipette tips to aspirate the required volume of each DNB and use wide bore tips to mix.

Chapter 5 Loading DNB

5.1 Preparing the post load plate and buffers

Two DNB loading protocols are listed below, please select the appropriate one depending on the sequencing kit used:

- Section 5.1.1 *Preparing the post load plate and buffer for FCL SE35/SE50/SE100/PE100, App-A FCL PE100, stLFR FCL PE100 sequencing on Page 43.*
- Section 5.1.2 *Post load plate and buffers preparation for FCL PE150, App-A FCL PE150 and App-D FCL PE150 sequencing on Page 440.*

5.1.1 Preparing the post load plate and buffer for FCL SE35/SE50/SE100/PE100, App-A FCL PE100, stLFR FCL PE100 sequencing

5.1.1.1 Thawing the post load plate

Perform the steps below:

1. Take out the Post Load Plate from the DNB Load Reagent Kit (For stLFR PE100, take out the Post Load Plate (stLFR) from the DNB Load Reagent Kit (stLFR)) and thaw in a water bath at room temperature for 2 hours.
2. Once the Post Load Plate is thoroughly thawed, place it in a 2 °C to 8 °C refrigerator until use. The Post Load Plate can also be thawed in 2 °C to 8 °C refrigerator overnight.
3. Gently invert the Post Load Plate to mix 5 times and then centrifuge for 1 minute before use.

5.1.1.2 Preparing the DNB loading reagents

Perform the steps below:

1. Take out the DNB Load Buffer II from the DNBSEQ-T7RS DNB Load Reagent Kit. If it is App-A FCL PE100, also take the App-A Insert primer 1 out from the High-Throughput Pair-End Sequencing Primer Kit (App-A). Thaw reagents at room temperature for approximately 0.5 hours.
2. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use. If crystal precipitation is found in DNB Load Buffer II, vigorously mix the reagent with 1–2 minutes of continuous vortex to re-dissolve the precipitate before use.

5.1.1.3 Preparing the 0.1 M NaOH reagent

Prepare 0.1 M NaOH according to the procedure described in 8.2 *Preparing wash reagents on Page 70*. Each Post Load plate requires at least 4 mL of 0.1 M NaOH.


5.1.1.4 Preparing DNB loading mixture

Perform the steps below:


1. Take a new 0.5 mL microfuge tube and add reagents following the table below.

Table 49 DNB loading mixture for FCL SE35/SE50/SE100/PE100, App-A FCL PE100 and stLFR FCL PE100

Adding order	Component	volume (μL)
1	DNB	270
2	DNB Load Buffer II	90
3	Make DNB Enzyme Mix II (LC)	1

 **NOTE** DNB in the above table refers to the pooled DNB in 4.5.2 *DNB pooling on Page 41*.

2. Gently pipette the DNB loading mix 5 to 8 times by using a wide bore tip.

 **NOTE**

- Do not centrifuge, vortex, vigorously pipette or shake the tube.
- It is important to point out that DNB loading mixture must be prepared fresh and used within 30 minutes.

5.1.2 Post load plate and buffers preparation for FCL PE150, App-A FCL PE150 and App-D FCL PE150 sequencing

5.1.2.1 Thawing the post load plate


Perform the steps below:

1. Take out the Rapid Post Load Plate from the DNB Rapid Load Reagent Kit and thaw in a water bath at room temperature for 2 hours.
2. Once the Rapid Post Load Plate is thoroughly thawed, place it in a 2 °C to 8 °C refrigerator until use. Rapid Post Load Plate can also be thawed in a 2 °C to 8 °C refrigerator overnight.
3. Gently invert the Rapid Post Load Plate to mix 5 times and then centrifuge for 1 minute before use.

5.1.2.2 Preparing the DNB loading reagents

Perform the steps below:

1. Take out the DNB Load Buffer IV from the DNB Rapid Load Reagent Kit.

-  **NOTE**
- If it is App-A FCL PE150, also take out the App-A Insert primer 1 from the High-Throughput Pair-End Sequencing Primer Kit (App-A).
 - If it is App-D FCL PE150, also take the App-D Insert primer 1 out from the High-Throughput Sequencing Primer Kit (App-D).

2. Thaw reagents in a water bath at room temperature for approximately 0.5 hours. After thawing, mix reagents by using a vortex mixer for 5 seconds, centrifuge briefly and place on ice until use.

5.1.2.3 Preparing the 0.1 M NaOH reagent

Prepare 0.1 M NaOH according to the procedure described in 8.2 *Preparing wash reagents on Page 70*. Each Post Load plate requires at least 4 mL of 0.1 M NaOH.


5.1.2.4 Preparing DNB loading mixture

Perform the steps below:


1. Take a new 0.5 mL microfuge tube and add reagents following the table below.

Table 50 DNB loading mixture for FCL PE150, App-A FCL PE150 and App-D FCL PE150

Adding order	Component	volume (μL)
1	DNB	300
2	DNB Load Buffer IV	150

-  **NOTE** DNB in the above table refers to the pooled DNB in 4.5.2 *DNB pooling on Page 41*.

2. Gently pipette the DNB loading mix 5 to 8 times by using a wide bore tip.

-  **NOTE**
- Do not centrifuge, vortex, vigorously pipette or shake the tube.
 - DNB loading mixture must be prepared fresh and used it within 30 minutes.

5.2 Preparing the sequencing flow cell

Perform the steps below:

1. Take out the flow cell from the DNBSEQ-T7RS sequencing flow cell box.
2. Unwrap the outer package.
3. Take out the flow cell from the inner package and inspect if the flow cell is intact.
4. Clean the back of the flow cell using dust remover.

5.3 DNB loading

Perform the steps below:

1. When starting the MGIDL-T7RS, the compartment doors need to be closed.
2. Start the MGIDL-T7RS program. Enter the user name **user** and password **123**, tap **Log in** to enter the main interface.




Figure 1 MGIDL-T7RS main interface

3. Tap on **A** or **B** to continue the operation, see the figure below:




Figure 2 MGIDL-T7RS selection interface

4. Tap on **Loading** and enter the information input interface, see the figure below:




Figure 3 MGIDL-T7RS information input interface

5. Open the loading compartment door.
6. Tap on the text box behind **DNB ID**, enter the DNB information into the text box.

NOTE Use only numbers or letters or a combination of numbers and letters for DNB ID.

7. Place the 0.5 mL micro tube containing DNB loading mix into the DNB tube hole (see *Figure 6 on Page 49*), the screen will prompt that the DNB tube is loaded.

8. Align the post load plate to the RFID scanning area and the ID information will appear in the text box.




Figure 4 The RFID scanning area of post load plate

9. Remove the seal of the post load plate and add 4 mL of 0.1 M NaOH into well No.11.




Figure 5 Adding 4 mL of 0.1 M NaOH into well No.11

- If the sequencing kit is App-A FCL PE100 or App-A FCL PE150, use a pipette to completely remove all the reagent in well No.1, then add 2 ml of App-A Insert primer 1 from High-Throughput Pair-End Sequencing Primer Kit (App-A).
- If the sequencing kit is for App-D FCL PE150, use a pipette to completely remove all the reagent in well No.1, then add 2 ml of App-D Insert primer 1 from High-Throughput Sequencing Primer Kit (App-D).

10. Place the prepared post load plate on the plate tray of MGIDL-T7RS. The screen will prompt that the post load plate is loaded.




Figure 6 Post-loading plate placement diagram

11. Align the flow cell to the RFID scanning area and ID information will appear in the text box.



Figure 7 Scanning the Flow cell ID

12. Hold both sides of the flow cell, upwardly align the locating bulge on the flow cell to the locating groove on the flow cell stage. Gently press down the edges of the flow cell, see the figure below:




Figure 8 Flow cell locating

NOTE Make sure that all the four rubber sealing rings are on the four corners of the flow cell.

13. Press the flow cell attachment button on the flow cell stage to ensure that the flow cell is securely seated and held on the stage. The screen will prompt that the flow cell is loaded.




Figure 9 Flow cell loaded

- NOTE**
- Remove the dust on both sides of the flow cell with a gas dust remover.
 - Do not press or touch the glass cover of the flow cell to avoid flow cell damage or fingerprints and impurities left on the glass surface.
 - Do not move the flow cell after installing the flow cell onto the stage, or it may cause the sealing gaskets to misalign with holes of the fluidics line.
 - If flow cell attachment fails, gently wipe the back of the flow cell and flow cell stage with microfiber clean wiper moistened with 75% ethanol, then clean with a gas dust remover.

14. Close the loading compartment door.

15. Tap the **Start** button and select **Yes** as shown in the Figure 10 on Page 51. Flow cell loading starts as shown in the Figure 11 on Page 51.




Figure 10 MGIDL-T7RS loading confirmation dialog box




Figure 11 MGIDL-T7RS flow cell loading interface

16. The process take around 2 hours. When the screen is shown as in the figure below, the flow cell loading is completed.




Figure 12 MGIDL-T7RS flow cell loading complete status window

17. Press the flow cell attachment button and remove the loaded flow cell from the stage. The flow cell is now ready for sequencing.



NOTE ● If sequencing cannot be performed immediately, put the loaded flow cell in a clean zip bag and store it at 4 °C until use.

- The maximum storage time for loaded flow cell is 48 hours.

18. When the loading is completed, install the washing flow cell onto the flow cell stage and press the flow cell attachment button. Tap **Confirm** button as shown in *Figure 12 on Page 51*.
19. Tap **Post-wash** and select **Yes** to start MGIDL-T7RS wash (see the figures below), which will take around 20 minutes.




Figure 13 MGIDL-T7RS post-wash interface




Figure 14 MGIDL-T7RS post-wash confirmation interface

20. MGIDL-T7RS wash starts, see the figure below:




Figure 15 MGIDL-T7RS wash interface

21. When the screen is shown as the figure below, the wash is completed. Tap **Finish** and another flow cell loading can be performed.




Figure 16 MGIDL-T7RS wash complete status window

Chapter 6 Preparation before sequencing

6.1 Preparing the sequencing reagent cartridge




Figure 17 Sequencing cartridge wells

Perform the steps below:

1. Take out the Sequencing Reagent Cartridge from the High-throughput Sequencing Kit. Thaw in a water bath at room temperature for 4 to 5 hours.
2. After being thoroughly thawed, store it at 2 °C to 8 °C until use (or thaw cartridge in 2 °C to 8 °C fridge 24 hours in advance).
3. Shake the cartridge vigorously in all directions 10 to 20 times to mix well.

NOTE It is a normal phenomenon that dark green crystal appears in well No.1, which is crystallization of raw materials of the reagent in this well. When the cartridge is thawed, mix the reagents in the cartridge well and the crystals will dissolve. Sequencing quality will not be affected.

4. Take out the dNTPs Mix IV (or dNTPs Mix V) and dNTPs Mix II from the sequencing kit and thaw them at room temperature.
5. After thawing, invert the dNTPs Mix IV (or dNTPs Mix V) and dNTPs Mix II 4 to 6 times. Centrifuge briefly and place on ice until use.

6. Process the primer according to different situations:
 - For Dual barcode PE sequencing
 - a. Take out the 1 μ M AD153 Barcode Primer 3 from CPAS Barcode Primer 3 Reagent Kit.
 - b. After thawing at room temperature, vortex for 5 seconds. Centrifuge briefly and place on ice until use.
 - For App-A PE sequencing
 - a. Take out the App-A Insert Primer 2, 1 μ M App-A MDA primer, 1 μ M App-A Barcode Primer 2 from High-Throughput Pair-End Sequencing Primer Kit (App-A).
 - b. Take out the 1 μ M App-A Barcode Primer 3 (just for Dual barcode App-A PE sequencing) from High-Throughput Barcode Primer 3 Reagent Kit (App-A).
 - c. After thawing at room temperature, vortex for 5 seconds. Centrifuge briefly and place on ice until use.
 - For App-D PE sequencing
 - a. Take out the App-D Insert Primer 2, 1 μ M App-D MDA primer, 1 μ M App-D Barcode Primer 2, 1 μ M App-D Barcode Primer 3 (just for Dual barcode App-D PE sequencing) from High-Throughput Sequencing Primer Kit (App-D).
 - b. After thawing at room temperature, vortex for 5 seconds. Centrifuge briefly and place on ice until use.
7. Open the kit cover and wipe any water condensation with lint-free paper. Spray 75% ethanol on the surface of the cartridge seal and clean the seal with lint-free paper.
8. Pierce the seal at the edge of well No.9 and No.10 to make a hole around 2 cm in diameter by using a 1 mL sterile tip.
9. Take out the Sequencing Enzyme Mix from High-throughput Sequencing Kit. Invert the Sequencing Enzyme Mix 4 to 6 times and place on ice until use.

10. Use a pipette with the appropriate volume range and add dNTPs Mix IV (or dNTPs Mix V) and Sequencing Enzyme Mix into well No.9 according to the table below.

Table 51 Sequencing cartridge well No.9 reagent adding

Sequencing kit	dNTPs mix IV volume (mL)	dNTPs mix V volume (mL)	Sequencing enzyme mix volume (mL)
FCL SE35	1.7	/	1.7
FCL SE50	2.0	/	2.0
FCL SE100	3.0	/	3.0
FCL PE100	/	3.9	3.9
FCL PE150	/	5.3	5.3
App-A FCL PE100	/	3.9	3.9
App-A FCL PE150	/	5.3	5.3
App-D FCL PE150	/	5.3	5.3
stLFR FCL PE100	5.4	/	5.4

11. Use a pipette with the appropriate volume range and add dNTPs Mix II and Sequencing Enzyme Mix into well No.10 following the table below:


Table 52 Sequencing cartridge well No.10 reagent adding

Sequencing kit	dNTPs mix II volume (mL)	Sequencing enzyme mix volume (mL)
FCL SE35	4.5	1.5
FCL SE50	5.4	1.8
FCL SE100	8.1	2.7
FCL PE100	10.2	3.4
FCL PE150	13.8	4.6
App-A FCL PE100	10.2	3.4
App-A FCL PE150	13.8	4.6
App-D FCL PE150	13.8	4.6
stLFR FCL PE100	14.7	4.9

12. Seal the loading wells of well No.9 and No.10 with the transparent sealing film.

 **NOTE** Do not cover the center of the well to avoid blocking the sampling needle.

13. Place the cartridge horizontally on the table, and hold both sides of the cartridge with both hands. Shake it clockwise 10 to 20 times, and then counterclockwise 10 to 20 times, make sure reagents are fully mixed. If the sequencing run is SE35/SE50/SE100, the preparation for the reagent cartridge is finished,

 **NOTE** If the sequencing run is SE35/SE50/SE100, the preparation for the reagent cartridge is finished.

14. Perform the following steps according to different situations:

- For PE sequencing
 - a. Pierce the seal of well No.8 by using a 1 mL sterile tip.
 - b. Add 600 μ L MDA Enzyme Mix to the MDA Reagent tube with a 1 mL pipette.
 - c. Invert the tube 4 to 6 times to mix the reagents.
 - d. Add the mixture to well No.8. When adding the mixture, make sure there are no bubbles at the bottom of the tube.





Figure 18 MDA mixing


 **NOTE** When using MDA Enzyme Mix, do not touch the wall of the tube to prevent influencing the enzyme activity.

- For Dual barcode PE sequencing
 - a. Pierce the seal of well No.3 by using a 1 mL sterile tip.
 - b. Add 3.5 mL of 1 μ M AD153 Barcode Primer 3 into well No.3 with a 1 mL pipette.

- For App-A PE sequencing
 - a. Pierce the seals of well No.3, No.4, No.6 and No.13.
 - b. Adding the reagents by using the appropriate pipette according to the table below:

Table 53 Reagents for App-A dual barcode


Primer working solution	Well	Volume (mL)
1 μ M App-A Barcode Primer 2	No.4	3.5
1 μ M App-A MDA Primer	No.6	4.2
1 μ M App-A Insert Primer 2	No.13	4.2
1 μ M App-A Barcode Primer 3	No.3	3.5

 **NOTE** App-A barcode primer 3 is Just for dual barcode App-A PE sequencing.

- For App-D PE sequencing
 - a. Pierce the seals of well No.3, No.4, No.6 and No.13.
 - b. Adding the reagents by using the appropriate pipette according to the table below:

Table 54 Reagents for App-D dual barcode

Primer working solution	Well	Volume (mL)
1 μ M App-D Barcode Primer 2	No.4	3.5
1 μ M App-D MDA Primer	No.6	4.2
1 μ M App-D Insert Primer 2	No.13	4.2
1 μ M App-D Barcode Primer 3	No.3	3.5

 **NOTE** App-D barcode primer 3 is Just for dual barcode App-D PE sequencing.

15. Close the sequencing cartridge cover.

6.2 Preparing the washing cartridge

Perform the steps below:

1. Shake the cartridge clockwise 5 to 10 times, and then counter clockwise 5 to 10 times to ensure the reagents are fully mixed.

2. Spray 75% ethanol on the surface of the cartridge seal and clean the seal with lint-free paper. Pierce either of the well No.2 by using a 1 mL sterile tip.




Figure 19 Washing cartridge

3. Add 45 mL of 0.1 M NaOH into well No.2 through the pierce by using an electronic pipette. Refer to 8.2 *Preparing wash reagents* on Page 70 for the preparation of 0.1 M NaOH.





Figure 20 Washing cartridge added 0.1 M NaOH

6.3 Filling the pure water container

Fill the pure water container with laboratory-grade water according to the table below:

Table 55 Pure water consumption (L)

Sequencing kit	1 flow cell	2 flow cell	3 flow cell	4 flow cell
FCL SE35	1.0	2.0	3.0	4.0
FCL SE50	1.0	2.0	3.0	4.0
FCL SE100	1.5	3.0	4.5	6.0
FCL PE100	3.0	6.0	9.0	12.0
FCL PE150	4.5	9.0	13.5	18.0
App-A FCL PE100	3.0	6.0	9.0	12.0
App-A FCL PE150	4.5	9.0	13.5	18.0
App-D FCL PE150	4.5	9.0	13.5	18.0
stLFR FCL PE100	3.5	7.0	10.5	14.0


-  **NOTE**
- Check whether the water in the pure water container is sufficient. If the pure water is insufficient, the sequencing will fail. Replenish pure water in time, and pay attention to opening the air vent of the pure water container.
 - The pure water will be used in sequencing so it must be kept clean. Renew the pure water in the pure water container on a weekly basis.
 - Before refilling the pure water container, empty the container and spray 75% ethanol on the inner surface of the container lid and the surface of the pure water tube. Wipe and clean the surfaces with new microfiber clean wipers. Rinse the container with fresh pure water 3 times.
 - Refer to *H-020-000157-00 DNBSEQ-T7RS Genetic Sequencer User Manual* for the preparation of the water container.

Chapter 7 Sequencing

7.1 Loading the reagent cartridge

Perform the steps below:

1. Open the reagent compartment door and clean the inner walls with a microfiber clean wiper or lint-free paper moistened with laboratory-grade water. Keep the compartment clean and dry.

 **NOTE** Be careful not to be scratched by the sampling needle above when cleaning the inner walls of the compartment.

2. Place the sequencing cartridge into the sequencing cartridge compartment and place the washing cartridge into the washing cartridge compartment.




Figure 21 Loading the reagent cartridge

3. Close the doors of both sequencing cartridge compartment and washing cartridge compartment, and then close the door of the reagent compartment.

7.2 Entering sequencing interface

Enter the user name **user** and password **123**, tap **Log in** to enter the main interface.




Figure 22 DNBSEQ-T7RS main interface

7.3 Loading the flow cell

Perform the steps below:

1. Select A/B/C/D respectively according to sequencing demand. Tap **Sequencing** and select **New run** (see the figure below).




Figure 23 DNBSEQ-T7RS selection interface

2. Clean the loaded flow cell with a dust remover to ensure no visible dust on the surface and back of the flow cell. Put the flow cell on the flow cell drive, and tap the flow cell drive control button to withdraw the flow cell drive.




Figure 24 The flow cell drive

7.4 Sequencing parameters

Perform the steps below:

1. Align the sequencing cartridge, washing cartridge and flow cell respectively to the RFID scanning area, the ID information will automatically display in the corresponding text box. If the reader fails, information can be entered manually.

5/21/2020 9:52:02 AM

A | B | C | D |

Sequencing cartridge ID 100001635451912260940 ✓

Washing cartridge ID 100001636052003310242 ✓

Flow cell ID E100003005 ✓

Recipe PE150+10 1-128

☒ Split barcode

Advanced settings ⌵

⌂ A: Preparing

◀ Previous ▶ Next

Figure 25 DNBSEQ-T7RS sequencing parameters

2. Tap ▼ behind **Recipe**. Select a appropriate sequencing recipe from the list. If Customize recipe is required, select **Customize** in the drop-down menu to enter the interface as shown in the figure below:

Figure 26 Customize a recipe

- The rules for filling in the **Customize a recipe** interface are as follows:
 - When name a sequencing recipe, use only letters, numbers, “+”, “_” and “-”.
 - Duplicate name check will be performed to ensure that each sequencing recipe name is unique. i.e. a new recipe name must not be the same with an existing recipe.
 - Enter numbers only to the read length of Read1, Read2, Barcode and DualBarcode.
 - Multiple segments of dark reaction cycles can be set in Read1 and Read2. Use “,” to separate each segment and the dark reaction cycles of each segment are presented in the format of “number” and “number-number”.
- Example:
 - The Read1 read length is 100 cycles and the Read2 read length is 100 cycles.
 - Barcode read length is 10 cycles and Dual Barcode read length is 10 cycles.
 - In the 100 cycles of Read1, the 20th to 30th cycles and the 50th to 60th cycles need to perform dark reactions. In the 100 cycles of Read2, the 20th to 30th cycles need to perform dark reactions.
 - Name this recipe as “PE100+10+10+Dark”.

- Fill the **Customize a recipe** interface as shown in the figure below:

Figure 27 Example

NOTE For stLFR FCL PE100, the Read1 read length is 100 cycles and the Read2 read length is 100 cycles, Barcode read length is 42 cycles and DualBarcode read length is 10 cycles.

3. Tap ▼ in the red box of the figure below and select the corresponding barcode sequence. If customized barcode sequence is required, selected the inputted barcode sequence. Select whether split barcode and Dual barcode (stLFR FCL PE100 just select split barcode)

Figure 28 Set the barcode sequence

4. Tap on the **Advanced settings** to enter the interface as shown in the figure below. Users can select whether the primer is **custom primers** and whether to perform **Auto wash**.

The screenshot shows the 'Advanced settings' screen for DNBSEQ-T7RS. At the top, there are tabs labeled A, B, C, and D. Below the tabs, the following settings are visible:

- Sequencing cartridge ID: 100001635451912260940 (with a green checkmark)
- Washing cartridge ID: 100001636052003310242 (with a green checkmark)
- Flow cell ID: E100003005 (with a green checkmark)
- Recipe: PE150+10 (with a dropdown arrow)
- Split barcode: ☒ (checked)
- Advanced settings: ☐ (unchecked)
- Custom primers: ☐ Yes ☒ No
- Auto wash: ☒ Yes ☐ No

At the bottom, there are 'Previous' and 'Next' buttons. A status indicator at the bottom left shows 'A: Preparing'.

Figure 29 DNBSEQ-T7RS advanced settings

NOTE Custom Primers refers to App-A/ App-D and others requiring primer replacement. stLFR sequencing does not belong to Custom primers.

7.5 Reviewing parameters

Tap **Next** to review the parameters and ensure that all information is correct, see the figure below taking PE150 as example:

The screenshot shows the 'Review' screen for DNBSEQ-T7RS. It displays a table of parameters and their values:

Item	Description
User name	user
Sequencing cartridge ID	100001635451912260940
Washing cartridge ID	100001635451912260941
Flow cell ID	E100003005
Recipe	PE150+10
Custom primers	No
Cycle	310
Read1	150
Auto wash	Yes

At the bottom, there are 'Previous' and 'Start' buttons. A status indicator at the bottom left shows 'A: Preparing'.

Figure 30 Reviewing information

7.6 Starting sequencing

Perform the steps below:

1. After confirming that all the information is correct, tap **Start** and select **Yes** (see the figure below):




Figure 31 Confirming sequencing interface

2. When the following screen appears, the sequencing is started.





Figure 32 Sequencing starts interface

3. During sequencing, tap  and the selection of **Auto wash** can be changed as shown in the figure below:




Figure 33 Disable auto wash

4. When the screen appears as shown in the figure below, the sequencing and wash process for this run are completed.




Figure 34 DNBSEQ-T7RS sequencing complete interface


Chapter 8 Device maintenance

8.1 Terminology and definition

Table 56 Wash methods

Wash type	Description
MGIDL-T7RS automatic wash	When the loading is completed, replace the flow cell with a used flow cell and tap Wash . The loader will automatically perform the wash without the need to change the post load plate.
DNBSEQ-T7RS automatic wash	Select Yes for Auto wash, the system will automatically perform a wash after each sequencing run.
MGIDL-T7RS manual wash	<p>Perform a wash manually under the following conditions:</p> <ul style="list-style-type: none">● The device is used for the first time● The device has not been used for 7 days or longer● Impurities are found in the device or flow cell● Replacing the tubing, sampling needles, or other accessories exposed to the reagents.
DNBSEQ-T7RS manual wash	<p>Perform a wash manually under the following conditions:</p> <ul style="list-style-type: none">● The device is used for the first time● The device has not been used for 7 days or longer● Impurities are found in the device or flow cell.● Replacing the tubing, sampling needles, or other accessories exposed to the reagents.

8.2 Preparing wash reagents

 **NOTE** The following Wash Reagents are stored at 4 °C and are valid for 28 days.

- Prepare Wash Reagent I (1 M NaCl+0.05% Tween-20) following the table below:

Table 57 Wash reagent I: 1 M NaCl+0.05% Tween

Reagent	Volume
5 M NaCl solution	200 mL
100% Tween-20	0.5 mL
Laboratory-grade water	799.5 mL

- Prepare Wash Reagent II (0.1 M NaOH) following the table below:

Table 58 Wash reagent II: 0.1 M NaOH

Reagent	Volume
2 M NaOH solution	50 mL
Laboratory-grade water	950 mL

8.3 Washing cartridge

- An empty washing cartridge and washing flow cell for a full wash are provided together with the device.
- The washing plate and washing cartridge must be cleaned before refilled with fresh washing reagents. Replace the washing plate and washing cartridge after three months of continuous use.
- Flow cells from previous runs can be used as washing flow cells. Each flow cell can be used for 3 times.
- Prepare the MGIDL-T7RS washing plate: take a clean and empty post-load plate, add 4 mL of 0.1 M NaOH into well No.11, 4 mL of 1 M Wash Reagent I (1 M NaCl+0.05% Tween-20) into well No.10, 4 mL of laboratory-grade water into well No.9 and 20 mL of laboratory-grade water into well No.12.
- Prepare DNBSEQ-T7RS washing cartridge 1: A clean and empty sequencing cartridge.
- Prepare DNBSEQ-T7RS washing cartridge 2: Take a clean and empty washing cartridge, add 45 mL of Wash Reagent II (0.1 M NaOH) into either of the well No.2, and 45 mL of Washing Reagent I (1 M NaCl+0.05% Tween-20) into either of the well No.3.

8.4 Wash procedures

8.4.1 MGIDL-T7RS manual wash

Perform the steps below:

1. Enter the program.
2. Enter the password **123**, tap **Log in** to enter the main interface.
3. Select the side that needs to be washed.
4. Open the loading compartment door.
5. Place the washing plate filled with wash reagents into the side that needs to be washed. Close the compartment door.
6. Press the flow cell attachment button and wait until the negative pressure is released. Remove the flow cell from the stage. Skip this step if no flow cell is on the stage.
7. Take out the washing flow cell and place it on the flow cell stage. Press the flow cell attachment button and press down the flow cell to ensure the flow cell is securely attached to the stage.
8. Tap **Wash** and select **Yes** to begin the MGIDL-T7RS wash, which will take around 20 minutes.

8.4.2 DNBSEQ-T7RS manual wash

Perform the steps below:

1. Make sure the pure water container is filled with at least 4.5 L of laboratory-grade water before performing the wash.
2. Enter the program. Enter the user name **user** and password **123**, tap **Log in** to enter the main interface
3. Tap **Wash**.
4. Install a used flow cell from a previous run. Press the flow cell drive control button again to withdraw the flow cell drive.
5. Place the clean and empty DNBSEQ-T7RS washing cartridge 1 into the low-temperature compartment on the side that needs to perform a wash, and then close the low-temperature compartment door.
6. Place the DNBSEQ-T7RS washing cartridge 2 filled with wash reagents into the room-temperature compartment on the side that needs to perform the wash, and then close the room-temperature compartment door and the reagent compartment door.

7. Tap **Start** and select **Yes** to begin the DNBSEQ-T7RS manual wash, which will take around 40 minutes.

Chapter 9 Troubleshooting

9.1 Low DNB concentration

When DNB concentration is lower than 8 ng/ μ L, perform the following steps:

1. Check whether the kit has expired.
2. Check whether the library meets the requirements.
3. If DNB concentration still does not meet the requirements after a new sample preparation, please contact the engineer.

9.2 Abnormal negative pressure

When the negative pressure value is shown in red, the negative pressure is abnormal, perform the following steps:

1. Gently wipe the stage surface with a damp lint-free paper or a lint-free cloth and blow the stage with a dust remover and ensure no dust is left.
2. Blow the back of the Flow Cell with a dust remover to ensure no dust is left.
3. If these solutions cannot solve the problem, please contact the engineer.

9.3 Bubbles

9.3.1 Bubbles in MGIDL-T7RS

When bubbles exist in MGIDL-T7RS, perform the following steps:

1. Check the rubber sealing ring to make sure it is in the right position.
2. Check the DNB loading cartridge to make sure enough reagent is in each well
3. Replace the used Flow Cell and inspect the pump.

4. If the problem persists, please contact the field service engineer.

9.3.2 Bubbles in DNBSEQ-T7RS

When bubbles exist in DNBSEQ-T7RS, perform the following steps:

1. Check the water container to make sure water is enough
2. Check the water tube in the water container to make sure it inserts to almost the bottom of the container
3. Check the reagent needles to make sure they can immerse fully into the reagent cartridge. Otherwise, restart the sequencing software.
4. If after restarts, the problem still persists after a restart, please contact the field service engineer.

9.4 Impurities

When impurities exist in the flow cell, perform the following steps:

1. Perform a manual wash on MGIDL-T7RS and DNBSEQ-T7RS.
2. If there is still no improvement after manual wash, follow 8.2 *Preparing wash reagents on Page 70* to reconstitute wash reagents, and perform manual wash again on MGIDL-T7RS and DNBSEQ-T7RS.
3. If there is still no improvement, please contact the field service engineer.

9.5 Pump fails

If liquids cannot be pumped into the flow cell, or large bubbles appear in the flow cell, perform the following steps:

1. Check if the pure water volume is sufficient.
2. When it happens in the MGIDL-T7RS and DNBSEQ-T7RS:
 - 1) remove the flow cell, check if there are impurities on the sealing gasket and remove any dust with a dust remover.
 - 2) Place the flow cell following the instruction and start the pump again.
3. Check if the sampling needles can move properly.

4. If the sampling needles cannot move properly, restart the sequencing software.
5. If the problem persists, please contact the field service engineer.

9.6 Reagent kit storage

- If the kit has been thawed (including dNTPs) but cannot be used within 24 hours, it can be frozen and thawed once only.
- If the kit has been thawed (including dNTPs) but cannot be used immediately, store it at 4 °C and it is strongly recommended to use it within 24 hours. A thawed kit stored at 4 °C less than 7 days still can be used, but it may compromise the sequencing quality. We do not recommend using a thawed kit stored at 4 °C more than 7 days.
- If dNTPs and enzyme have been added into the cartridge, i.e. the cartridge has been prepared but cannot be used immediately, store it at 4 °C and use it within 24 hours.
- If dNTPs and enzyme have been added into the cartridge, i.e. the cartridge has been prepared and the needles have punctured the seal but the cartridge cannot be used in time, the cartridge must be covered with foil or plastic wrap. Store the kit at 4 °C and use it within 24 hours.

Chapter 10 Warnings and precautions


- This product is restricted for research use only, please read the user manual carefully before use.
- Make sure that you are familiar with the SOP and of all the laboratory apparatus to be used.
- Avoid direct skin and eye contact with any samples and reagents. Don't swallow. Please wash with plenty of water immediately and go to the hospital when this happens.
- All the samples and waste materials should be disposed of according to relevant laws and regulations.

Appendix 1 Qubit ssDNA assay kit

 **NOTE** Be careful not to create bubbles.

Perform the steps below:

1. Prepare the Qubit working solution.
 - 1) Diluting the Qubit ssDNA Reagent 1:199 in Qubit ssDNA Buffer. Each sample requires 190 μL of Qubit working solution.
 - 2) Mix by vortexing 2 to 3 seconds.

 **NOTE**

- Each sample DNB quantification requires the preparation of a 200 μL Qubit Working solution.
- Two additional 200 μL Qubit Working solutions are required to build the standard curve.

2. Prepare the required number (N+2) of 0.5 mL tubes for standards and samples (N).
3. Prepare the standard tubes and sample tubes to be tested according to the table below.

/	S1 (μL)	S2 (μL)	D1 (μL)	D2 (μL)	D3 (μL)
working solution	190	190	198	198	198
S1 (0 ng/ μL)	10	/	/	/	/
S2 (20 ng/ μL)	/	10	/	/	/
Sample	/	/	2	2	2
Total	200	200	200	200	200

4. Allow all tubes to incubate at room temperature for 2 minutes.
5. On the Home screen of the Qubit 3.0 Fluorometer, press DNA, then select ssDNA as the assay type. The “Read standards” screen is displayed. Press Read Standards to proceed.
6. Insert the tube containing S1 into the sample chamber, close the lid, then press Read standard.
7. When the reading is complete, put S2 into the sample chamber, close the lid, then press Read standard.
8. When the reading is complete, remove S2. The instrument displays the results on the Read standard screen.

9. Press Run samples. On the assay screen, select the sample volume and units, Insert a sample tube into the sample chamber, close the lid, then press Read tube. When the reading is complete, remove the sample tube. Repeat until all samples have been read.

Appendix 2 Manufacturer

Manufacturer	MGI Tech Co., Ltd. / Wuhan MGI Tech Co., Ltd.
Address	Main Building and Second floor of No.11 Building, Beishan Industrial Zone, Yantian District, Shenzhen, 518083, Guangdong, China
	Building 24, Stage 3.1, BioLake Accelerator, No.388, 2nd Gaoxin Road, East Lake High-Tech Development Zone, 430075, Wuhan, P.R.China
	Building B13, No.818, Gaoxin Avenue, East Lake High-Tech Development Zone, 430075, Wuhan, P.R.China
Service hotline	(+86)4000-966-988
E-mail	MGI-service@genomics.cn
Website	en.mgi-tech.com

---This page is intentionally left blank.---