



**Minimum**  
10 mins

**32**  
Samples

Perform 16/32 samples  
extraction in approximately 10 mins



AUTOMATED NUCLEIC ACID EXTRACTOR |

# MGISP-NE32



## MGISP-NE32

MGISP-NE32 is an automated nucleic acid extractor adopting magnetic rod technology. Together with nucleic acid extraction kit, MGISP-NE32 is able to extract and purify nucleic acid from 16/32 samples in approximately 10 mins.

### Technical Principle

The nucleic acid extraction and purification is performed in four automated steps – lysis, binding, washing and elution. Using magnetic rods and disposable tips, MGISP-NE32 is intended for automated moving and processing of magnetic particles in a 96-well plate.



#### 8-inch Built-in Touch Screen

Create, edit and delete program

#### Reagent and Consumable

96-well pre-packed plate  
Disposable tip

#### Flexible

- Adjustable mixing module with multiple modes & gears, compatible with different magnetic particles and reagents.

#### Efficient

- Moving and processing magnetic particles, extracting and purifying nucleic acid from 16/32 samples in approximately 10 mins.

#### Safe

- UV lamp, 96-well pre-packed plate and disposable tip, avoiding cross-contamination by minimizing contact with sample.

#### Stable

- Specially designed heating block for 96-well plate, ensuring even heating during the experiment.




#### Temperature Control Module

Temperature Control, Lysis/Elution well	From +5°C above ambient temperature to 120°C
Temperature Precision	±1°C

#### Mixing Module

Adjustable mixing module with multiple modes & gears

## Performance



- Test method**

MGISP-NE32 VS Liquid handling method

- Evaluation standard**

RT-PCR

- Test samples**


Simulated samples of throat swab containing novel corona virus (RNA), dilute them in different concentrations

- Test results**

Ct value of fluorescent quantitative PCR is almost equal, that demonstrated the yield of virus RNA is coincident between MGISP-NE32 and liquid handling method. Meanwhile, the reproducibility is great, the CV of Ct value is less than 2.5%.

## ONE-STOP PLATFORM

Provide instrument, consumables and kits for extraction



### MGISP-NE32

- Automated Nucleic Acid Extractor
- Automated extraction step
- Moving and processing magnetic particles
- Nucleic acid extraction from 16/32 samples in approximately 10 mins.



### Nucleic Acid Extraction Kit

- 96-well pre-packed plate
- Good extraction stability and reproducibility
- Compatible with Automated Nucleic Acid Extractor
- Adapted to most downstream applications: PCR, RT-PCR, Sequencing, etc.

# MGISP-NE32

Automated Nucleic Acid Extractor

Perform 16/32 samples extraction in approximately 10 mins

## ORDER INFORMATION

### Equipment

Cat. No.	Product name	Specification
950-000020-00	MGISP-NE32RS Automated Nucleic Acid Extractor	EA, CE RUO
950-000013-00	MGISP-NE32 Automated Nucleic Acid Extractor	EA, CE IVD

### Reagent

Cat. No.	Product name	Specification	Feature
1000023774	MGEasy Nucleic Acid Extraction Kit (OP02-32)	32 preps, RUO	35 mins, Standard Version
1000022606	Nucleic Acid Extraction Kit (OP02-32)	32 preps, CE IVD	35 mins, Standard Version
1000023938	MGEasy Nucleic Acid Extraction Kit (VDR03P-32)	32 preps, RUO	10 mins, Fast Version
1000023937	Nucleic Acid Extraction Kit (VDR03P-32)	32 preps, CE IVD	10 mins, Fast Version
940-000113-00	MGEasy FFPE DNA Extraction Prepacked Kit	32 preps, RUO	High efficiency

## CONTACT US

### MGI Tech Co., Ltd

Building 11, Beishan Industrial Zone, Yantian District, Shenzhen

✉ MGI-service@mgi-tech.com 🌐 en.mgi-tech.com ☎ +86 4000-688-114

#### Copyright Disclaimer

The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

