


The background of the slide is a dark blue gradient with several 3D rendered SARS-CoV-2 virus particles. The particles are spherical with a textured surface and numerous spike-like protrusions. One large particle is on the right side, and several smaller ones are scattered across the left and center.

MGI SARS-CoV-2

Automated Extraction Solutions

MGI has been at the forefront of developing and deploying SARS-CoV-2 testing solutions from day one. We've enabled accurate and reliable virus detection through high throughput automated extraction solutions in more than 80 countries. Improving accuracy and throughput of testing is our top priority.

SARS-CoV-2 Testing Workflow



Swabs & VTM

SAMPLE COLLECTION



● ● Sterile ● ● Safe

Product details

Product Name

Disposable Sampling Kit

Specification

50 PCS/BOX

Product Components

A pcs including a swab and a preservation solution

Tip Material

Nylon flocking

Regulatory Approvals

NMPA; CE
[Check with MGI sales for local compliance]

MGISTP-7000



Convenient Sample Loading

Airtight sample tubes in a standard rack reduces exposure to samples. System can run up to 7000 samples per day.

No Manual Handling

One Epson robot arm, four recap and de-cap modules, four barcode scanners, and two channel-independent pipette heads. Barcode information can be automatically stored and transferred to any LIMS system.

Highly Versatile


Two pipetting channels allow free configuration for volume range and disposable tips. Samples can be transferred into ANSI/SBS 96-well microplates rapidly without manual pipetting errors, supporting a wide variety of automation laboratory equipment.

Unmatched Safety

Re-cap sample tubes after liquid transfer to minimize sample exposure and prevent cross contamination. All samples will be returned to their original location. Provide HEPA-filtration and UV light to ensure a safer operation environment.



High Productivity




*Universal fit. Excellent fit on 5mL and 10mL Transport Tubes, with Screw cap.

Efficient, safe, and reliable



10,000 samples per day:



Product Specifications

Functions	Integrates multiple functions, including decapping, recapping, barcode identification, automated liquid transfer, HEPA-filtration and negative pressure protection	Pipetting Precision	Pipetting Method	Based on the principle of air displacement, with disposable tips
Throughput	192 samples/40 min; 7000 samples/day		Volume Range	10 μ L-1000 μ L
Barcode Scanner	1D & 2D Scanner Sensors		CV	10 μ L: <5%; 100 μ L: <2%
Tube Compatibility	Universal fit. Excellent fit on 5 mL and 10 mL Transport Tubes, with Screw cap		Accuracy	1000 μ L: <1%
Dimensions	1470mm(W) * 960mm(D) * 2100mm(H)		Weight	10 μ L: \pm 10%; 100 μ L: \pm 5%
				1000 μ L: \pm 2%

Ordering Information

Equipment

Cat. No.	Product name	Specification
900-000336-00	MGISTP-7000RS Sample Transfer Processing System	CE RUO
900-000334-00	MGISTP-7000RS Sample Transfer Processing System	CE IVD

Consumables

Cat. No.	Product name	Specification
1000023969	1000ul black conductive tips, sterile, filter	16 boxes/case
1000023970	1000ul transparent tips, sterile, filter	16 boxes/case
014-000101-00	MGISTP-7000 sample rack A	EA

MGIEasy Nucleic Acid Extraction Kit



- Magnetic beads-based automated extraction
- Higher specificity for viral SARS-Cov-2 nucleic acids

EXTRACTION



Workflow:

Lysis/Bind

Sample lysis, nucleic acid binding with magnetic particles.

Wash

Magnetic particles washing, removing contaminants and inhibitors.

Elute

Nucleic Acid separation with magnetic particles and elution in elution buffer.

MGISP-960 High-throughput Extraction

- 192 samples in 80 min
- High recovery and reproducibility
- Upgradable to sequencing library prep capability
- Compatible with throat swabs
- Complete set of supplies



EXTRACTION



RT-PCR PREP

MGISP-100B Medium Throughput and Flexible



- 8, 16, 24, or 32 samples in 40-80 min
- Safe with UV disinfection and positive pressure HEPA system
- Can be upgraded to sequencing library prep capabilities
- All supplies readily available

EXTRACTION



MGISP-NE32 Medium Throughput and Economical

- 1-32 samples in approximately 9 min
- Based on magnetic rod technique, moves magnetic particles instead of liquids; collection efficiency > 98%
- Ready 96-well pre-packed plate and disposable tip
- Lowest per sample extraction cost
- 8-inch built-in touch screen



EXTRACTION




MGISP-NE384

NUCLEIC ACID EXTRACTION



MGISP-NE384 is a high-throughput automated nucleic acid extractor adopting magnetic rod technology. Together with nucleic acid extraction kits, MGISP-NE384 is able to extract and purify nucleic acid from 96/192/288/384 samples. HEPA filtration system and UV lamp ensure a safe working environment for lab technicians.



UV Lamp

More than 100 000 μ W.s/cm²

HEPA Filtration System

99.99% @ 3 μ m

23-inch Built-in Touch Screen, Windows software

User-friendly Interface
Create, edit and delete program

4x96 Mixing & Magnetic Module

Positioning Precision of Robotic Arm: ± 0.1 mm
Z axis Independent Control, 1-4 modules running choice

Temperature Control Module

Temperature Range

From +5 $^{\circ}$ C above ambient temperature to 115 $^{\circ}$ C

Temperature Precision

$\leq \pm 1^{\circ}$ C

Flexible

- Mixing and magnetic module Z axis independent control, 1-4 modules running choice for 96/192/288/384 samples.

Safe

- HEPA filtration system and UV lamp, ISO Class 5, avoiding cross-contamination and ensuring a safe working environment.

Efficient

- Moving and processing magnetic particles, extracting and purifying nucleic acid from 384 samples.

Stable

- Reliable hardware design, 96-well plate reagent and disposable tip, minimizing operational errors.

SARS-CoV-2 Extraction using MGISP-960


- Test method**

Compare automated extraction using RT-PCR with MGEasy Nucleic Acid and Qiagen Virus RNA Extraction Kits with fluorescent PCR to quantify virus RNA.


- Test results**

Lower Ct values using the MGEasy Extraction Kit demonstrated a higher viral RNA yield compared to Qiagen Extraction Kit.


- Comparison of fluorescent PCR peak diagram of two nucleic acid extraction reagent kits :**




Simulated sample of saliva 1



Simulated sample of saliva 2




Simulated sample of saliva 3



Simulated sample of throat swab

In the above diagrams, the blue curve refers to the fluorescent PCR test result of MGEasy Nucleic Acid Extraction Kit + MGISP-960, and the green curve refers to fluorescent PCR test result of competitor's Q Virus RNA Extraction Reagent Kit + manual.

SARS-CoV-2 Extraction using MGISP-100B



- Test method**

Compare manual vs. automation using MGEasy Nucleic Acid Extraction Kit and RT-PCR


- Test samples**

Simulated samples of throat swab containing novel corona virus (RNA), dilute them in different concentrations

- Test results**

Data demonstrates a higher viral RNA yield and better CV using MGISP-100B

SARS-CoV-2 Extraction using MGISP-NE32



- Test method**

Compare manual vs. automation using RT-PCR


- Test samples**

Simulated samples of throat swab containing novel corona virus (RNA), dilute them in different concentrations

- Test results**

Comparable Ct values using both methods demonstrated the similar viral RNA yields between MGISP-NE32 and manual extraction. The CV for Ct value using MGISP-NE32 is less than 2.5%

SARS-CoV-2 Extraction using MGISP-NE384



- Test Method**

MGISP-NE384 VS Liquid handling method

- Evaluation Standard**

RT-PCR


- Test Samples**

Simulated samples of throat swab containing virus(RNA), dilute them in different concentrations

- Test Results**

Ct value of fluorescent quantitative PCR is almost equally, that demonstrated the yield of virus RNA is coincident between MGISP-NE384 and liquid handling method. Meanwhile, Ct value reproducibility is great.

High-throughput Nucleic Acid Testing Laboratory Total Solution



MGISP-NE384

- Automated nucleic acid extraction
- Moving and processing magnetic particles
- 1-4 modules running choice for 96/192/288/384 samples



Nucleic Acid Extraction Kit

- Good extraction stability and reproducibility
- Compatible with Automated Nucleic Acid Extractor
- Adapted to most downstream applications: PCR, RT-PCR, Sequencing, etc.

ORDERING INFORMATION

MGISP-960

Equipment

Cat. No.	Product name	Specification
900-000154-00	High-throughput Automated Sample Preparation System MGISP-960RS, V9	EA, CE RUO
900-000165-00	High-throughput Automated Sample Preparation System MGISP-960, V9	EA, CE IVD

Reagents

Cat. No.	Product name	Specification
1000020261	MGIEasy Magnetic Beads Virus DNA/RNA Extraction	1728 preps, RUO
1000020471	MGIEasy Magnetic Beads Virus DNA/RNA Extraction	96 preps, RUO

Consumables

Cat. No.	Product name	Specification
1000000723	250µL automated filter tips	10 boxes/case
1000004644	1.3mL U-bottom deep-well plate	2 plates/bag
1000012059	Hard-shell thin-wall 96-well skirted PCR plates, white shell/clear well	10 plates/bag

MGISP-100B

Equipment

Cat. No.	Product name	Specification
900-000286-00	MGISP-100BRS Automated Nucleic Acid Extraction and Purification System	EA, CE RUO
900-000285-00	MGISP-100B Automated Nucleic Acid Extraction and Purification System	EA, CE IVD

Reagents

Cat. No.	Product name	Specification
1000020471	MGIEasy Nucleic Acid Extraction Kit	96 preps, CE RUO
1000021042	Nucleic Acid Extraction Kit	96 preps, CE IVD

Consumables

Cat. No.	Product name	Specification
1000000723	250µL automated filter tips	10 boxes/case
1000004644	1.3mL U-bottom deep-well plate	2 plates/bag
1000023283	0.2mL PCR 8-Strips with flat caps	125 strips/box

MGISP-NE32

Equipment

Cat. No.	Product name	Specification
950-000020-00	MGISP-NE32RS Automated Nucleic Acid Extractor	EA, CE RUO
950-000013-00	MGISP-NE32 Automated Nucleic Acid Extractor	EA, CE IVD

Reagents

Cat. No.	Product name	Specification
1000023774	Nucleic Acid Extraction Kit	32 preps, CE RUO
1000022606	Nucleic Acid Extraction Kit	32 preps, CE IVD

MGISP-NE384

Equipment

Cat. No.	Product name	Specification
900-000358-00	MGISP-NE384RS High-throughput Automated Nucleic Acid Extractor	EA, CE RUO
900-000359-00	MGISP-NE384 High-throughput Automated Nucleic Acid Extractor	EA, CE IVD

CONTACT US

MGI Tech Co., Ltd

Building 11, Beishan Industrial Zone, Yantian District, Shenzhen

✉ MGI-service@mgi-tech.com 🌐 en.mgi-tech.com ☎ 4000-966-988

Copyright Disclaimer

The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

