


# ┌ Packages for Population Genomics ┐

Ultra-high Throughput with Flexible and Customizable Automation



## Main Features

- ✓ Unique packages for population genomics
- ✓ Achievable for 50k high-depth WGS per year per set
- ✓ Flexible customization for the full workflow from 50k to 1 million genomes sequenced annually
- ✓ Automation from sample to report

## © DNBSEQ-T10 × 4RS Introduction

DNBSEQ-T10 × 4RS is an ultra high-throughput sequencer based on DNBSEQ™ sequencing technology. It is designed to meet the requirement of population genomics sequencing market.

Read length	PE100	PE150
Effective reads**	360 B	256 B
Average data output	72 Tb	76.8 Tb
Run time***	96 h	108 h
Data quality****	Q30 ≥ 85%	Q30 ≥ 85%

\* The sequencing depth is 30 × ;

\*\* The sequencing data was measured by running 8 sequencing slides;

\*\*\* Run time includes DNB loading and sequencing. FASTQ file generation does not take up sequencer time;

\*\*\*\* ≥Q30 base ratio obtained by sequencing WGS standard library using DNBSEQ-T10 × 4RS

## Sequencing Specification

One DNBSEQ-T10 × 4RS supports the operation of 8 sequencing slides simultaneously, producing up to 20Tb data per day (about 200 whole human genomes sequencing data\*). A single set of DNBSEQ-T10 × 4RS can produce more than 50,000 WGS per year.



The dip-immersion biochemistry, which ensures more uniform biochemical reaction of sequencing slide with larger size, as well as significantly saves the cost of reagents.



Sequencing slide  
DNB loading density is up to 4 million/mm<sup>2</sup>, and one slide can load 70G DNB spots.

## Innovative Sequencing Technologies

Different from ordinary sequencing slide and closed reaction systems used in most platforms, DNBSEQ-T10 × 4RS uses dip-immersion biochemistry and open systems firstly to achieve the best balance among sequencing read length, throughput, data quality and cost.


## Sequencing Data Quality

	PE100		
Sample	NA12878	NA24694	NA24695
Clean data rate	99.97	99.98	99.98
Clean Q30(%)	91.565	91.625	91.62
Mapping rate(%)	99.81	99.81	99.81
Duplicate rate(%)	1.03	0.83	0.88
Mismatch rate(%)	0.3	0.3	0.3
Average sequencing depth(x)	30.66	30.76	30.71
Coverage(%)	99.2	99.89	99.21
Coverage at least 20x(%)	94.39	90.91	94.39
SNP_Precision	0.9994	0.9991	0.9992
SNP_Sensitivity	0.9908	0.9929	0.9927
INDEL_Precision	0.9901	0.9918	0.9917
INDEL_Sensitivity	0.9807	0.9803	0.9802

Human reference genomes NA12878, NA24694 and NA24695 were sequenced using DNBSEQ-T10 × 4RS. The bioinformatic analysis of 30 × WGS data showed that the accuracy and sensitivity of SNPs exceeded 99%, and the accuracy and sensitivity of Indel also exceeded 98%, reaching industry-leading level.

- © Sample — Human Cell Line
- © Library Prep Kit — MGIEasy FS PCR-Free DNA Library Prep Set
- © Sequencing — PE100
- © Data analysis — MegaBOLT


# ©Total Package – Customizable and Automated



\* Calculation is based on human WGS sequenced at 30× depth.

© The sequencing laboratory can be customized with an annual output of  $\geq 50,000$  WGS

3D Layout of Sequencing Laboratory Total Package



Demo Planning for Annual Sequencing of 50,000 and 500,000 WGS\*

	Zone	Biobank	Library center		Sequencing center	Data center
			Library preparation	DNB preparation		
	<b>Main Function</b>	Whole blood separated into plasma and buffy coat for gDNA extraction	From gDNA to library construction	From library to make DNB	DNB loading and sequencing	Write FQ, data analysis, laboratory management system
<b>50,000 WGS per year</b>	<b>Area(m<sup>2</sup>)**</b>	50	50	50	100	20
	<b>Key Equipment</b>	MGISP-960	MGISP-96XL	MGISP-96XL	DNBSEQ-T10 × 4RS	ZLIMS/ZTRON/MegaBOLT Pro
<b>500,000 WGS per year</b>	<b>Area(m<sup>2</sup>)</b>	100	150	150	500	50
	<b>Key Equipment</b>	MGI Lab-S MGISP-960	MGI Lab-L	MGISP-96XL	DNBSEQ-T10 × 4RS	ZLIMS Elite/ZTRON/MegaBOLT Pro

\* Recommended equipment quantity, configuration and site area for annual sequencing capacity of 50,000~500,000 WGS.

\*\* The recommended area is only for core functional zone. The actual area should be determined by the on-site conditions.

MGI Tech Co., Ltd.

Building 11, Beishan Industrial Zone, Yantian District, Shenzhen, CHINA

www.mgi-tech.com

MGI-service@mgi-tech.com

Twitter: [https://twitter.com/MGI\\_BGI](https://twitter.com/MGI_BGI)

LinkedIn: <https://www.linkedin.com/company/mgi-bgi>

+86-4000-688-114

Copyright Disclaimer: The copyright of this brochure is solely owned by MGI Tech Co., Ltd. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating or otherwise) without the prior written permission by MGI Tech Co., Ltd. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

\*Unless otherwise informed, this StandardMPS sequencing reagent is not available in Germany, UK, Sweden, and Switzerland.

Version: August 2023 | MGPA0402202

Information in this brochure is updated to [08/23/2023] and only for your reference. In no event shall the brochure be regarded as warranty or commitment made by MGI Tech Co., Ltd. All rights and obligations shall subject to final executed agreement.

Website

