


# STOmics

## Stereo-CITE Product Note





## Simultaneously characterizes hundreds of proteins along with expressed RNA to facilitate multi-omics study, cell phenotyping and biologically relevant mechanisms.


### 01 Stereo-CITE Proteo-Transcriptomics solution for Protein & RNA co-analysis

Combined with Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE-seq), STOmics Stereo-CITE Proteo-Transcriptomics Set is intended for simultaneous detection of the whole transcriptome and high-plex protein on the same tissue section. Built upon DNA Nanoball (DNB) technology, STOmics Stereo-CITE Proteo-Transcriptomics Set enables a “tissue-to-data” solution through *in situ* capture of the whole transcriptome, at nanoscale resolution and centimeter-sized Field of View (FOV). The Stereo-seq Chip T (poly-T-based chip) is loaded with capture probes containing spatial coordinate information. Through a series of biochemical processes, the probes can capture mRNA molecules and antibody-derived tags (ADTs) *in situ* within the tissue, perform cDNA synthesis, and obtain transcriptome plus multi-protein spatial distribution information of the entire tissue through sequencing and a complementary visualization platform.

The advanced bioinformatics analysis tools now perform a comprehensive analysis of proteome and transcriptome. This feature enables researchers to profile spatial expression maps of both the whole transcriptome and multi-plex proteins on the same section at a large field-of-view. Stereo-CITE facilitates the study of tissue heterogeneity and associated biological significance in complex (disease) tissues, leading to a broader understanding of spatial omics studies.

### 02 Application highlights

- **Spatial *in situ* multi-omics**  
Unbiased spatial profiling of whole transcriptome and hundreds of protein markers on the same tissue section.
- **High Resolution**  
Subcellular resolution for both transcriptome and protein analysis.
- **High-plex protein detection**  
High detection efficiency of 100+ plex proteins with antibody combinations of your choice from in-house validated vendor.
- **Spatial protein profiling via high-throughput sequencing**  
No auto-fluorescence interference or antigen instability caused by multiple rounds of detection.



## Applications of Stereo-CITE in the multi-omics study of human paracancer lymph nodes

### 03 Characterizing spatial gene expression profile on histology images





Figure 1.  
 A. Human paracancer lymph nodes spatial gene expression profile at Bin20 (~10 $\mu$ m)  
 B. Human paracancer lymph nodes spatial protein expression profile at Bin20 (~10 $\mu$ m).  
 C. Pseudo-color image of 4 selected protein expressions.


### 04 Stereo-CITE recapitulated the organizational structure of human lymph nodes

The figure below shows the comparison of the spatial co-clustering results with the H&E-stained image of adjacent sections. It was found that the multi-omics spatial co-clustering results perfectly recapitulated the tissue structure of the lymph nodes. It was observed that the mantle area was composed of mainly naive B cells, the germinal center was composed of mainly cycling B cells, the capsule was composed of mainly endothelial cells, and the subcapsular sinus was composed of mainly mast cells, T cells, B cells, and platelets, which is highly consistent with prior biological and histological knowledge.






Co-clustering proves to be an effective approach for the integration of multiple omics analyses. Specifically, the method enables a clear separation of the red and yellow clusters, which remain indistinguishable through RNA and protein-only clustering approaches. Stereo-CITE exemplifies the advantages of combining diverse omics data to achieve a more comprehensive understanding of biological systems.




**05 Joint analysis of RNA and protein information provides more accurate spatial clustering**

Protein and transcriptional information contribute to the clustering results differently. Statistical diagram demonstrated that proteomic information plays a key role in grouping cancer metastatic related cells, endothelial cells, mast cells and others, where transcriptome information is the main contribution in the cell type identification of B\_plasma and Mast groups.



Proteome and transcriptome information contribute to joint clustering results differently



## Stereo-seq's multi-omics solutions bring greater possibilities and less effort to gaining research insights

### 06 Product Reagents

Stereo-seq Permeabilization Set for Chip-on-a-slide				
Part Number	Product	Specification	Description	
211SP118	Stereo-seq Permeabilization Kit	111KP118	8 RXN	For determining permeabilization parameters to optimize mRNA capture
	Stereo-seq Chip P Slide (1cm*1cm)	210CP118	8 EA	
	STOmics Accessory Kit	1000033700	5 PCS	
Stereo-CITE Proteo-Transcriptomics Set				
Part Number	Product	Specification	Description	
211PT114	Stereo-seq Transcriptomics T Kit	111KT114	4 RXN	For generating a spatially-resolved 3' mRNA library and ADT (protein) library from the same biological tissue sections
	Stereo-seq Chip T Slide (1cm*1cm)	210CT114	4 EA	
	STOmics Accessory Kit	1000033700	5 PCS	
	Stereo-seq Protein Assisted Kit	212KA114	4 RXN	
Stereo-seq Library Preparation Kit				
Part Number	Product	Specification	Description	
111KL114	Stereo-seq Library Preparation Kit	4 RXN	For constructing STOmics Library	
Stereo-seq 16 Barcode Amplification Kit, 16 RXN				
Part Number	Product	Specification	Description	
111KB016	Stereo-seq 16 Barcode Amplification Kit, 16 RXN	16 RXN	16 barcode kit for preparing for and constructing of Stereo-CITE libraries.	
Stereo-seq PCR Adaptor				
Part Number	Product	Specification	Description	
301AUX001	Stereo-seq PCR Adaptor	2 EA	Compatible with PCR thermal cycler as a heating unit	

User manuals and documentations: <https://en.stomics.tech/resources/sop>


### 07 Intuitive Software Pipelines and Suits

**ImageStudio** image processing software, **SAW** (Stereo-seq Analysis Workflow) and **StereoMap** visualization platform are offered free of charge to users and enables them to discover spatial biology knowledge with multiplexed tissue images.

#### Image processing

SAW (Automatic)

ImageStudio (Manual)



#### Reads processing & Quantification


SAW

		MIDCount	ExonCount	CellID
gene1	x1 y1	2	2	0
	x1 y2	3	1	0
gene2	x2 y1	1	0	1
	x1 y1	3	2	0
	x1 y2	2	2	0
	x2 y2	5	3	2

#### Visualization & Additional analysis

StereoMap

Stereopy & Other community developed tools



Learn more: <https://en.stomics.tech/BioinfoTools>





**STOmics** | **MGI**

Reach out to us to learn more:

Version : June 2024

Information in this brochure is only for your reference.  
In no event shall the brochure be regarded as  
warranty or commitment.

**Web: [en.mgi-tech.com](http://en.mgi-tech.com)**

**Email: [MGI-service@mgi-tech.com](mailto:MGI-service@mgi-tech.com)**

