MGI

MGI Non-invasive Prenatal Paternity Testing Package Quick Operation Guide

Part No.: H-020-001063-00Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples) Fill the

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

Overview

The quick operation guide is intended to guide you to use DNBelab-D4RS Pa-SNPs Library Prep Set to perform paternity test.

The material information is as follows:

Туре	Name	Remark		
	DNBelab-D4RS Digital sample preparation system	Cat. No.: 900-000822-00		
	DNBSEQ-G99RS Genetical Sequencer	Cat. No.: 900-000607-00		
Equipment	Pipette (10 μL/20 μL/200 μL)	General laboratory brand		
	Mini vortex mixer	General laboratory brand		
	Mini centrifuge	General laboratory brand		
Software	FGID Forensic DNA Analysis System	Cat. No.: FGI20230002		
	DNBelab-D4RS Pa-SNPs Library Prep Set	Cat. No.: 940-002477-00		
Reagent	DNBSEQ-G99RS High-throughput Sequencing Reagent Set (FCL SE100/PE50)	Cat. No.: 940-001268-00		
	Filtered tip	General laboratory brand		
	Тір	General laboratory brand		
Consumable	0.2 mL PCR tube	General laboratory brand		
	Wide-bore tip	BIOFOUNT, BI-200KX-H or AXYGEN T-205-WB-C		

Preparing reagents

- 1. Take out magnetic beads (A567-Pa-BE) from DNBelab-D4RS Pa-SNPs Genotyping Kit (Box2) (Cat. No.: 940-002108-00) and equilibrate it to room temperature. Mix other reagents by using the mini vortex mixer, centrifuge them briefly and place at normal temperature for later use.
- 2. Take out all the reagents from DNBelab-D4RS Pa-SNPs Genotyping Kit (Box1) (Cat. No.: 940-002106-00), thaw them on ice, and tap the tube eight times (twice for C6-Pa-DNB 2). Centrifuge the reagents briefly and place them on ice for later use.
- 3. Use one of the following groups of BC:

Group	Reagent	Barcode list (1-128)	Group	Reagent	Barcode list (1-128)
	B5-Pa-BC 1	1		B5-Pa-BC 9	9
1	B6-Pa-BC 2	2	7	B6-Pa-BC 10	10
I	B7-Pa-BC 3	3	5	B7-Pa-BC 11	11
	B8-Pa-BC 4	4		B8-Pa-BC 12	12
	B5-Pa-BC 5	5		B5-Pa-BC 13	13
2	B6-Pa-BC 6	6		B6-Pa-BC 14	14
Z	B7-Pa-BC 7	7	4	B7-Pa-BC 15	15
	B8-Pa-BC 8	8		B8-Pa-BC 16	16

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples Fill the samp

Fill the sample preparation cartridge Ma

Make DNBs Perform sequencing

Create and perform an analysis task

- It is recommended to physically separate working areas in the laboratory. It is recommended to prepare samples and prepare sample preparation cartridge by using the filtered tips in the pre-PCR area and prepare DNBs by using the DNBelab-D4RS in the post-PCR area. By doing this, the contamination can be prevented and the accuracy of the results can be guaranteed.
- Do not place C6-Pa-DNB 2 at room temperature. Do not touch the tube wall for a prolonged time.
- Protect C7-Pa-QU away from light.

Preparing samples

🕜 Tips

- It is recommended that the concentration of gDNA samples should be greater than or equal to 0.2 ng/µL and the concentration of cfDNA samples should be greater than or equal to 0.4 ng/µL.
- It is recommended to extract a fresh cfDNA before the library preparation.
- It is recommended to store the cfDNA at -20 °C or -80 °C and use it within a week. If the cfDNA is stored more than a week, it is recommended to reextract the cfDNA by using the plasma stored at -80 °C or the fresh plasma.
- 1. Separately aspirate and dispense 4 ng cfDNA of fetus and 2 ng gDNA of father and mother into three 0.2 mL PCR tubes, and add A123-Pa-EB until the total volume reaches 10 µL. Use a new tip when aspirating reagents every time.
- 2. Sample mixture: add 20 μ L A567-Pa-BE into the 10 μ L sample in the step 1. Mix the mixture three times by using the mini vortex mixer (3 seconds each), centrifuge the mixture briefly to collect the mixture at the bottom of the tube, and place it at room temperature for later use.

Filling the sample preparation cartridge

1. Place the instruction card onto the sample preparation cartridge.



2. Insert the funnel into well SF/F3.



3. Uncap the sealing fluid tube, and empty the sealing fluid into the opening of the funnel at one time.



MGI Non-invasive Prenatal Paternity Testing Package Quick Operation Guide

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples Fill the

Fill the sample preparation cartridge

Make DNBs P

Perform sequencing Create and perform

Create and perform an analysis task

- 4. Remove the funnel when the sealing fluid has spread throughout the cartridge. Ensure that no bubbles exist in the cartridge before going to the next step.
- 5. Add the reagents into the cartridge according to the table below. Use new tips to add the reagents every time.

- The loading order for the reagent is well A, well B, and well C.
- Insert the tip vertically into the well, tilt it about 30° towards the intermediate reaction zone, then dispense the liquid.



Avoid generating bubbles.

Name	Cap color	Volume	Well	Storage temperature
A123-Pa-EB		30 µL	A1, A2, A3	2 ℃ to 8 ℃
A4-Pa-DNB 1		30 µL	A4	-25 ℃ to -15 ℃
A567-Pa-BE		30 µL	A5, A6, A7	2 ℃ to 8 ℃
B1-Pa-DNBB	0	15 µL	B1	-25 ℃ to -15 ℃
B2-Pa-PCR 1		15 µL	B2	-25 ℃ to -15 ℃
B4-Pa-ST	0	15 µL	B4	-25 ℃ to -15 ℃

Name	Cap color	Volume	Well	Storage temperature
B5-Pa-BC 1				
B5-Pa-BC 5		<i>.</i> .	55	25.00.015.00
B5-Pa-BC 9		6 µL	В5	-25 °C to -15 °C
B5-Pa-BC 13	-			
B6-Pa-BC 2				
B6-Pa-BC 6		<i>.</i> .	5.0	25.02.1.15.02
B6-Pa-BC 10		6 µL	B6	-25 °C to -15 °C
B6-Pa-BC 14	-			
B7-Pa-BC 3				
B7-Pa-BC 7				
B7-Pa-BC 11		6 µL	В7	-25 °C to -15 °C
B7-Pa-BC 15	-			
B8-Pa-BC 4				
B8-Pa-BC 8		<u> </u>	50	25.00.00
B8-Pa-BC 12		6 µL	88	-25 °C to -15 °C
B8-Pa-BC 16	-			
C6-Pa-DNB 2	0	6 µL	C6	-25 °C to -15 °C
C7-Pa-QU	0	6 µL	C7	2 ℃ to 8 ℃

6. Separately add 30 μL of different sample mixtures into F1, F2, F3, and F4 wells.

For Research Use Only. Not for use in diagnostic procedures.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents P

Prepare samples Fill the sample preparation cartridge

rtridge Make DNBs

e DNBs Perform sequencing

Create and perform an analysis task



Tips

The order of samples mixtures in wells F1-F4 is father, mother, and fetus.

 If there are more than one father's sample, add the sample mixtures into F1 to F4 according to the figure below.



Mother 1 Father 2

 If there is only one father's sample, add the sample mixtures into F1 to F4 according to the figure below.



Fetus Mother 1

- 7. Ensure that all the reagents and samples have been added into preparation cartridge and covered by the sealing fluid.
- 8. Remove the instruction card.

Making DNBs

1. Turn the power switch on the back of the DNBelab-D4RS Digital sample preparation system to the **ON** position. The device starts self-check.

After completion, the login interface is displayed.

- 2. Input the username (*admin*) and default password (*123456*) in the login interface and tap **Login** to go to the main interface.
- 3. Tap (B) to start the workflow.
- 4. Scan the QR code on the packaging of the Prep set by using the scanner to input the information. After completion, tap () to go to the next step.



5. Scan the QR code on the package of the sample preparation cartridge by using the scanner to input the information. After completion, tap (>) to go to the next step.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.





6. Tap the **Sequencing platform** list, and select **DNBSEQ-G99**. After completion, tap (>) to go to the next step.



7. Fill in with sample name, barcode, and sample attribute. After completion, tap (>) to go to the next step, and tap (>) again to go to the next step. The compartment door and cartridge rack open automatically.



8. Place the prepared sample preparation cartridge onto the cartridge rack according to the on-screen instruction. After completion, tap (>>) to go to the next step.

() ()) () ()	0224/09/11 10:43:20
Please install the cartridg	ge.

- 9. Ensure that the cartridge is installed. Tap **Yes**. The Cartridge loading interface is displayed, and the cartridge rack withdraws automatically.
- 10. When you are prompted that the cartridge loading is completed, manually close the compartment door. After completion, tap () to go to the next step.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.





11. Check the information. Ensure that the information is correct, and tap $\widehat{(}^{\!\!\!\text{Run}\!\!\!\!}$.



During running, the interface displays the remaining time and completion time.



12. After library preparation is completed, tap \bigcirc . A library preparation report is displayed.



The report displays the DNB concentration and recommended input volume. The DNB concentration should be greater than or equal to 8 ng/ μ L.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples Fill the s

Fill the sample preparation cartridge

Make DNBs Perform sequencing

g Create and perform an analysis task

			Summar	y report	
СН	Sample ID	Barcode	Attrs	Result(ng/µL)	Input volume(µL)
CH1	Demo-F1	1	F1	16.16	1.86
CH2	Demo-F2	2	F2	17.04	1.76
СНЗ	Demo-M	3	M1	17.00	1.18
CH4	Demo-C	4	Ftl	9.59	8.34

- 13. Tap (>) to go to the next step. The compartment door and cartridge rack open automatically.
- 14. Remove the sample preparation cartridge, and place it on a level and clean surface for later use.



15. Aspirate DNBs:

- (1) Adjust the pipette range to 50 μ L.
- ⁽²⁾ The wells D1 to D4 correspond to the well F1 to F4. Insert four wide-bore pipette tips vertically into the wells D1-D4 of the sample preparation cartridge, separately aspirate the DNBs from the wells, and transfer them into four new PCR tubes respectively. Ensure that the tips remain tightly closed when inserted into the wells.

- 3 Remove the upper oil by using new 200 µL or 10 µL sharp pipette tips until the transparent water phase (the DNB product) in the middle exposes.
- ④ Take out a piece of DNB separation paper and place it on the laboratory bench.
- (5) Adjust the pipette range to 10 μ L, insert a new 10 μ L sharp pipette tip vertically into the DNBs, and aspirate 10 μ L DNBs. During aspirating, lower the tip when the DNB level declines. Avoid aspirating the oil phase. Use a new tip when aspirating the DNBs every time.
- ⁽⁶⁾ Transfer the DNBs onto the DNB separation paper. If there is no wet spot on the paper, proceed to step (7). If there is wet spot on the paper, transfer the DNBs to other clean position of the paper by using a new wide-bore tip until no obvious wet spot is left on the paper.
- 7 Aspirate the DNBs without oil and transfer them into four new PCR tubes (do not transfer them into one tube) by using four new 10 µL sharp pipette tip.
- $(\ensuremath{\$})$ According to the mix volume on the interface of DNBelab-D4RS, aspirate the corresponding volume of DNBs, transfer them into a new 1.5 mL centrifuge tube, and gently mix the DNBs twice by using a new wide-bore pipette tip. The mixture volume should be greater than or equal to 7.5 µL. If the volume is less than 7.5 µL, increase the volume of DNBs in equal proportion. The DNB mixture and the remaining DNBs can be stored in a 4 °C refrigerator.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples) Fill the

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

Tips

- If you fail to aspirate the DNBs, dispense the oil into F well, change a new wide-bore tip and aspirate again.
- The DNBs can be stored in 4 °C refrigerator for later use (use it within 48 hours).
- 16. Tap \bigcirc > Yes. The cartridge rack withdraws automatically.
- 17. After completion, manually close the compartment door. Tap (>) to return to the main interface.



18. If you do not continue to use, tap (III) > Shut down, and turn the power switch on the back of the device to the OFF position.

Performing sequencing

Preparing flow cell

Place the flow cell at room temperature for 30 minutes to 24 hours before loading the DNBs.

Tips

Do not open the package at this moment.

Preparing the sequencing cartridge

- Take out the sequencing reagent cartridge and place it in water bath at room temperature for 2 hours, or pre-thaw it in the refrigerator at 2 °C to 8 °C one day ahead. Then, place the sequencing cartridge at 2 °C to 8 °C until use.
- 2. Invert the cartridge for 5 times. Tear off the package and wipe the condensate water on the package and wells with a dust-free cloth.
- 3. Use the Puncher to pierce the M1, M2, M3, M4 wells of the cartridge with the pre-mixed reagents.





• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples) Fill the

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

4. Hold the A and B sides of the cartridge with both hands by following the marks A and B on the cartridge, and shake it vertically and horizontally for 10 to 20 times to mix the reagents thoroughly.



5. Pierce the MDA well on the cartridge with a clean 1 mL pipette tip.

The prepared cartridge can be placed at 4 $\,\,^{\rm o}{\rm C}$ until use and should be used within 24 hours.

Sequencing

- 1. Tap _____ in the upper-right corner of the main interface of the control software of the sequencer. Enter the user name (*user*) and password (*123*). Tap Login. The main interface returns to view.
- 2. Tap **Sequence** in the A or B operation area in idle status. If you need to perform sequencing on both the stages A and B, tap **Sequence A&B**.
- 3. After tapping **Sequence**, perform one of the following operations:
 - If the waste container compartment door opens automatically, place an empty waste container into the compartment according

to the on-screen prompts, and gently press the compartment door to close it. The system automatically starts the check before sequencing.

- If the waste container compartment door does not open, the system automatically starts to check.
- 4. After the check is completed, tap **Next** to enter the sequencing parameters.

Select **Sequencing only** and **No** in the **BBS** box by default. Input the DNB ID in the **DNB ID** box.

- 5. Select **SE50+10** in the **Recipe** box. Then select the barcode list **1-128** in the box next to the **Recipe** box. If there is no such recipe or the SoftwareVersion of the device has been updated to 1.7.0.674 or above from a previous version, delete the existing recipe and customize a recipe. For detailed operation, refer to the Q18 of *FAQ on Page 16*.
- 6. Select **Yes** for both **Split Barcode** and **Auto Wash** in the advanced settings.
- 7. Tap Next. The auto-sliding screen moves up automatically.
- 8. Push the sequencing cartridge into the reagent compartment. The system will automatically scan the cartridge ID and display it in the **Sequencing Cartridge ID** box.

Tips

If the automatic scanning fails, you can manually enter the ID.

- 9. Tap **Prime** and then **Yes** in the pop-up window. The system automatically starts priming.
- 10. After priming, the auto-sliding screen moves up automatically.

MGI

MGI Non-invasive Prenatal Paternity Testing Package Quick Operation Guide

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples Fill the s

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

Loading DNBs

- 1. Take out DNB Load Buffer II and TE Buffer place it on ice for 30 minutes until it thaws. Mix it thoroughly by using the vortex mixer, centrifuge it for 5 seconds, and place it on ice until use.
- 2. Mix 7 µL DNB load buffer II and 1 µL Make DNB Enzyme Mix II (LC), pipette gently to mix well to prepare as the DNB loading mix solution.
- 3. Add the DNB loading mix solution into the mixed DNBs prepared in *Making DNBs on Page 4* (the volume should be 7.5 to 21 μ L) according to the table below. The total volume of the mixture should be greater than 10 μ L.

	Volume (µL)					
Component	7.5 ≤ V < 21	V≥ 21				
DNB loading mix solution	8×(V/21)	8				
Mixed DNBs	V	21				
Total volume	29×(V/21)	29				

- 4. Gently mix the mixture for 5 to 8 times by using the wide-bore nonfiltered pipette tip. Place the mixture at 4 °C until use.
- 5. Take out the flow cell from the inner package and inspect it to ensure that the flow cell is intact.
- 6. Aspirate 10 μ L of the DNB loading mixture by using the 200 μ L, nonfiltered pipette tip and insert the tip into the inlet.

🕜 Tips

- The DNB loading mixture should be prepared freshly.
- Do not pipette vigorously, or shake the tube.

- 7. Fix the tip with one hand, press the tip ejector on the pipette to unload the tip with the other hand, and observe the liquid level in the tip:
 - If the liquid level drops automatically, the DNB loading mixture will automatically flow into the flow cell.
 - If the liquid level does not drop, perform the following steps:
 - a. Leave the tip with DNB loading mixture in the inlet.
 - b. Adjust the aspirate volume to 2 μL and take a new 200 μL non-filtered pipette tip.
 - c. Hold the new empty pipette tip with one hand and gently insert it into the outlet while pressing the button down with the other hand.
 - d. Gently release the button and remove the tip in the outlet after the liquid level of the tip at the inlet drops.
- 8. When the liquid level in the pipette tip stops dropping, remove the pipette tip in the inlet. Turn the flow cell upside down.

Loading flow cell and reviewing sequencing information

- 1. Insert the prepared flow cell into the flow cell compartment in the direction indicated by the arrow on the flow cell. The system will scan the flow cell ID automatically. If automatic scanning fails, you can manually input the ID.
- 2. Tap **Next**. The auto-sliding screen moves down automatically, and the review interface is displayed.

MGI

MGI Non-invasive Prenatal Paternity Testing Package Quick Operation Guide

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples) Fill the s

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

3. In the review interface, review all items in the review interface. After checking that all information is correct, tap **Sequence**, and then select **Yes** to start sequencing.

The Sequencing interface of the control software displays real-time sequencing progress and step.

4. After the sequencing is completed, tap **Finish** to end the sequencing workflow. The auto-sliding screen automatically moves up and the waste container compartment door is automatically opened. Check the summaryReport, and ensure that the TotalReads(M) is greater than or equal to 80 M, the Q30(%) is greater than or equal to 90%, the SplitRate is greater than or equal to 90%.

Device maintenance

If **Auto Wash** is selected when setting the sequencing parameters, the device starts an automatic wash after the sequencing is completed. For details, refer to *DNBSEQ-G99RS* High-throughput Sequencing Set User Manual.

Creating and performing an analysis task

Selecting analysis product

1. Launch the Chrome browser, set the language to English, type the following IP address in the address bar, and press **Enter** to go to the login interface of the FGID:

192.168.1.3:8080

2. Select the default username to go to the main interface.

3. Select **FGID** to go to the system.

Inputting the sample information

Click **Sample** and select **Sample List** to go to the sample list window.

	Sample List										halos_ma	nana
erensic DNA Analysis System						Name	Received Date					
FGID 🗸	Submission Date	e 🔹 Submission Date	Search Refres									
	Add	import Export	Delete Template Down	nload								
Sample	No.	Sample Id $\$	Associated Project	Submitter	Client Institution	Received Date 🗘	Submission Date	Laboratory	Sample Name	Opera	ion	
Sample List	0 1	240429-LY20240								Edit	Delete	ð
E Project	2	240429-LY20240								Edit	Delete	ð
Task	3	240429-LY20240								Edit	Delete	8
Analysis	. 4	240429-XYZ-071								Edit	Delete	
, analysis	5	240429-YJJ-0419								Edit	Delete	,
Report	6	240429-9947A-1								Edit	Delete	,
Applications v	7	240429-9947A-1								Edit	Delete	
	8	240429-9947A-1								Edit	Delete	,
	9	240429-9948-15								Edit	Delete	
	10	240429-9948-15								Edit	Delete	,
	11	240429-9948-15								Edit	Delete	
	12	240429-9948-10								Edit	Delete	,
	<											>

Select one of the following methods to import sample information:

- Method 1: Manually inputting
 - 1 Click Add and input sample information in the Add Sample window.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents	Prepare samples	Fill the sample preparation cartridge	Make DNBs	Perform sequencing	Create and perform an analysis task
------------------	-----------------	---------------------------------------	-----------	--------------------	-------------------------------------

* Sample Id:	Demo_F	Inspection Person:	
Delegate Unit:		Test Address:	
Accept Date:		Test Date:	
Patient Name:		Gender:	Please select V
Id Card:		Birthday:	
Age:		Project Name:	

- 2 Click **Submit** or **Submit** and **Continue** to continue inputting sample.
- Method 2: Inputting information by Excel
 - ① Click **Template Download** to download the sample template.
 - 2 Input the information into the template and click **Save** after completion.

Tips

A field with a red asterisk (*) is required, and other fields are optional $(Age \ge 1, the Id Card number should include less than or equal to 18 characters.).$

3 Click Import, click Upload in the Import window, select the *.xlsx* document that has been filled in samples information, and click Confirm.

* Order Form Template	FGID Sample V
Project Name	
* Import File	Upload

Creating an analysis task

Create an analysis task according to the following steps:

remain DNA Academia Sentem	Ratab	No			hale Ctatus	ER Create T	ima Croate Tima	Controls	Defeash		
FCID	Daton	I NO.	Plow C		ask Status V		me • Create time	Search	Reifesti		
roib v	A	dd Task	Import Task	Template Down	load						
	No.	В	atch No. 🗘	Flow Cell ID $\ \hat{\bigcirc}$	Sequencer ID	Platform	Progress	Start Time 🗘	End Time 🗘	Task Status Task	Errs Operation
Sample ^	1	R	240510000531	FT100025863	R1100700220011	DNBSEQ-G99				Started	Redo
Sample List	2	R	240510000530	FT100025863	R1100700220011	DNBSEQ-G99				Started	Redo
Project	3	R	240510000529	FT100048021	R12345678	DNBSEQ-G99	-			Started	Redo
= Troject	4	R	240510000528	FT100048021	R12345678	DNBSEQ-G99	-			Started	Redo
Task	5	R	240510000527	FT100048021	R12345678	DNBSEQ-G99	-			Started	Redo
Analysis	6	R	240509000526	FT100025863	R1100700220011	DNBSEQ-G99				Started	Redo
Report	7	R	240509000525	FT100025863	R1100700220011	DNBSEQ-G99				Started	Redo
Hoport	8	R	240509000524	FT100048601	R1100700220011	DNBSEQ-G99	-			Started	Redo
Applications ~	9	R	240509000523	FT100048544	R1100700220011	DNBSEQ-G99	-			Started	Redo
	10	R	240509000522	FT100025863	R1100700220011	DNBSEQ-G99				Await	Stop Delete Redo
	11	R	240509000521	FT100025863	R1100700220011	DNBSEQ-G99				Started	Redo
	12	R	240508000520	FT100047970	R1100700220011	DNBSEQ-G99				Started	Redo
	13	R	240508000519	FT100049584	R1100700220011	DNBSEQ-G99	-			Started	Redo
	14	R	240508000518	FT100048601	R1100700220011	DNBSEQ-G99				Started	Redo
	15	R	240507000517	FT100051477	R12345678	DNBSEQ-G99				Started	Redo
	16	R	240507000516	FT100050760	R12345678	DNBSEQ-G99				Started	Redo

- 1. Click **Task > Template Download** to download the task template.
- 2. Input **Cutoff** in the H1 cell of the template.

NG

MGI Non-invasive Prenatal Paternity Testing Package Quick Operation Guide

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples Fill the s

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

- 3. Input the following information into the **FGID_task_template**:
 - ① Input L01 in the Lane ID column.
 - 2 Input DNB ID in the **DNB ID** column and the DNB ID in an analysis task should be consistent.
 - ③ Select No for Is Dual Barcode column.
 - ④ Input barcode ID in the barcode ID1 column and leave barcode ID2 column blank.
 - 5 Select MGEasy Pa-SNPs Genotyping Kit for the Kit Name column.
 - 6 Input data amount for the sample (unit: M) in the Cutoff column. Input 5 for the parent sample which indicates that cut 5 M reads for analysis. Input 25 for fetal sample and input None if cutoff is not required.

Tips

- The sample ID should be consistent with the sample ID imported into the FGID system (including the upper- and lower-case letters), and one barcode can only correspond to one sample ID within a flow cell.
- It is recommended to submit one analysis task for one sample ID to avoid sample confusion during the analysis.

	A	В	C	D	E	F	G	Н
1	*Sample ID	*Lane ID	*DNB ID	*Is Double Barcode	*barcode ID1	barcode ID2 (If [Is Double Barcode] is yes, barcode ID2 is required)	•Kit Name	Cutoff
2	Demo_F	L01	20240430	No	1		MGIEasy Pa-SNPs Genotyping Kit	5
3	Demo_M	L01	20240430	No	2		MGIEasy Pa-SNPs Genotyping Kit	5
4	Demo_C	L01	20240430	No	3		MGIEasy Pa-SNPs Genotyping Kit	25
5								
6								
7								
8								

- 4. Click Import Task in the Task interface.
- 5. Select Sequencer ID, input Flow Cell ID, and click Upload.
- 6. Select the task list, click **Preview** to preview the imported task details. Ensure that the information is correct and click **Confirm**.

For Research Use Only. Not for use in diagnostic procedures.

When the task status is in **Started**, the task is submitted successfully.

Starting analysis

After the sequencing is completed, the analysis system automatically starts the basic analysis task.

Viewing basic data analysis results

In the FGID system main interface, click **Analysis** to view the progress of basic data analysis.

D -														
-	Batch No.	Plow Cell ID		Start Time +	Start Time	Search	Retresh							
	No.	Batch No. 0		Flow Cell ID	Total Of Sample	SubBatch Count	Tinished Ratio			Start Time 0	End Time (Chip Q	C Files
	1	R240827000815		· PT100020749	4	1			0/1 (0%)	2024-05-27 30:17:58			Down	load
	No.	Sub Datch No. Kit	Name		Current Step	Sample C	reat Status	Error Mesage				Start Time	End Time	Opr
173.	1	R24052700051501 38	GIEany Pa-SNIPs Gen	otyping Kit	cut	4	Antiyong					2024-08-27 10:17:58		544
	2	R240823000811		~ PT100011848	1	1			1/1 (100%)	2024-08-23 14:45-07	2024-08-23 15:20:54		Down	load
_	Ne.	Sab Bands No. Kit	Name		Carreat Step	Sample C	read Status	Error Mesage				Start Time	End Time	Ope
		#340#32000#1101 \k	(Thurs Do 1930). Com	staries Vit			Patient					2014 /8 21 14 45 27	MAK 48 13 15 16 16	
				or the set										
ns v	3	R240823000810		> PT100038317	3	1			1/1 (100%)	2024-08-23 13:49:55	2024-08-23 13:57:23		Down	load
	4	R240823000809		> PT100038317	4	1	-		1/1 (100%)	2024-08-23 13:06:50	2024-08-23 13:32:53		Down	load
	5	R240823000808		> PT100011545	7	1	-		1/1 (0099)	2024-08-23 13:05:00	2024-05-23 13:30:24		Down	load
	6	R240823000807		> PT100020781	3	1	-		1/1 (100%)	2024-08-23 09:58:18	2024-08-23 30:26:09		Down	load
	7	R240823000806		> FT100025970	1	1			1/1 (100%)	2024-08-23 09:58:15	2034-08-23 32:06:11		Down	load
	8	R240823000805		> V350067820	3	1			1/1 (100%)	2024-08-23 09:58:26	2034-08-23 30:52:08		Down	load
	9	R240823000804		> PT100001546	3	1			1/1 (1099)	2034-08-23 09:58:09	2024-08-23 32:34:38		Down	load
	10	R240825000803		> PT100021129	1	1			1/1 (100%)	2024-08-23 09:57:56	2024-08-23 32:34:39		Down	load
	11	R240823000802		> PT100007998	3	1			1/1 (100%)	2024-08-23 09:02:54	2024-05-23 09:47:08		Down	load
	12	R240823000801		> PT100025864	1	7			00010 T/7	2024-08-23 09:01:06	2024-08-23 10:06:09		Down	load
	13	R240822000800		> PT100021895	3	1	-		1/1 (100%)	2034-08-22 17:47:34	2024-08-22 18:13:36		Down	load
	14	R240822000799		> PT100045514	4	1			1/1 (100%)	2034-08-22 17:39:39	2024-08-22 18:09:36		Down	load
	15	R240822000798		> PT100021456	4	1	-		1/1 (000%)	2024-08-22 17:38:04	2034-05-22 18:06:36		Down	load
	16	R240822000796		> PT100021895	3	1			1/1 (100%)	2024-08-22 16:43:16	2034-08-22 17:16:06		Down	load
	17	R240822000795		> PT100001546	3	1			1/1 (100%)	2034-08-22 16:39:30	2034-05-22 17:00:36		Down	load
	18	R240822000794		> FT100020781	3	1			1/1 (100%)	2024-08-22 35:05:01	2034-08-22 16:20:35		Down	load
	19	R240822000795		> PT100020781	3	1			1/1 (100%)	2024-08-22 14:25:21	2024-05-22 14:41:05		Down	load
	20	R240822000792		> FT100020789	8	1			1/1 (200%)	2024-08-22 10:56-28	2024-08-22 11:30:05		Down	load
	21	R240822000791		> PT100020721	3	1			0/1 (210)	2024-08-22 10:50:58	2034-06-22 32-51-35			
				-									200	

- 1. Click > in front of **Flow Cell ID** to check the progress or error reason of this analysis task.
 - Status shows the current analysis task status. If it displays Finished, it indicates that basic data analysis is completed and you can proceed to the next step of paternity analysis.

Part No.: H-020-001063-00
Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples) Fill the sam

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

- Current step indicates the analysis step of the current task. Click the status below Current step to view the analysis flowchart and current progress.
- 2. After the task is completed, click **Batch No.** to view the basic data analysis results for individual samples.

Click **Export** to export the result files.

Submitting paternity analysis tasks

1. Click **FGID** > **Applications** > **Forensic analysis** > **View More** > **NIPPT** to go to the paternity analysis application.





2. Select Gene Locus Type > Pa-SNPs Genotyping Kit Analysis parameters.

GXD 📑	Tweese Analysis 1 Treeses Enalysis	adı
FCID V	O Forware Liketifación belete Tangie	
Sample ~	Naglodan 1077 v	
lask Analyzis	 *Analysis of PS-SPF Complete Leaders *Data Complete 10 	
pplications		
Forensic Analysis		
Identification Result		
	Catool Next	

8G

MGI Non-invasive Prenatal Paternity Testing Package **Quick Operation Guide**

• Part No.: H-020-001063-00 • Version: 1.0

• Release date: Oct. 2024 ©MGI All rights reserved.

Prepare reagents

Prepare samples

Fill the sample preparation cartridge

Make DNBs Perform sequencing

Create and perform an analysis task

3. Click **Next** to go to the paternity analysis sample selection interface. Select samples which need to perform paternity analysis:

FGID .	Forensic Analysis / Forensic A	nalysis												admin-
FGID V			Form	(1)- neie Identification							Se	(2) lect Sample		
(h. name	Unselected							s	ielected					
G Task	Demo	Project Name		Sample Name		learch C I	Tafwah	A	Application : NEPPT Gene La	cua Type : Pa-SNPa (Genetyping Kit Analysis			
CO Machine	No. Sample Id	Associated Project	Batch No.	Marker Type	Sample Name	Client Institution	Received Date		No. Identity Requirement	Sample ID	Sample Name	Sample Type	Marker Type	Operation
	🖬 1 Demo_M		9240823000810	NPCR						Denzo_M		Alleged Mother	U NPCR	
E Report	2 Demo_C		R240823000810	NPCR.					1 Trio Paternity 🧠	Demo_C		Alleged Child	· NPCR	Delete
S. Applications	3 Demo_F		9240823000810	NPCR.						Denao_F		Alleged Father	- NPCR	
								>						
	and a state	000												
							Previous	8	Submit					

- ① Enter the corresponding sample information in the search bar for retrieval. Check the samples that need paternity analysis (ensure to check Sample ID and Batch No.), and click \rightarrow to add them to the Selected list.
- 2 Select the sample identity in the **Sample Type** column.
- ③ Click **Submit** to submit the paternity analysis task.
- 4. View paternity analysis results

		o Re-generale Report										
Na.	Batch No.	Sample ID	Relationship	Identity Status	Identity Result	Document Status	Gens Locus Type	Application	Additional Parameter	Document Update	Operator	Operation
1	FGID3409200001	Demo_F,Demo_M,Demo_C	Alleged Father, Alleged Mother, Allege	Finish	Inclusion	Finish	NFCR	NIPPT	Trio Paternity	2024-09-20 13:40	Admin	Data Download Preview Delete

- ① Click FGID > Applications > Identification Result to go to the analysis results interface.
- 2 Click **Preview** to view the analysis report.

3 Click **Data Download** to download the paternity analysis result files.

The analysis result file is a compressed file that includes a report document in .docx format and a result folder. The result folder contains the following contents:

- a. r0.paternity_results.txt: Indicates paternity analysis SNPs information, including parental SNP genotypes, fetal four-base depth, and PI results of each effective SNP.
- b. r0.paternity_results.xlsx: The content is same as the r0.paternity results.txt but in .xlsx format.
- c. r1.result.txt: It is the paternity analysis result file, referenced in the Word report.

d. r2.analysis.xlsx: It includes paternity analysis indicators

🖓 Tips

- To ensure the reliability of the results, check the r2.analysis.xlsx in the result folder for each paternity analysis result to confirm if any indicators are abnormal. For more detailed information, refer to the Q16 of FAO on Page 16.
- Paternity analysis results include Inclusion, Exclusion, Inconclusive, Inclusion-WrongMother, and WrongSample. If the result is WrongSample, it may indicate that the selected sample role identity is incorrect, and it is recommended to resubmit the paternity analysis task. For more detailed information, refer to the Q16 of FAQ on Page 16.

FAQ

Q1 What kind of laboratory environment is required for this experiment? How to maintain the experiment environment daily?

It is recommended to separately perform operations in the normal pre-PCR and post-PCR area and clean the laboratory regularly. Perform preparations in the pre-PCR area

- A before preparing libraries to reduce aerosol contamination. If you cannot separately perform operations, perform preparations in a clean bench.
- Q2 Are there any restrictions on the mother's samples?
 - This product is only suitable for cfDNA samples extracted from plasma of 8 weeks and above gestation, with a fetal concentration greater than or equal to 2%.
- This product is not suitable for twins or multiple pregnancy samples.
 - This product is not suitable for samples from mothers with tumor diseases, pre-eclampsia/eclampsia during pregnancy, or those who have received blood transfusions, bone marrow or organ transplants, or stem cell therapy.
- Q3 Are there any restrictions on the father's samples?
- A The extracted DNA concentration must be $\ge 0.2 \text{ ng/}\mu\text{L}$.

Q4 Are there recommended kits for DNA extraction?

• For gDNA extraction, it is recommended to use the MGIEasy genomic DNA extraction kit (Brand: MGI, Cat. No.: 1000010524, suitable for saliva, blood, and oral swabs), and for trace sample genomic DNA extraction (Brand: TIANGEN, Cat. No.:

- A DP316, suitable for samples with low DNA content like hair, semen, dry blood spots, chewing gum, and nails).
 - For cfDNA extraction, it is recommended to use the HiPure Circulating DNA Midi Kit (Brand: MAGEN, Cat. No.: D3182-03A).

Q5 Is the input amount for gDNA limited to 2 ng and cfDNA to 4 ng?

If the input amount is too large, it may lead to significant differences in sample output data amounts. If you need to increase the input amount for library preparation, it is

- A recommended to perform sequencing for only one family on a single DNBSEQ-G99 FCL flow cell.
- Q6 Can the sample preparation cartridge be inverted after loading the sealing liquid?

After loading the sealing liquid, the sample preparation cartridge should not be inverted or excessively tilted. Because this will cause the sealing liquid to leak, reduce the

A volume and fail to prepare library. The preparation cartridge cannot be used further. When transferring the sample preparation cartridge, keep it as horizontal as possible and do not tilt it.

- Q7 What is the size range of the target amplification fragments of the experiment?
- A The target insert fragment size range of the experiment is 60 bp to 90 bp and the target amplification fragment size of 144 bp to 174 bp.
- Q8 Can the library be sequenced with the NIPT library? How to do it?
- They can be sequenced together, but use different groups of barcode, make DNBs
- A separately, and mix the two kinds of DNBs in corresponding mass according to the required data amount.
- Q9 Is it necessary to remove the upper oil phase?
- A The upper oil phase of the DNB should be removed. Otherwise, the residual oil will affect sequencing quality.
- Q10 If the data come from other machine or the data are generated previously, they should be analyzed by using the server, and how to operate?

You can copy the fq files to the hard drive and transfer them to the server by using Filezilla. Perform the following steps:

- 1. Press down the server top cover, then lift it.
- 2. Insert the hard drive into the left USB port.
- Double-click the Filezilla icon on the desktop to open the software, and click File > Site Manager in the upper-left corner, or directly click the site manager shortcut. Select the Upload account in My Sites, and click Connect to log in to the Uploader account.
- 4. Input /mnt/das/Data in the right remote site, and select the sequencer ID (which may not match the actual sequence ID).
- A 5. Select the flow cell to upload in the left local site, with a folder directory of *flow cell number/L01*.
 - 6. Select the flow cell number to upload, right-click the flow cell number, and select **Upload**.
 - In the right remote site under the uploaded sequencer ID directory, click Info > Upload go to the path, right-click in the blank area, and select Create New File to create a flag file named as chip number_L01.txt.
 - 8. After uploading the fq file, you can submit the analysis task.

🕜 Tips

The flag file must be created; otherwise, the system will not start the analysis for that flow cell.

Part No.: H-020-001063-00
Version: 1.0

- Q11 If automatic analysis fails after sequencing, how to restart the analysis?
 - 1. Log in to the Uploader account.

Α

- 2. Enter /mnt/das/Data in the right remote site, select the folder corresponding to the sequencer number, and check if the flow cell file exists.
- 3. If the corresponding flow cell file exists, check if the task submission form
- is correct, or go to the system interface, click **Task** and the Batch No. of the corresponding task number to confirm if the task submission information is correct. After confirming that the information is correct, click **Task**, then click **Redo** in Operation to restart the analysis.
 - 4. If the corresponding flow cell file does not exist, refer to steps 3-8 in Q10 to upload the flow cell to the server.
- Q12 What is the minimum data amount for a single family, and how many families can be sequenced on each G99 flow cell?
- A single family requires a minimum of 3M Reads for gDNA and 15M Reads for cfDNA. Up to 3 families of sample can be sequenced on each G99 flow cell.
- Q13 What files are included in the exported basic data analysis results?

The exported basic data analysis results include infoQC.txt, r0.QualityControl.xlsx, r1.CompareSummary.xlsx, and the genotyping results and vcf files for each sample.

- infoQC.txt: It refers to basic data analysis results for the samples, including MappedRate, TargetOnMap, Coverage ≥ 100X, Coverage ≥ 1000X, Uniformity (0.1), and other data quality control results.
- r0.QualityControl.xlsx: The content is same as the infoQC.txt but in .xlsx format.
- A r1.CompareSummary.xlsx: It refers to the SNP comparison results between samples. When sample names contain "-", the first "-" is used as the group marker, and characters before "-" that are consistent are considered as the same group for SNP comparison. For example, for a-1-1, a-2-1, a_1-1, and a_1-2, a-1-1 and a-2-1 will be compared because they have "a-" which is regarded as the grouping marker. a_1-1 and a_1-2 will be compared because they have "a_1-" which is regarded as the grouping marker.
 - The folder named after the sample contains the genotyping results and vcf files.
- Q14 What are the principles for selecting SNP points for analysis? What is the origin of the frequency information? How to replace it with local population allele frequencies?

The SNP frequencies used come from the 1000 Genomes Project which utilizes the global frequency data. Each selected SNP has an allele frequency higher than 0.35

A with no linkage between SNPs. To replace it with local population frequencies, you can consult the technical support, provide the frequency information and let the technical support replace it.

Q15 Why are the basic data analysis results empty?

Α

In this case, the sample may lack of fq file. Check the corresponding barcode of the sample in the submitted task. To view the task, check the task submission form,

- A or click **Task** and the Batch No. of the corresponding task number to view the task submission information.
- Q16 What are the determination indicators for paternity analysis results?
 - Inclusion: The indicator is CPI > 1.0e+04 and Z-Score ≥ 3 , indicating that the likelihood of the alleged father being related to the fetus is > 99.999999%. When the result is Inclusion, check if Err(SNPs) is less than or equal to 1.00%. If the value is greater than 1.00%, contamination may exist in the experiment, causing unreliable results, and library repreparation is recommended.
 - Exclusion: The indicator is CPI < 1.0e-04 and Z-Score < 10, indicating that the likelihood of the alleged father being related to the fetus is < 0.00000001%.
 - Inclusion-WrongMother: The indicator is CPI < 1.0e-04 and Z-Score ≥ 10, indicating that the mother of this family may not be the biological mother of the fetus. When the result is Inclusion-WrongMother, check sample concentration, library concentration, sample storage time, and conditions.
 - WrongSample: The indicator is Fetal(SNPs) \geq 100% or Err(SNPs) \geq 10%. When Fetal(SNPs) is greater than or equal to 100%, it may indicate that the father's sample was analyzed as the mother. When Err(SNPs) is greater than or equal to 10%, it may indicate that the father's sample was analyzed as the fetus.
 - Inconclusive Range: The indicator is CPI \geq 1.0e-04 and \leq 1.0e+04 and situations that cannot be determined as the above four results. The indicator results may be resulted from a low number of SNPs available for analysis due to contamination in the sample or poor sample quality. You can check if the SNP values in PI(SNPs) are low.
- Q17 Why is the FGID_sample_template or FGID_task_template downloaded from the FGID system blank?
- A This may occur due to a caching delay in the system when downloading files. Please wait until the download is complete before opening the downloaded file.

Q18	How to create a new SE50+10 sequencing recipe?

Α

If the device's Software Version has been updated from a lower version to 1.7.0.674 or above, and the ISW Version is V1.3.0.568 or above, you need to delete the original sequencing recipe and customize it. Perform the following steps:

1. Tap the () > Settings on the upper-right corner of the main interface G99 main

interface, and select **Recipe** to go to the Sequencing recipe settings interface. In this interface, you can create, delete a sequencing recipe and adjust the display order of sequencing recipes.

- 2. Tap Create to go to the Create Recipe window.
- 3. Input SE50+10 in the box behind the Recipe name, 50 in the box under the Read1, 0 in the box under read2, 10 in the box under Dulbarcode, and 0 in the box under Barcode.
- 4. Tap Save. The Recipe interface returns to view.
- 5. Tap **Close**. The main interface returns to view.