
MGIEasy Small RNA Library Prep Kit V2.0

Features



Introduction

Small RNAs are an endogenous non-coding single-stranded RNA with a length of approximately 20-24 nucleotides and includes microR-NAs, siRNA and piRNA. They play a variety of important regulatory roles in cells such as inducing gene silencing, including RNA interference (RNAi), translation inhibition, heterochromatin formation and others which regulate various cellular physiological processes.

MGIEasy Small RNA Library Prep Kit provides an efficient workflow for generating libraries suitable for MGI high-throughput sequencing platforms with a broad range of input amounts (10 ng - 1 µg) of total RNA. Libraries prepared using the MGIEasy Small RNA Library Prep Kit can be used to detect known and unknown small RNAs and their target genes from human, animal and plant sample types.


Product Specification		
Assay Time	~13 hours	
Hands-On Time	~50 min	
Input Quantity	10 ng - 1 μg of total RNA	
Sample types	Tissues	
Species Compatibility	Human, animals, plants, such as mouse, rice and Arabidopsis	
Applications	Small RNA Sequencing	
Platform Compatibility	BGISEQ-500, MGISEQ-2000, DNBSEQ-G400*	
Recommended Read Length	SE50	
Recommended sequencing data per sample	25 M raw reads	

For Research Use Only. Not for use in diagnostic procedures.

Data Performance

High Reproducibility

MGIEasy Small RNA Library Prep Kit is compatible with a range of total RNA input (10 ng -1 µg) and exhibits excellent quality sequencing data and miRNA detection. The protocol exhibits high reproducibility and correlation across range of sample inputs (Fig. 2a). Sequencing data from samples of two input amounts (500 ng and 1 ug) were highly correlated when compared using Spearman and Pearson correlation coefficients (Fig. 1b).







Fig. 1a Total number of detected miRNAs from different RNA inputs.

Fig.1b Gene expression level reproducibility in different RNA input.

Fig. 1 Libraries were prepared from Human Brain Total RNA (HBR) using the MGIEasy Small RNA Library Prep Kit, and sequenced on MGI sequencer at SE50 read-length. Approximately 25M raw reads were collected per sample for analysis.

Compatible with various species

MGIEasy Small RNA Library Prep Kit can be widely applied to human, animal, plant and plasma sample. To demonstrate this human, mouse and rice total RNA samples were used to prepare sequencing libraries. The sequenced small RNA library data (Fig. 2) shows that the MGIEasy Small RNA Library Prep Kit is compatible with different species and has high technical reproducibility.





Libraries were prepared from HBR, mouse and rice total RNA using the MGIEasy Small RNA Library Prep Kit, and sequenced on MGI sequencer at SE50 read-length. Approximately 25M raw reads were collected per sample for analysis.

High sensitivity

A comparison was performed between the MGIEasy Small RNA Library Prep Kit sequenced on an MGI sequencer and other comparable Small RNA kits sequenced on a different platform. The results in Fig. 4 demonstrate that the libraries constructed using MGIEasy Small RNA Library Prep Kit detect more miRNA with the same amount of data (Fig. 3).




Fig. 3 Number of detected miRNA using different sequencing platforms.

The MGI libraries were generated from Human Brain Total RNA using the MGIEasy Small RNA Library Prep Kit, and sequenced on a MGI sequencer at SE50 read-length. The other libraries were prepared with Human Brain Total RNA using others kits and sequenced on "N" platform at SE50 read-length. Approximately 25M raw reads were collected per sample for analysis.

Summary

MGIEasy Small RNA Library Prep Kit can efficiently prepare Small RNA libraries from human, animal and plant samples with a range (10 ng -1 µg) of total RNA input. This enables excellent quality sequencing data and outstanding sensitivity to detect microRNAs with high technical reproducibility. MGIEasy Small RNA Library Prep Kit provides two options of size selection based either on beads or gel purification to meet the needs of most species types.

Ordering information			
Product	Configuration	Catalog No.	
MGIEasy Small RNA Library Prep Kit	24 RXN	1000005269	
MGISEQ-2000RS High-Throughput Sequencing Set (Small RNA)	SE 50	1000006138	

Contact Us

MGI Tech Co., Ltd Add : Building11, Beishan Industrial Zone, Yantian District, Shenzhen, CHINA 518083 Email: MGI-service@mgi-tech.com Website: en.mgi-tech.com Tel: 4000-688-114 Version: November 2022 | MGPD111810200-06



5

https://www.linkedin.com/company/mgi-bgi



https://twitter.com/MGI_BGI

Copyright © 2018 MGI, Inc. All Rights Reserved. The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating, or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd.. and their respective producers.