



# An Ultra-High Throughput Sequencer

DNBSEQ-T20×2

For Research Use Only.



Say Hello to

**DNBSEQ-T20x2**





# An Ultra-High Throughput Sequencer

DNBSEQ-T20×2 is a fully automated sequencer comprising all components of a genetic sequencer integrated by robotics to process six slides per run. DNBSEQ-T20×2 achieves remarkable throughput far greater than existing "high-throughput" sequencers, and can drive down sequencing cost to a sub \$100 price per genome as a trailblazer.

**>50,000 WGS\*/Year**

\*30× Genome Coverage

## Throughput Improvement

# Up to 10.5×

PE100



4 Tb

# 42 Tb

output per run

# Up to 4.5×

PE150\*



16 Tb


# 72 Tb

output per run

\*PE150 is estimated to be released in 2023 Q3.

# Trailblazing Path To Sub Hundred-dollar Genome

DNBSEQ-T20×2 is designed to meet the most challenging sequencing scenarios. It aims to empower large population genome projects around the world, and revolutionize healthcare with OMICS technologies. DNBSEQ-T20×2 will be **the perfect choice** to bring sequencing efficiency and scale to the next level.



# Product Components

## 6-axis Robotic Arm

Transfers slides between different modules

## Power Module

Supplies electricity and compressed air power

## BIT Module

Generates FASTQ data or performs advanced bioinformatics analysis

## Biochemical Reaction Module

Performs biochemical reactions on slides

## Reagent Module

Stores unreacted reagents

## Storage Module

Stores slides and reacted reagents

## Optical Imaging Module 1

Captures labeled fluorescent signal

## Optical Control Module 1

Controls each component of optical imaging module 2 and transfers fluorescent signal to CAL file

## Optical Imaging Module 2

Captures labeled fluorescent signal

## Optical Control Module 2

Controls each component of optical imaging module 2 and transfers fluorescent signal to CAL file

# Six Super-sized Slides Provide Excellent Performance

DNBSEQ-T20×2 supports **different sequencing applications to be run on different slides at the same time**. Glass covers have been removed from the DNBSEQ-T20×2 slides to facilitate the dip-immersion mode for biochemistry reactions. This allows reagents to be re-circulated for reuse, while achieving high uniformity and stability of reactions.

**PE100**

Read Length

**35 B**

Effective Reads/Slide

**42 Tb**

Data Output/Run

**2.5 d<sup>++</sup>**

Sequencing Time

**≥85**

Q30 (%)<sup>+</sup>

**PE150\***

Read Length

**40 B**

Effective Reads/Slide

**72 Tb**

Data Output/Run

**3.5 d<sup>++</sup>**

Sequencing Time

**≥80**


Q30 (%)<sup>+</sup>

\*PE150 is estimated to be released in 2023 Q3.

<sup>+</sup>The percentage of base above Q30 is the average of an internal standard library over the entire run.

The actual performance is affected by factors such as sample type, library quality and type, and insert fragment length.

<sup>++</sup>The run time includes slide loading, sequencing, and generation of Cal. file. Cal. file is a binary file format generated by DNBSEQ™ platform basecall software.



Slides



# Innovative Dip-Immersion Biochemistry

DNBSEQ-T20×2 adopts an innovative **dip-immersion biochemical technology**. In contrast to the flow cell, the dip-immersion biochemical reaction is carried out by consecutively immersing multiple sequencing slides (no cover) in the same reagent container, enabling reagent reuse and uniform biochemistry on wafer-sized sequencing arrays.



# Powered by DNBSEQ™

Leveraging the core advantages of **DNBSEQ™** combined with **highly accurate base-calling and copy number correction algorithms**, DNBSEQ-T20×2 achieves outstanding sequencing performance and exceptional data quality with low duplication rates and ultra-high effective data output.



INCREASED ACCURACY



DECREASED DUPLICATES




REDUCED INDEX HOPPING

DNA single strand  
circularization

DNB Making

Patterned Array

DNB Loading  
(DNA nanoballs)



# Shape The Future Of Ultra-High Throughput Sequencing

DNBSEQ-T20×2 covers all high-throughput sequencing needs in scientific research and clinical research, including: WGS, WGBS, stLFR, single-cell sequencing, stereo-seq, and etc.

\*PE150 is estimated to be released in 2023 Q3.

# Estimated Sample Throughput For Key Applications

Application	Sample Throughput	Reads/Sample	Data/Sample	Depth/Sample	Application filed
WGS	~ 336/Run	500 M	100 Gb	30X	Scientific research
cWGS	~ 240/Run	700 M	140Gb	40X	Clinical research
stLFR WGS	~ 168/Run	750 M	150 Gb	30X	Scientific research
WGBS	~ 168/Run	550 M	110 Gb	30X	Scientific research
Single-cell Sequencing	~ 192/Run	800 M 10,000 cells/sample	/	/	Scientific research
Stereo-seq	~ 144/Run	1000 M	/	/	Scientific research

\*All sample throughputs are estimated by PE100.


\*\*The stLFR Library Prep Kit is Only for Sale in the United States, and for Research Use Only.

# Start Sequencing Immediately

We offer a **complete package** with the DNBSEQ-T20×2, including sample preparation reagents and equipment, automated library preparation systems, and ultra-high throughput sequencers. In addition, a range of tools such as massive data management systems and processing modules are available to users.



# The One-stop Package



# Full Portfolio of DNBSEQ™ Platform





# Hardware Specification

Parameter	Specification
Models	Genetic Sequencer DNBSEQ-T20×2RS
Dimensions	Biochemical Reaction Module: 1750 mm (L)* 1220 mm (W) * 2000 mm (H) DNBSEQ-T20×2RS: 4200 mm (L)* 4800 mm (W) * 2000 mm (H)
Net Weight (without BIT module)	DNBSEQ-T20×2RS: ~2743 kg
Touch Screen	Touch Screen Type: LCD Touch Screen Size: 21.5 inch Touch Screen Resolution: 1920×1080
Power	Rated Voltage: 3~, 380/400 V Rated Frequency: 50 Hz Rated Power: 20000 VA
Operating Environment Requirements	Temperature: 19 °C~22 °C Relative Humidity: 30% RH~70% RH Atmospheric Pressure: 90 kPa~106 kPa
Shipping/Stroage Environment Requirements	Temperature: 15 °C~30 °C Relative Humidity: 30% RH~70% RH Atmospheric Pressure: 90 kPa~106 kPa

\*The maximum sound pressure of DNBSEQ-T20×2 is 75 dB (A).

\*\*The shell protection grade of DNBSEQ-T20×2 is IPX0.

# Ordering Information

Cat.No.	Product Name
900-000544-00	Genetic Sequencer DNBSEQ-T20x2RS (CHN-RUO)
900-000543-00	Genetic Sequencer DNBSEQ-T20x2RS (CE-RUO)
900-000667-00	Genetic Sequencer DNBSEQ-T20x2RS (US-RUO)
900-000546-00	DNB Loader MGIDL-T20ARS (CHN-RUO)
900-000546-00	DNB Loader MGIDL-T20ARS (CE-RUO)
900-000668-00	DNB Loader MGIDL-T20ARS (US-RUO)
940-000350-00	DNBSEQ-T20x2RS High-throughput Sequencing Set


\* Unless otherwise informed, StandardMPS and CoolMPS sequencing reagents, and sequencers for use with such reagents are not available in Germany, Spain, UK, Sweden, Belgium, Italy, Finland, Czech Republic, Switzerland, Portugal, Austria and Romania. Unless otherwise informed, StandardMPS sequencing reagents, and sequencers for use with such reagents are not available in Hong Kong.


MGI Tech Co., Ltd

Building 11, Beishan Industrial Zone,  
Yantian District, Shenzhen. CHINA 518083

Version: February 2023

Information in this brochure is updated to [02/24/2023] and only for your reference. In no event shall the brochure be regarded as warranty or commitment made by MGI Tech Co., Ltd. All rights and obligations shall be subject to the final executed agreement.

 [en.mgi-tech.com](https://en.mgi-tech.com)

 [MGI-service@mgi-tech.com](mailto:MGI-service@mgi-tech.com)

