

» PFI

One click Pathogen Fast Identification Server



» Product Introduction

PFI is a solution that combines the functions of software and hardware. It involves a large microorganism database with analyzing function, and is able to manage the whole process including inputting samples and generating reports via ZLIMS. The entire system is loaded on the server of Bioinformatics analysis accelerator.

PFI can rapidly, accurately and comprehensively identify microorganisms from the original samples at a metagenome/meta-transcriptome level, and automatically generate identification reports and analytical results. It can also provide assistance and references for the diagnosis and treatment of infectious diseases.


» Database

The PFI system collects more than 27,000 microbial genomic sequences (the information of species that have multiple reference genomes is also included) for the database to provide rapid and precise detection.

Microbial classification	Species	Genus
Bacteria	8335+	1925+
Archaea	361+	130+
Viruses	9461+	2606+
Fungi	9031+	1050+
Parasites	642+	361+

» Workflow

- Format of the input data : FASTQ



» Total analysis time

No. Reads/sample	No. Samples/time	Read lengths	Total analysis time (hour)
100M	1	SE50	0.5
	1	SE100	0.85
	1	PE50	0.92
	1	PE100	1.53
100M	5	SE50	1.18
	5	SE100	1.77
	5	PE50	3.23
	5	PE100	2


» The server

Highly-stable and reliable HP Z8 G4 workstation



- CPU: Intel Xeon Gold 6240
- SSD: 2TB+256GB
- Memory: 128GB DDR4
- HDD: 30TB 7200 SATA 3.5 inch

» Our advantages



Simplified operation process

Web system, one-click operation, easy to use, no need of Linux command and programming.

With Zlms to enable the one-stop automatic operation (sequencing and analyzing); totally automatic, no manual calculating process is needed.



Large database

Host database: reference sequences of human beings, pigs, goats, sheep, mice, rats, carp, domestic geese, chicken, ducks, cattle, cats, dogs, rabbits and other common animals; automatically filter those sequences or customized according to customer needs.

Microbial database: entire NCBI RefSeq genome sequences, including bacteria (including archaea), viruses, fungi, and parasites.



Fully functional analysis

Able to undertake analysis on RNA meta-transcriptome sequencing and DNA metagenome sequencing independently, RNA meta-transcriptome sequencing and DNA metagenome sequencing at the same time.

Solve the problem of background microorganisms' interference, identify the samples more precisely, find the active pathogenic microorganisms.



A wide range of samples

- Including environmental and clinical samples
- Compatible with samples from various sources such as blood of human or animals, alveolar lavage fluid, intestines, etc.

Pathogen Fast Identification helps understand the diversity of samples by sequencing all nucleotides from both host and microbes. As this method does not require preliminary knowledge of pathogenic microbial genomes, it can identify unknown pathogens in infectious disease. Importantly, the unique identification technique supports to develop strategies to control and prevent human and animal infectious diseases.

> Overview

> Solution

> Result

Complete Genome Identify Sketch Map

Human Enterovirus C96

3 nt gap 13 nt gap 24 nt gap

0k 2k 4k 6k 7.4k ← Genome

Human Poliovirus 1

34 nt gap

0k 2k 4k 6k 7.4k ← Genome

Human coxsackievirus A24

0k 2k 4k 6k 7.4k ← Genome

Metagenomics sequence

ABI 3730 sequence

Fig.1 Complete genome identify sketch map. A total of 10 contigs were assembled after metagenomic sequencing (yellow-highlighted bars). Sequencing was performed with an ABI 3730 genetic analyzer.

the sequence was verified, and the gaps were filled with overlap for the whole genome (green-highlighted bars) (Color figure online).

> Paper

A novel Enterovirus 96 circulating in China causes hand, foot, and mouth disease
published on Virus Genes on February 7th, 2017



Case 2 Identification of co-infections

> Overview

A 26-year-old woman developed a mild respiratory illness on Jan. 28, 2017, but symptoms progressed to recurrent fever, cough, chills, expectoration, slight hemoptysis, muscle and joint pain in the following days. On February 3, the patient was hospitalized due to worsening symptoms of cyanotic lips, fever of 39.3 ° C, heart rate of 144 beats/ min and diagnosed with severe pneumonia with ARDS. She was treated with antibiotics and antiviral therapies and then discharged on February 17. Blood and respiratory secretions were collected during her hospitalization for pathogen testing. The screening results from bacteria and fungi culture-based test, G-test and GM-test were all negative. In addition, HIV, HBV, influenza viruses, SARS-CoV, MERS-CoV and other coronaviruses were negative by ELISA and/or (RT-)PCR assays.

> Solution

Pulmonary secretions from the patient were collected on the first day of hospitalization and analyzed using metagenomic sequencing to determine the cause of infection.

> Result

Two respiratory viruses, HRV and HBoV were identified in high abundance using sequencing and confirmed by specific (RT-)qPCR assays and a report generated to diagnose acute co-infection of HBoV1 and HRV-C.

In this case, metagenomic sequencing showed a significant advantage in detection of the causative agents of severe illness over traditional methods such as culture, ELISA, PCR, etc. because prior knowledge was not assumed or required.




Figure A. Viral reads distribution.

> Paper

In a case published on Journal of Infection in March 2018 of an adult suffering from severe pneumonia.



Case 3 Diagnosis of rare pathogen

> Overview


On August 20, 2018, A 42-year-old man was hospitalized after presenting with symptoms of severe head-ache, fever of 38.4 ° C, and elevated CSF leukocyte and protein levels. 24 hours post-presentation, the patient spoke incoherently, had breathing difficulties, became comatose and was subsequently transferred to ICU. Further examination by CT scan showed hydrocephalus and brain edema. Four days later, a culture from cerebrospinal fluid samples showed negative results for bacteria and fungi, therefore, to identify the pathogen, the sample was further analyzed using high-throughput sequencing on August 31. Results were reported to clinicians 2 days later.

> Solution

The cerebrospinal fluid sample collected from the patient was analyzed on the MGI sequencing platform using the Pathogen Fast Identification workflow.

> Result

Interestingly the results revealed a low level (0.0266%) of sequencing reads that identified as *Naegleria fowleri*, a rare amebic pathogen that can cause primary amebic meningoencephalitis (PAM).



PAM caused by *Naegleria fowleri* infection is extremely rare in China but almost always fatal. The patient went to the Songkran Festival prior to the onset of illness and may have come into contact with sewage. In this case, traditional methods failed to detect the pathogen, however, the MGI high-throughput sequencing platform successfully identified the rare pathogen.

Case 4 Public health issue

> Overview

A 45-year-old male returning to China from Angola showed symptoms of Rift Valley fever including fever (38.8 ° C), chills, headache, arthralgia, anorexia and enervation on July 13 and was hospitalized for treatment. BGI assisted the Entry-exit Inspection and Quarantine of China to obtain a whole genome sequence of Rift Valley fever virus from the individual using NGS technology. As a result, BGI helped identify, quarantine and treat the individual and prevent a local outbreak of RVF in China.

> Solution


RVFV isolation and culture identification were done in biosafety lab of Guangdong Inspection and Quarantine Technology Center. BGI laboratory performed high-throughput sequencing of the sample to gather genomic information about RVFV.

> Result

Alignment of the full genome sequence of the RVFV isolate (named RVFV-Beijing strain) revealed 100% identity of three gene segments and 98% homologous to RVFV Kakamas isolate in South Africa.

> Paper

Isolation and phylogenetic study of Rift Valley fever virus from the first imported case to China
Published on VIROLOGICA SINICA in June 2017



VIROLOGICA SINICA 2017, 32 (3): 253-256
DOI: 10.1007/s12250-017-3949-2

LETTER

Isolation and phylogenetic study of Rift Valley fever virus from the first imported case to China

Dear Editor,

Rift Valley fever (RVF) is an anthroponosis caused by Rift Valley fever virus (RVFV). RVFV belongs to the *Phlebovirus* genus in the family *Bunyaviridae*, which is circulating among ruminants. Human infection with RVFV is generally asymptomatic; however, minority of patients develop severe RVF diseases like encephalitis or hemorrhagic complications, with a fatality rate up to 50%. RVF is mostly distributed in eastern and southern Africa, including Kenya, Zimbabwe, Zambia, Namibia and Somalia (Liljequist et al., 2011; Shivanya et al., 2016). Infected cases were also reported in Egypt, Saudi Arabia and Yemen (Hegans and Malins, 2011). Large outbreaks of RVF occurred in Kenya, Madagascar and South Africa with rapid spread (from four September to October 2015). In recent years, the import and export trade of animals increased remarkably along with the development of transportation and tourism, which raise the risk of RVFV transmission. The epidemic in Niger caused wide spread concern since the first case in August, 2016. With regard to China, RVFV had not been reported until July 23, 2016, when a Chinese worker returning from Angola was diagnosed with RVFV infection and was confirmed to be the first imported case of RVF in China (World Health Organization, 2016).

Here, we describe the RVFV isolation in serum sample and its complete genome from the first imported case to China. The 45-year-old male returning to China from Angola showed symptoms of fever (38.8 °C), chills, headache, arthralgia, anorexia and enervation on July 13 and returned to China for treatment on July 21 through the port of Beijing Airport. Because of yellow fever outbreak in Angola, differential diagnosis of more than ten pathogens in the lab were carried out, including yellow fever virus, dengue virus, Chikungunya virus and Psittacosis. The real-time fluorescent RT-PCR was performed with specific primers (RVFV-GAAATCTCTT-GAGACAC-ATGGAT; RVFV-TTCACCTCTCTT-TTCACCTCTT) and primers (RVFV-FAMCACA-GTCCACACGGCCCTTACATT-BHQ1) to detect RVFV, and the serum sample was determined to be RVFV positive.

RVFV isolation and culture identification were done in Biosafety Level 3 lab of Guangdong Inspection and Quarantine Technology Center. Vero and BHK-21 cells were used for virus isolation from serum according to the reference (Baron et al., 2014). The pre-treated serum sample was inoculated into monolayer Vero and BHK-21 cells in six-well plate. Cells were cultured at 37 °C in 5% CO₂ for up to 7 days. Obvious cytopathic effect (CPE) was observed on Vero and BHK-21 cells 72 h post-inoculation. Vero cells showed CPE of swelling, particle increasing, disruption and necrosis, with BHK-21 turning round, shrinkage and rupture (Figure 1a). The supernatant of infected Vero and BHK-21 cells was collected 5 days post inoculation and virus RNA was extracted. The real-time fluorescent RT-PCR and RT-PCR were performed to detect RVFV proliferation in cells (Figure 1b, c). RT-PCR showed two amplified bands of 564 bp for L segment (RVFV-LP: GGATCATATCTGAAATCTCTTCTCTT; RVFV-RP: CAAGTCAACAAGATCTCTTCTCTT) and 835 bp for S segment (RVFV-SFP: CTGACATCTGGGATTGGAG; RVFV-RP: CCTCTGTGTTCTCTCTCTT) (Figure 1c).

The viral RNA of RVFV extracted from cell culture supernatant was subjected to next generation sequencing according to the in-house SOP of BGI library construction for BGISEQ500 (BGI-Shenzhen, China). About 13.63 Gb short reads were obtained and analyzed using Kraken system (Wood and Salzberg, 2014), most of which are from host cells or unassigned. To obtain the genome of RVFV, all short reads were first mapped to whole genomes of *Phlebovirus* retrieved from GenBank (until Aug., 2016) with a match rate threshold at 0.9 according to the coverage mapped by short reads. Then the consensus genome sequences were called from the dominant bases at each nucleotide site. Full length RVFV genome sequences could be obtained with an average depth at 1.5 × from 100 M SE100 reads, which were submitted to GenBank under accession No. of KX632066, KX632067 and KX632068, corresponding to L, M and S segment respectively. The length of L segment is 6402 bp with the region of 17 nt-629 nt including the RNA polymerase (L protein). M segment contains 3872 bp and

Case 5 Diagnosis of animal disease

> Overview

A group of goats, infected with unknown pathogens developed scabby lesions around their lips, muzzle, and in their mouth. To efficiently control the unknown infection, throat swab samples from the affected goats were tested using MGI Pathogen Fast Identification system with a report generated within two days.

> Solution

The throat swab samples from goats were processed using the MGI Pathogen Fast Identification system with automatic DNA and RNA sample extraction.

> Result

A large proportion of sequencing reads in both DNA (64.2%) and RNA (44%) samples mapped directly to the Orf virus which associates closely with the clinical symptoms presented.

Table 1 The microorganism identification of DNA and RNA sample from goat swab

Rank	name of pathogen	DNA		RNA	
		reads number	relative abundance	reads number	relative abundance
1	Orf_virus	358593	64.20%	11311	44%
2	Pseudocowpox_virus	26658	4.80%	650	2.50%
3	Bacillus_subtilis	3130	0.60%	518	2.00%
4	Pseudomonas_aeruginosa	2011	0.40%	217	0.80%
5	Staphylococcus_aureus	1158	0.20%	220	0.90%

Direct comparison of the obtained sequencing reads to reference genome of Orf virus genome showed 86.7% identity, 87.6% average coverage and 200X depth. (see Figure below)




Fig.1 Alignment linear graph of assembled sequence and viral genome




Fig.2 Average coverage of viral genome




Fig.3 Average sequencing depth of viral genome

To verify the result, a traditional PCR assay was then performed and it confirmed the accuracy of the result. The MGI sequencing technology is highly accurate for pathogen identification and capable of rapid diagnosis and treatment of animal diseases.

» MGI Global Presence

• China



Changchun R&D centre	Shenzhen R & D Center Production Center Training Center Logistics Centre After-sale Center	Wuhan R & D Center Production Center Training Center Logistics Centre After-sale Center Quality Center Qualification Center	Hong Kong Logistics Centre After-sale Center
--------------------------------	--	---	---

• Latvia



Riga Qualification Center After-sale Center Training Center Logistics Centre

• UAE



Dubai After-sale Center

• Japan



Kobe After-sale Center

• US



San Jose R & D Center Production center After-sale Center

» Ordering Information

Product	Part No.
Fast identification of pathogenic infections (with server)	510-000164-00

MGI Tech Co., Ltd.

Service & Support

MGI has accumulated rich experience in gene sequencing with an excellent team of scientists and engineers, who are committed to providing comprehensive technical support in each section: from the installation, testing and operation, training, maintenance to subsequent upgrades, as well as the laboratory system construction, experiment scheme design and sequencing data analysis. You will experience an unprecedented journey of sequencing.

Contact Us

Address.: Building 11, Beishan Industrial Zone, Yantian District, Shenzhen, China 518083

Email: MGI-service@mgi-tech.com

Website: en.mgi-tech.com

Tel: 4000-966-988

<https://www.linkedin.com/company/mgi-bgi>

@MGI_BGI



Wechat

For Research Use Only

Copyright Disclaimer: The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating, or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

Version: January 2021 | MGPE2003002