



Part No.: H-020-000720-00

# MGIEasy

## Total RNA Extraction Set

### Instructions for Use

Version: 4.0

Leading Life Science Innovation

Address: Building 24, Stage 3.1, BioLake Accelerator, No.388, 2nd Gaoxin Road,  
East Lake High-Tech Development Zone, 430075, Wuhan, P.R. China  
E-mail: [MGI-service@mgi-tech.com](mailto:MGI-service@mgi-tech.com)  
Website: [en.mgi-tech.com](http://en.mgi-tech.com)

Research Use  
Only

**Wuhan MGI Tech Co., Ltd.**

---

## About the instructions for use

This instructions for use is applicable to MGIEasy Total RNA Extraction Set. The version of the instructions for use is 4.0 and the set version is 1.0.

This instructions for use and the information contained herein are proprietary to Wuhan MGI Tech Co., Ltd. (hereinafter referred to as MGI), and are intended solely for the contractual use of its customers in the use of the products described herein and for no other purpose. This instructions for use and its contents shall not be reprinted, reproduced, modified, distributed, or disclosed to others, in whole or in part, without prior written consent from MGI.

MGI makes no commitment to this instructions for use, including (but not limited to) any special commercial purpose and any reasonable implied warranties. MGI has taken measures to ensure the correctness of this instructions for use. However, MGI is not responsible for any missing parts of the instructions for use and reserves the right to revise the instructions for use and modify the set, so as to improve the reliability, performance or design.

Figures in this instructions for use are for illustrative purpose only. The content may be slightly different from the set. For the most up-to-date details, refer to the set purchased.

MGIEasy™ is the trademark of MGI, or its subsidiaries in China and/or other countries. Other names and trademarks mentioned in this instructions for use may be the property of their respective companies and subsidiaries.

©2023-2024 Wuhan MGI Tech Co., Ltd. All rights reserved.

---

# Revision history

Version	Date	Description
4.0	July 4, 2024	<ul style="list-style-type: none"><li>Updated the volume of Wash Buffer II</li><li>Added the sample requirements and operations for MGISP-NE32RS and MGISP-NEXRS extraction</li></ul>
3.0	December 25, 2023	Updated the operations
2.0	October 27, 2023	Updated instructions and tips for 4.1, 4.2, and 4.3
1.0	May 15, 2023	Initial release

# Contents

<b>Chapter 1 Introduction</b>	<b>1</b>
1.1 Product name	1
1.2 Specifications	1
1.3 Intended use	1
1.4 Working principle	1
1.5 Main components	1
<b>Chapter 2 Applicable device</b>	<b>3</b>
<b>Chapter 3 Sample requirements</b>	<b>3</b>
3.1 Applicable sample	3
3.2 Sample amount requirements	3
3.3 Sample storage	4
3.4 Sample transportation	4
3.5 Sample safety	4
<b>Chapter 4 Operation</b>	<b>4</b>
4.1 Preparing materials	4
4.2 Pretreating samples	5
4.2.1 Human whole blood	6
4.2.2 Animal tissue	6
4.2.3 Cell	7
4.2.4 Bacteria	7
4.2.5 Yeast	8
4.3 Extracting the nucleic acids	9
4.3.1 Extracting the nucleic acids manually	9
4.3.2 Extracting the nucleic acids automatically on MGISP-960RS	11
4.3.3 Extracting the nucleic acids automatically on MGISP-NE384RS	15

.....	
4.3.4 Extracting the nucleic acids automatically on MGISP-NE32RS	20
.....	
4.3.5 Extracting the nucleic acids automatically on MGISP-NEXRS	26

---

<b>Chapter 5 Warnings and precautions</b>	<b>34</b>
---	-----------

---

<b>Appendix 1 Manufacturer information</b>	<b>34</b>
--	-----------

# Chapter 1 Introduction

## 1.1 Product name

MGIEasy Total RNA Extraction Set

## 1.2 Specifications

Set name	Model	Component	Cat. No.	Specification
MGIEasy Total RNA Extraction Set Cat. No.: 940-000880-00	MRT96	MGIEasy Total RNA Extraction Kit	940-000877-00	96 RXN/Kit
		DNase I	940-000879-00	
MGIEasy Total RNA Extraction Set Cat. No.: 940-000875-00	MRT384	MGIEasy Total RNA Extraction Kit	940-000878-00	384 RXN/Kit
		DNase I	940-000876-00	

## 1.3 Intended use

This set is used to extract high-quality and high-purity total RNA from cell, animal tissue (fresh or frozen), blood and other samples.

## 1.4 Working principle

By using this product, salt ions with high concentration lyse and release RNA from the animal cell, animal tissue (fresh or frozen at -80 °C) and blood samples. The released RNA is then captured by magnetic beads and washed by specific wash buffer to remove proteins, salt and other impurities. After being dried, the RNA in magnetic beads is eluted by elution buffer and high-purity total RNA is obtained.

## 1.5 Main components



**Tips** To avoid frequent freeze-thaw cycles, Buffer RDD can be stored at room temperature not exceeding 30 °C.

**Table 1 MGIEasy Total RNA Extraction Set (MRT96) Cat. No.: 940-000880-00**

Kit name	Component	Specification	Storage condition	Validity period	Transportation condition
MGIEasy Total RNA Extraction Kit Cat. No.: 940-000877-00	Buffer LY	29 mL/tube×1	2 °C to 30 °C	12 months	2 °C to 30 °C
	Buffer WB I	81 mL/tube×1			
	Buffer WB II	41 mL/tube×1			
	RNase Free Water	15 mL/tube×1			
	Proteinase K	2 mL/tube×1			
	Magnetic Beads T	6 mL/tube×1			
	Buffer LYR	168 mL/tube×2			
DNase I Cat. No.: 940-000879-00	DNase I	0.8 mL/tube×1	-25 °C to -15 °C		-25 °C to -15 °C
	Buffer RDD	15 mL/tube×1	-25 °C to 30 °C		-25 °C to 30 °C

**Table 2 MGIEasy Total RNA Extraction Set (MRT384) Cat. No.: 940-000875-00**

Kit name	Component	Specification	Storage condition	Validity period	Transportation condition
MGIEasy Total RNA Extraction Kit Cat. No.: 940-000878-00	Buffer LY	116 mL/tube×1	2 °C to 30 °C	12 months	2 °C to 30 °C
	Buffer WB I	323 mL/tube×1			
	Buffer WB II	162 mL/tube×1			
	RNase Free Water	60 mL/tube×1			
	Proteinase K	8 mL/tube×1			
	Magnetic Beads T	24 mL/tube×1			
	Buffer LYR	672 mL/tube×2			
DNase I Cat. No.: 940-000876-00	DNase I	0.8 mL/tube×4	-25 °C to -15 °C		-25 °C to -15 °C
	Buffer RDD	61 mL/tube×1	-25 °C to 30 °C		-25 °C to 30 °C

## Chapter 2 Applicable device

- MGISP-960RS High-throughput Automated Sample Preparation System
- MGISP-NE384RS Automated Nucleic Acid Extractor
- MGISP-NE32RS Automated Nucleic Acid Extractor
- MGISP-NEXRS Automated Nucleic Acid Extractor

## Chapter 3 Sample requirements

### 3.1 Applicable sample

This product is applicable to samples from cultured eukaryotic cell, solid tissue and blood of human or animals, and prokaryotic cell such as G+ bacteria and G- bacteria.

### 3.2 Sample amount requirements

**Table 3 Required sample volume for different extraction methods**

	Human whole blood	Animal tissue		Cell	Bacteria	Yeast
		Tissue of liver, spleen and kidney	Tissue of heart and lungs			
Manual extraction	100 $\mu$ L to 200 $\mu$ L	1 mg to 20 mg	5 mg to 15 mg	$1 \times 10^5$ to $5 \times 10^6$	$5 \times 10^7$ to $5 \times 10^9$	$5 \times 10^7$
MGISP-960RS	100 $\mu$ L to 200 $\mu$ L	1 mg to 30 mg	2 mg to 5 mg	$1 \times 10^5$ to $2.5 \times 10^6$	$5 \times 10^7$ to $5 \times 10^9$	$5 \times 10^7$
MGISP-NE384RS	100 $\mu$ L to 200 $\mu$ L	5 mg to 30 mg	5 mg to 20 mg	$1 \times 10^5$ to $1 \times 10^6$	$5 \times 10^7$ to $5 \times 10^9$	$5 \times 10^7$
MGISP-NE32RS	100 $\mu$ L to 200 $\mu$ L	1 mg to 20 mg	1 mg to 5 mg	$1 \times 10^5$ to $1 \times 10^6$	$5 \times 10^7$ to $5 \times 10^9$	$5 \times 10^7$
MGISP-NEXRS	100 $\mu$ L to 200 $\mu$ L	5 mg to 30 mg	5 mg to 20 mg	$1 \times 10^5$ to $1 \times 10^6$	$5 \times 10^7$ to $5 \times 10^9$	$5 \times 10^7$



**Tips** The amount of human whole blood sample should be adjusted according to the white blood cells count in the blood.



---

### 3.3 Sample storage

- For human whole blood samples, it is recommended to use EDTA-anticoagulated whole blood or sodium citrate-anticoagulated whole blood.
- For tissue samples, it is recommended to use fresh tissue or frozen tissue at -80 °C for up to 3 months.
- For cell or bacteria samples, it is recommended to use fresh precipitate cell or bacteria without culture medium or use frozen cell or bacteria at -80 °C for up to 6 months.
- Do not freeze and thaw frozen samples frequently. Otherwise, the RNA quality may decrease.
- Please thaw and mix the frozen samples thoroughly before use.
- Before use, take out all components in the reagent set, equilibrate to room temperature (10 °C to 30 °C) and mix them thoroughly before adding to wells. If solid object appears, heat the reagent at 50 °C to redissolve it, which does not affect the reagent's extraction effect.

---

### 3.4 Sample transportation

Use the dry ice for transportation for up to 7 days. During transportation, avoid frequent freeze-thaw cycles.

---

### 3.5 Sample safety

- All samples are regarded potentially infectious.
- All samples should be extracted after being inactivated according to relevant national regulations.

---

## Chapter 4 Operation

---

### 4.1 Preparing materials

Prepare the following materials:

**Table 4 User-supplied materials**


Type	Item	Description
Equipment	Mini centrifuge	With a speed no less than 10000 rpm
	Plate centrifuge	None

Type	Item	Description
Equipment	Tubular centrifuge	None
	Vortex mixer	None
	Thermomixer compact	It can be replaced by a water bath
	1.5 mL magnetic rack	None
	Pipette	1 mL/200 $\mu$ L/20 $\mu$ L
	Grinding media	3 mm, zirconia, RNase-free
	Grinding mill	-10 $^{\circ}$ C , low-temperature
Reagent	Absolute ethanol	Analytically pure
	Lysozyme	Recommended brand: TIANGEN
	PBS solution	Recommended brand: Sangon Biotech
	DEPC	Recommended brand: Sangon Biotech
	1xTE Buffer	PH 8.0. Recommended brand: Sangon Biotech
	TritonX-100	Recommended brand: Aladdin
	Lysostaphin	Recommended brand: Coolaber
	$\beta$ -Mercaptoethanol	Recommended brand: Aladdin
Consumables	Pipette tips	None
	Centrifuge tube	<ul style="list-style-type: none"> <li>50 mL/1.5 mL/0.5 mL</li> <li>DNase-free, RNase-free</li> </ul>
	Tips	1 mL/200 $\mu$ L/20 $\mu$ L
	Glass beads, acid-washed	Recommended brand: Magen

## 4.2 Pretreating samples


It is necessary to pretreat samples before nucleic acid extraction.

## 4.2.1 Human whole blood

 **Tips** The human whole blood here represents the EDTA-anticoagulated whole blood or sodium citrate-anticoagulated whole blood. It is recommended to extract nucleic acids directly from the human whole blood. If you need to store the human whole blood for a long time, pretreat it whose volume should be more than 200  $\mu\text{L}$ , and store it at  $-80\text{ }^{\circ}\text{C}$  for up to 1 month. Blood is complicated in structure, so during blood pretreatment, RNA from blood sample may degrade, but the purity is not affected.

Perform the following steps:

1. Use a new centrifuge tube. Add 200  $\mu\text{L}$  to 1000  $\mu\text{L}$  of fresh human whole blood and 1 $\times$  Buffer LYR into the tube. The adding volume of 1 $\times$  Buffer LYR is 5 times that of fresh human whole blood. Incubate the tube on ice for 10 to 15 minutes.

 **Tips**

- To mix thoroughly, the volume of the mixture of blood and 1 $\times$  Buffer LYR cannot exceed three-fourths of the tube height.
- During incubation, when the mixture is translucent, the red blood cell is lysed.
- The incubation time could be extended to 20 minutes if necessary.

2. Place the tube into a centrifuge with a speed of 2100 rpm (about 400  $\times g$ ), centrifuge it at  $4\text{ }^{\circ}\text{C}$  for 10 minutes, and remove the supernatant completely.
3. Add 1 $\times$  Buffer LYR whose volume is twice that of human whole blood into the centrifuge tube and resuspend the cell.
4. Place the tube into a centrifuge with a speed of 2100 rpm (about 400  $\times g$ ), centrifuge it at  $4\text{ }^{\circ}\text{C}$  for 10 minutes, and remove the supernatant completely.

## 4.2.2 Animal tissue

Perform the following steps:

1. Use a new 1.5 mL centrifuge tube. Add 1 mg to 20 mg of fresh or frozen tissue at  $-80\text{ }^{\circ}\text{C}$  into the tube.
2. Add 100  $\mu\text{L}$  to 500  $\mu\text{L}$  of Buffer LY and 2 to 5 RNase-free zirconia beads into the tube. Place the tube into an electronic homogenizer and start reaction with a frequency of 70 Hz at  $4\text{ }^{\circ}\text{C}$  for 1 minute.

 **Tips** After the sample is grinded by the Buffer LY, a large amount of bubbles will appear. Therefore, it is recommended to reserve partial absolute ethanol during lysis buffer preparation. Add the reserved absolute ethanol into the centrifuge tube containing pretreated samples to remove bubbles and transfer the samples into the plate for sample according to the requirement of total volume.

3. Aspirate the supernatant slowly for extraction.

## 4.2.3 Cell



**Tips** The following methods are applicable to extract 10 µg to 30 µg of total RNA from  $1 \times 10^6$  cultured eukaryotic cells.

### 4.2.3.1 Cell suspension for collection

Perform the following steps:

1. Add cells to a new 1.5 mL centrifuge tube and estimate the cell count.
2. Collect the cells to the tube and centrifuge it in a centrifuge at 300 ×g for 5 minutes.
3. Remove the supernatant of culture medium.

### 4.2.3.2 Trypsin treatment

Perform the following steps:

1. Add the cells to a new 1.5 mL centrifuge tube, estimate the cell count, and remove the culture medium.
2. Add the PBS solution to wash cells and remove the PBS solution.
3. Add the PBS solution containing 0.1% to 0.25% trypsin into the tube.
4. When the cells detach from the wall of the tube, add the culture medium with serum to inactivate trypsin. Transfer the cell solution into a RNase-free centrifuge tube and centrifuge it in a centrifuge at 300 ×g for 5 minutes.
5. Collect cell pellets and remove the supernatant.



**Tips** When collecting cells, completely remove the cell culture medium. Otherwise, the cell may not be lysed completely, which affects the combination between RNA and magnetic beads, and even the RNA yield.


## 4.2.4 Bacteria

Perform the following steps:


1. Use a new 1.5 mL centrifuge tube. Estimate the bacteria number and collect the bacteria to the tube and centrifuge it in a centrifuge at 500 ×g for 5 minutes.
2. Remove the supernatant of culture medium.
3. Add 100 µL of lysozyme solution into the tube. The concentration and dissolution method are as follows:

**Table 5 Lysozyme concentration and dissolution methods for different types of bacteria**

Bacteria type	Concentration of lysozyme	Dissolution method
G- bacteria	1 mg/mL	1xTE Buffer
Most G+ bacteria	15 mg/mL	1xTE Buffer with 1.2% TritonX-100

-  **Tips**
- Most G+ bacteria here represents *Clostridium butyricum*, *Clostridia sporogenes* and *Listeria monocytogenes*.
  - For the extraction from *Staphylococcus aureus* and *Staphylococcus epidermidis*, mix 50 mM Tris-HCl (pH 7.5) and lysostaphin to prepare lysostaphin solution (2.4 U/ $\mu$ L). Store the solution at -20 °C and avoid frequent freeze-thaw cycle.

4. Vortex the tube to mix it thoroughly. Place the G- bacteria and most G+ bacteria at room temperature for 3 to 5 minutes and 5 to 10 minutes, respectively.

-  **Tips** For *Staphylococcus aureus* and *Staphylococcus epidermidis*, add 10  $\mu$ L of lysostaphin solution and 50  $\mu$ L to 90  $\mu$ L of PBS solution into the tube, vortex the tube for 10 seconds and incubate it for 30 minutes. At this time, extract the nucleic acids from the sample directly.


5. Add 300  $\mu$ L of Buffer LY and 20  $\mu$ L of Proteinase K into the tube. Mix the tube thoroughly and let stand for 5 minutes.
6. Add 60  $\mu$ L of Magnetic Beads T and 300  $\mu$ L of absolute ethanol into the tube. Mix the tube thoroughly and let stand for 8 minutes.

## 4.2.5 Yeast

-  **Tips** It is recommended that the number of yeasts should not exceed  $5 \times 10^7$ .

Perform the following steps:

1. Estimate the count of the bacteria, add the bacteria to a new 1.5 mL centrifuge tube and centrifuge the tube in a centrifuge at 500  $\times g$  for 5 minutes.
2. Slowly remove the culture medium.
3. Add the acid-washed glass beads into the tube until the volume reaches 100  $\mu$ L mark. Then, add 300  $\mu$ L of Buffer LY into the tube.
4. Place the tube into an electronic homogenizer and start reaction with a frequency of 70 Hz at 0 °C for 1 minute. Grind the sample for 2 times after 30 seconds.

-  **Tips** A large amount of bubbles may appear in the tube after grinding. Just add 300  $\mu$ L of absolute ethanol to remove the bubbles.

- Aspirate the supernatant for extraction. Do not touch the glass beads in the process.

## 4.3 Extracting the nucleic acids



- Tips**
- You can extract the nucleic acids manually or on automation devices. For automated nucleic acid extraction, ensure that you prepare applicable consumables.
  - During extraction on MGISP-960RS or MGISP-NE384RS, pretreat cell, bacteria, human whole blood and other samples, according to the requirement of 50  $\mu\text{L}$ /well, use an appropriate volume of 1 $\times$  PBS solution to disperse samples evenly, and use pipette to transfer samples into the 96-well plate for the lysis buffer.

### 4.3.1 Extracting the nucleic acids manually

Perform the following steps:


- Take out a 1.5 mL centrifuge tube and add reagents into the tube to prepare the lysis buffer for different samples according to table below:

**Table 6 Proportion for lysis buffer**

Sample type	Human whole blood	Animal tissue	Cell	Bacteria	Yeast
Magnetic Beads T	60 $\mu\text{L}$	60 $\mu\text{L}$	60 $\mu\text{L}$	60 $\mu\text{L}$	60 $\mu\text{L}$
Proteinase K	20 $\mu\text{L}$	20 $\mu\text{L}$	20 $\mu\text{L}$	20 $\mu\text{L}$	20 $\mu\text{L}$
Buffer LY	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$
Absolute ethanol	400 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$
$\beta$ -Mercaptoethanol	6 $\mu\text{L}$	/	/	/	/



- Tips**
- The prepared lysis buffer is required to be used within 30 minutes. If you want to prepare the buffer in advance, add Proteinase K before adding the buffer to the plate to avoid inactivating Proteinase K caused by long-time preparation.
  - During extraction from the human whole blood, add  $\beta$ -Mercaptoethanol, whose volume is 2% that of Buffer LY, to the tube containing Buffer LY. For example, for 1 mL of Buffer LY, you need to add 20  $\mu\text{L}$  of  $\beta$ -Mercaptoethanol. Use the reagent immediately after preparation. The prepared reagent can be stored at 4  $^{\circ}\text{C}$  for 1 month. If solid object appears, heat the reagent at 37  $^{\circ}\text{C}$  to redissolve it.

2. Add the lysis buffer into the centrifuge tube with pretreated samples. Vortex the tube to mix thoroughly and place at room temperature for 8 minutes during which vortex the tube 2 to 3 times with 5 seconds for each time.
3. Place the tube on the magnetic rack and let stand for 2 to 3 minutes during which slowly invert the tube on the magnetic rack to wash the Magnetic Beads T on the tube wall and cap. When Magnetic Beads T is adsorbed completely, use a pipette to remove the supernatant.
4. Add 700  $\mu\text{L}$  of Buffer WB I to the tube. Vortex the tube for 1 minute to mix thoroughly and place the tube on the magnetic rack for 1 minute during which slowly invert the tube on the magnetic rack. When Magnetic Beads T is adsorbed completely, use a pipette to remove the supernatant.
5. (Optional) For extraction from the human whole blood, repeat step 4 once. Add 700  $\mu\text{L}$  of Buffer WB II to the tube. Vortex the tube for 1 minute to mix thoroughly, briefly centrifuge and place the tube on the magnetic rack for 1 minute. When Magnetic Beads T is adsorbed completely, use a pipette to remove the supernatant.
6. Decap and dry the tube for 5 to 10 minutes.
7. Add 8  $\mu\text{L}$  of DNase I and 72  $\mu\text{L}$  of Buffer RDD into the tube. Vortex the tube for 1 minute to mix thoroughly and place the tube at room temperature for 15 minutes during which vortex the tube for 10 seconds every 5 minutes.
8. Add 700  $\mu\text{L}$  of Buffer WB I or Buffer WB II (for human whole blood extraction only) into the tube. Vortex the tube for 1 minute to mix thoroughly and let stand at room temperature for 3 minutes. Centrifuge the tube briefly and place on magnetic rack for 1 minute during which slowly invert the tube on the magnetic rack. When Magnetic Beads T is adsorbed completely, use a pipette to remove the supernatant.
9. Add 700  $\mu\text{L}$  of Buffer WB II to the tube. Vortex the tube for 1 minute to mix thoroughly, centrifuge briefly and place the tube on the magnetic rack for 1 minute. When Magnetic Beads T is adsorbed completely, use a pipette to remove the supernatant.
10. Repeat step 9 once.  
 **Tips** Skip this step for human whole blood extraction.
11. Decap and dry the tube for 5 to 10 minutes.
12. Add 80  $\mu\text{L}$  of RNase Free Water to the tube. Vortex the tube for 1 minute and let stand at room temperature for 5 minutes. Centrifuge briefly and place the tube on the magnetic rack. When Magnetic Beads T is adsorbed completely, aspirate the supernatant. The aspirated supernatant is the required product. Place the product at  $-80\text{ }^{\circ}\text{C}$  for storage.

## 4.3.2 Extracting the nucleic acids automatically on MGISP-960RS

### 4.3.2.1 Preparing consumables

According to the following table, prepare consumables for a workflow of automated extraction on MGISP-960RS and place them at room temperature until use:

**Table 7 Required automated consumables for extraction on MGISP-960RS**

Name	Brand	Cat. No.	Quantity
250 $\mu$ L automated filter tips	MGI	1000000723	7
2.2 mL V-bottom 96-well deep-well plate	MGI	1000008088	4 or 6 (for human whole blood extraction only)
1.3 mL U-bottom 96-well deep-well plate	MGI	1000004644	3 or 2 (for human whole blood extraction only)
Hard-shell thin-wall 96-well skirted PCR plates	MGI	1000012059	1

### 4.3.2.2 Preparing samples

You can extract 1 to 96 samples on MGISP-960RS.

Perform the following steps:

1. Ensure that samples to be extracted have been pretreated.
2. Add samples to the 96-well deep-well plate with the volume no more than 50  $\mu$ L for tissue, cell, bacterial or yeast samples, or 200  $\mu$ L for human whole blood extraction for each well, and pipette the samples 2 to 3 times. Ensure that no liquid exists on the wall of the well.



**Tips** When the sample volume is greater than recommended, magnetic beads will be clustered, for which a part of beads is aspirated during liquid transfer and loss is caused. For example, when the input of liver cell ranges between 10 mg and 30 mg, magnetic beads may be clustered. The reduction of beads affects the yield but it does not affect the purity.

3. Place the 96-well deep-well plate with samples on ice until use.

### 4.3.2.3 Preparing reagents

Perform the following steps:

1. Add absolute ethanol into the Buffer WB I according to the label.



2. Add absolute ethanol into the Buffer WB II according to the label.
3. Take out a 1.5 mL centrifuge tube and add reagents into the tube to prepare the lysis buffer for different samples according to table below:

**Table 8 Proportion for lysis buffer**

Sample type	Human whole blood	Animal tissue	Cell	Bacteria	Yeast
Magnetic Beads T	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L
Proteinase K	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L
Buffer LY	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L
Absolute ethanol	400 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L
$\beta$ -Mercaptoethanol	6 $\mu$ L	/	/	/	/



- Tips**
- The prepared lysis buffer is required to be used within 30 minutes. If you want to prepare the buffer in advance, add Proteinase K before adding the buffer to the plate to avoid inactivating Proteinase K caused by long-time preparation.
  - During extraction from the human whole blood, add  $\beta$ -Mercaptoethanol, whose volume is 2% that of Buffer LY, to the tube containing Buffer LY. For example, for 1 mL of Buffer LY, you need to add 20  $\mu$ L of  $\beta$ -Mercaptoethanol. Use the reagent immediately after preparation. The prepared reagent can be stored at 4  $^{\circ}$ C for 1 month. If solid object appears, heat the reagent at 37  $^{\circ}$ C to redissolve it.

4. Prepare three or two (only for human whole blood extraction) 1.3 mL U-bottom 96-well deep-well plates, three or six (only for whole blood extraction) 2.2 mL V-bottom 96-well deep-well plates and a hard-shell thin-wall 96-well skirted PCR plate. Add reagents according to the following table:






**Table 9 Required reagent volume for each plate**




Reagent name	Adding volume	Plate type
RNase Free Water	100 $\mu$ L/well	1.3 mL U-bottom 96-well deep-well plate
Buffer WB I	1500 $\mu$ L/well	2.2 mL V-bottom 96-well deep-well plate
Buffer WB II	1500 $\mu$ L/well	2.2 mL V-bottom 96-well deep-well plate

Reagent name	Adding volume	Plate type
Buffer WB II (only for human whole blood extraction)	800 $\mu$ L/well	2.2 mL V-bottom 96-well deep-well plate
DNase I + Buffer RDD	8 $\mu$ L/well + 72 $\mu$ L/well	1.3 mL U-bottom 96-well deep-well plate
Lysis buffer	680 $\mu$ L/well or 786 $\mu$ L/well (only for human whole blood extraction)	1.3 mL U-bottom 96-well deep-well plate or 2.2 mL V-bottom 96-well deep-well plate (only for human whole blood extraction)
Waste plate	/	2.2 mL V-bottom 96-well deep-well plate
Waste plate	/	2.2 mL V-bottom 96-well deep-well plate
Product plate	/	Hard-shell thin-wall 96-well skirted PCR plates

#### 4.3.2.4 Starting extraction

Perform the following steps:

1. Switch to the  position to power on the device.
2. Turn on the computer and the desktop appears. Double-tap  to run the software.
3. Select **User** and **Real**. Enter the password.
4. Tap **Login** to enter the main interface.
5. On the upper-right corner of the control software, tap  and select **WDesigner**. The home interface is displayed.
6. Ensure that the application file in the .wfex format has been prepared.
7. Tap  in the toolbar and find the file location in the pop-up window.
8. Select the file and tap **Open**, fill in the **Application** and **Project**, and tap **Confirm** to save the application file. Then this application file can be executed in the control software.
9. After the file is imported successfully, tap  in the toolbar.

10. Tap **Initialize** on the top of the interface to start initializing.  
You will be prompted after a successful initialization.
11. Tap  on the left of the interface, and select **Clean > Pre-clean > Start**.
12. Follow the on-screen instructions to complete operations and tap **Continue**.  
The UV lamp and air filter start working.  
 **CAUTION** Ultraviolet radiation is harmful to the human body, so do not open the door after the pre-clean starts.
13. Tap  > **Run Wizard** to enter the Run Wizard interface.
14. Tap the drop-down list of **Solution** and select **JB-A09-137 MGIEasy Total RNA Extraction Kit\_RV1.0\_SV1.0**. Tap the drop-down list of **Script**, select **JB-A09-137 MGIEasy Total RNA Extraction Kit\_RV1.0\_SV1.0.py** for cell, tissue, bacteria and yeast samples, and select **JB-A09-137 MGIEasy Total RNA Extraction Kit (Blood)\_RV1.0\_SV1.0.py** for the human whole blood sample. Place samples, reagents and consumables according to the following table and figure:

**Table 10 Sample, reagent, and consumable placement position**

Name	Position
250 µL automated filter tips	Pos1 to Pos7
RNase Free Water	Pos13
Buffer WB I	Pos14
Buffer WB II	Pos15
Buffer WB II (only for human whole blood)	Pos23
Lysis buffer	Pos20
DNase I + Buffer RDD	Pos21
Waste plate	Pos16
Waste plate	Pos18
RNA product (empty PCR plate)	Pos12




Figure 1 Plate map for tissue, cell, bacterial or yeast samples




Figure 2 Plate map for human whole blood sample

15. Tap **Start**. The extraction workflow starts. It takes roughly 1.5 hours.  
During the workflow, tap **Pause** to pause and tap **Resume** to resume the workflow if required.
16. At the end of the workflow, remove the RNA product from Pos12.  
If the product is not used immediately, seal and store it in a freezer at  $-80^{\circ}\text{C}$ .
17. Dispose of the used deep-well plates, PCR plates and waste bag.  
If no experiment is to be conducted on the rest of the day, clean the operation deck of the device according to *MGISP-100&MGISP-960 Cleaning Instructions*.

### 4.3.3 Extracting the nucleic acids automatically on MGISP-NE384RS

#### 4.3.3.1 Preparing consumables

According to the following table, prepare consumables for a workflow of automated extraction on MGISP-NE384RS and place them at room temperature until use:

**Table 11 Required automated consumables for extraction on MGISP-NE384RS**

Name	Brand	Cat. No.	Quantity
2.2 mL V-bottom 96-well deep-well plate	MGI	1000008088	96 preps: 6 or 7 (only for human whole blood)
			384 preps: 24 or 28 (only for human whole blood)
96-well tips comb	MGI	1000025661	1 (96 preps)
			4 (384 preps)

### 4.3.3.2 Preparing samples

You can extract 96 to 384 samples on MGISP-NE384RS.

Perform the following steps:

1. Ensure that samples to be extracted have been pretreated.
2. Add samples to the 96-well deep-well plate. Seal the plate and vortex it until the magnetic beads T is mixed thoroughly. Remove the seal and centrifuge the plate to ensure that no liquid exists on the wall of the well.



**Tips** When the sample volume is greater than recommended, magnetic beads will be clustered, for which a part of beads is remained in the tips comb during liquid adsorption and transfer. For example, when the input of liver cell ranges between 15 mg and 30 mg, magnetic beads may be clustered. The reduction of beads affects the yield but it does not affect the purity.

### 4.3.3.3 Preparing reagents

Perform the following steps:

1. Add absolute ethanol into the Buffer WB I according to the label.
2. Add absolute ethanol into the Buffer WB II according to the label.
3. Take out a 1.5 mL centrifuge tube and add reagents into the tube to prepare the lysis buffer for different samples according to table below:

**Table 12 Proportion for lysis buffer**

Sample type	Human whole blood	Animal tissue	Cell	Bacteria	Yeast
Magnetic Beads T	60 µL	60 µL	60 µL	60 µL	60 µL
Proteinase K	20 µL	20 µL	20 µL	20 µL	20 µL
Buffer LY	300 µL	300 µL	300 µL	300 µL	300 µL
Absolute ethanol	400 µL	300 µL	300 µL	300 µL	300 µL
β-Mercaptoethanol	6 µL	/	/	/	/



- Tips**
- The prepared lysis buffer is required to be used within 30 minutes. If you want to prepare the buffer in advance, add Proteinase K before adding the buffer to the plate to avoid inactivating Proteinase K caused by long-time preparation.
  - During extraction from the human whole blood, add  $\beta$ -Mercaptoethanol, whose volume is 2% that of Buffer LY, to the tube containing Buffer LY. For example, for 1 mL of Buffer LY, you need to add 20  $\mu$ L of  $\beta$ -Mercaptoethanol. Use the reagent immediately after preparation. The prepared reagent can be stored at 4 °C for 1 month. If solid object appears, heat the reagent at 37 °C to redissolve it.

4. Take out six or seven (only for human whole blood) 2.2 mL V-bottom 96-well deep-well plates and add reagents according to the following table:

**Table 13 Required reagent volume for each plate**

Reagent name	Adding volume
Lysis buffer and sample <b>Tips</b> For bacteria samples, you only need to add the pretreated samples.	680 $\mu$ L/well or 786 $\mu$ L/well (only for human whole blood)
Buffer WB I	700 $\mu$ L/well
Buffer WB I (only for human whole blood)	700 $\mu$ L/well
Buffer WB II	700 $\mu$ L/well
Buffer WB II	700 $\mu$ L/well
RNA product (containing RNase Free Water)	80 $\mu$ L/well

#### 4.3.3.4 Starting extraction


Perform the following steps:

1. Switch to the position to power on the device.
2. Turn on the computer and the desktop appears. Double-tap the icon of MGISP-NE384RS to run the software.
3. Select **User** and **Real**, and enter the password. Tap **Login** to enter the main interface.
4. Tap **Initialize** on the top of the interface to start initializing.  
You will be prompted after a successful initialization.
5. Tap **Process manage** > to import the script.
6. Tap > **Workflow**.

7. Select the script and place the sample, reagents and consumables according to different sample types.


- For cell, tissue, bacteria and yeast samples, tap the drop-down list of **Script**, select **MGIEasy Total RNA Extraction Kit\_V1.0**, and place the sample, reagents and consumables according to the following table:

**Table 14 Sample, reagent, and consumable placement position**

Reagent name	Position
Lysis buffer and sample  <b>Tips</b> For bacteria samples, you only need to add the pretreated samples.	Pos1
Buffer WB I	Pos2
Buffer WB II	Pos4
Buffer WB II	Pos5
RNA product (with RNase Free Water)	Pos6

- For the human whole blood sample, tap the drop-down list of **Script**, select **MGIEasy Total RNA Extraction Kit (Blood)\_V1.0** and place the sample, reagents and consumables according to the following table:

**Table 15 Sample, reagent, and consumable placement position**

Reagent name	Position
Lysis buffer and sample  <b>Tips</b> For bacteria samples, you only need to add the pretreated samples.	Pos1
Buffer WB I	Pos2
Buffer WB I (Used to replace the original one during human whole blood extraction)	Pos2
Buffer WB II	Pos3
Buffer WB II	Pos5
RNA product (with RNase Free Water)	Pos6

8. Tap **Run**. The device starts extraction according to the following table.



**CAUTION** Please add reagents manually according to the prompts within 15 minutes. Otherwise, RNA may degrade, which causes experimental failure.

During the workflow, tap **Pause** to pause and tap **Resume** to resume the workflow if required.

■ Extraction from tissue, cell, bacterial or yeast samples

Before starting step 3, a message box appears to confirm that you have placed the deep-well plate with 8 µL of DNase I and 72 µL of Buffer RDD in each well into Pos3. Tap **Confirm**. Step 3 starts.

Before starting step 4, a message box appears to confirm that you have added 700 µL of Buffer WB I into Pos3. Tap **Confirm**. Step 4 starts.

**Table 16 Automated extraction settings for tissue, cell, bacterial, or yeast samples**

Step No.	1	2	3	4	5	6	7	8
Step name	Lysis	Wash	Bind (manually)	Wash (manually)	Wash	Wash	Elution	Release
Position	1	2	3	3	4	5	6	2
Volume (µL)	680	700	80	780	700	700	80	700
Delay time (s)	0	0	300	20	0	0	300	0
Mix	True	True	True	True	True	True	True	True
Mix type	Normal	Normal	Normal	Normal	Normal	Normal	Normal	Normal
Mix rate	HighMiddle	HighMiddle	High	HighMiddle	HighMiddle	HighMiddle	High	High
Mix time (s)	400	60	600	180	60	60	180	10
Collect	True	True	False	True	True	True	True	False
Collect mode	Cycle	Cycle	/	Cycle	Cycle	Cycle	Cycle	/
Collect cycle	6	3	1	2	2	2	6	1
Collect time (s)	1	1	1	1	1	1	1	1
Dialog	False	True	True	False	False	False	False	True

■ Extraction from the human whole blood sample

Before starting step 3, a message box appears to confirm that you have placed a new plate containing Buffer WB I into Pos2. Tap **Confirm**. Step 3 starts.

Before starting step 5, a message box appears to confirm that you have placed the deep-well plate with 8 µL of DNase I and 72 µL of Buffer RDD in each well into Pos4. Tap **Confirm**. Step 5 starts.



Before starting step 6, a message box appears to confirm that you have added 700 µL of Buffer WB II into Pos4. Tap **Confirm**. Step 6 starts.

**Table 17 Automated extraction settings for the human whole blood sample**

Step No.	1	2	3	4	5	6	7	8	9
Step name	Lysis	Wash	Wash	Wash	Bind	Wash	Wash	Elution	Release
Position	1	2	2	3	4	4	5	6	2
Volume (µL)	786	700	700	700	80	780	700	80	700
Delay time (s)	0	0	10	0	300	20	0	300	0
Mix	True	True	True	True	True	True	True	True	True
Mix type	Normal	Normal	Normal	Normal	Normal	Normal	Normal	Normal	Normal
Mix time (s)	480	60	60	60	600	180	60	180	10
Collect	True	True	True	True	False	True	True	True	False
Collect mode	Cycle	Cycle	Cycle	Cycle	/	Cycle	Cycle	Cycle	/
Collect cycle	6	3	3	2	/	2	2	6	/
Collect time (s)	1	1	1	1	1	1	1	1	1
Dialog	False	True	False	True	True	False	False	False	True

9. After the program ends, transfer the 96-well tips comb to the medical waste bag.
10. Remove the 96-well plate from Pos6 and transfer the product to a new tube.  
If the product is not used immediately, seal it and store it in a freezer at -80 °C .

### 4.3.4 Extracting the nucleic acids automatically on MGISP-NE32RS

#### 4.3.4.1 Preparing consumables

According to the following table, prepare consumables for a workflow of automated extraction on MGISP-NE32RS and place them at room temperature until use:


**Table 18 Required automated consumables for extraction on MGISP-NE32RS**

Name	Brand	Cat. No.	Quantity
2.2 mL 96-well V-bottom deep-well plate (MGISP-NE32)	MGI	091-000444-00	2 or 4 (only for human whole blood)
Plastic Magnetic Wand Cover	MGI	1000022599	2

#### 4.3.4.2 Preparing samples

You can extract 1 to 32 samples on MGISP-NE32RS.

Ensure that samples to be extracted have been pretreated.

-  **Tips**
- For pretreated samples such as precipitated cells and bacteria, it is recommended to use 50  $\mu$ L to 100  $\mu$ L of 1x PBS to disperse the samples before transferring samples to the plate containing lysis buffer for extraction.
  - When the sample volume is greater than recommended, magnetic beads will be clustered, for which a part of beads is remained in the tips comb during liquid adsorption and transfer. For example, when the input of liver cell ranges between 15 mg and 30 mg, magnetic beads may be clustered. The reduction of beads affects the yield but it does not affect the purity.

#### 4.3.4.3 Preparing reagents

Perform the following steps:

1. Add absolute ethanol into the Buffer WB I according to the label, and mark it on the label. Mix the tube and store it at room temperature until use.
2. Add absolute ethanol into the Buffer WB II according to the label, and mark it on the label. Mix the tube and store it at room temperature until use.
3. Take out a 1.5 mL centrifuge tube and add reagents into the tube to prepare the lysis buffer for different samples according to table below:

**Table 19 Proportion for lysis buffer**

Sample type	Human whole blood	Animal tissue	Cell	Bacteria	Yeast
Magnetic Beads T	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L	60 $\mu$ L
Proteinase K	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L	20 $\mu$ L
Buffer LY	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L
Absolute ethanol	400 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L	300 $\mu$ L
$\beta$ -Mercaptoethanol	6 $\mu$ L	/	/	/	/



- Tips**
- The prepared lysis buffer is required to be used within 30 minutes. If you want to prepare the buffer in advance, add Proteinase K before adding the buffer to the plate to avoid inactivating Proteinase K caused by long-time preparation.
  - During extraction from the human whole blood, add  $\beta$ -Mercaptoethanol, whose volume is 2% that of Buffer LY, to the tube containing Buffer LY. For example, for 1 mL of Buffer LY, you need to add 20  $\mu$ L of  $\beta$ -Mercaptoethanol. Use the reagent immediately after preparation. The prepared reagent can be stored at 4 °C for 1 month. If solid object appears, heat the reagent at 37 °C to redissolve it.

4. According to different sample types, take out an appropriate number of 96-well deep-well plates and add the required reagents. To avoid confusion, designate the left side of the device as Position A and the right side as Position B.
  - For cell, tissue, bacterial or yeast samples, take out two 2.2 mL 96-well V-bottom deep-well plates and label them as "A1" and "B1", respectively. Add an appropriate amount of reagent according to the table below.

**Table 20 Required reagent volume for each plate (cell, tissue, bacterial or yeast samples)**

Plate name	Reagent name	Col. No.	Adding volume/well
A1/B1	Lysis buffer	1 and 7	680 $\mu$ L
	Buffer WB I	2 and 8	700 $\mu$ L
	DNase I	3 and 9	8 $\mu$ L
	Buffer RDD		72 $\mu$ L
	Buffer WB II	4 and 10	700 $\mu$ L
	RNase Free Water	5 and 11	80 $\mu$ L
	Buffer WB II	6 and 12	700 $\mu$ L

- For the human whole blood sample, take out four 2.2 mL 96-well V-bottom deep-well plates and label them as "A1", "B1", "A2", and "B2", respectively. Add an appropriate amount of reagent according to the table below.


**Table 21 Required reagent volume for each plate (human whole blood sample)**


Plate name	Reagent name	Col. No.	Adding volume/well
A1/B1	Lysis buffer	1 and 7	784 $\mu$ L
	Buffer WB I	2 and 8	700 $\mu$ L
	Buffer WB I	3 and 9	700 $\mu$ L
	Buffer WB II	4 and 10	700 $\mu$ L
A2/B2	DNase I	1 and 7	8 $\mu$ L
	Buffer RDD		72 $\mu$ L
	RNase Free Water	2 and 8	80 $\mu$ L
	Buffer WB II	5 and 11	700 $\mu$ L

#### 4.3.4.4 Starting extraction

Perform the following steps:

1. Power on MGISP-NE32RS, and the device starts self-test.  
Self-test takes approximately 10 seconds.
2. Configure a script by using either of the following methods:
  - In the main interface, tap **New** to enter the file editing interface. Edit the script according to step 5 below for different samples.
  - Insert the USB drive, and tap **Run** in the main interface. Select the script file from the USB drive, tap **Import** in the lower left corner, and you will be prompted that the script is successfully imported.
3. Tap **Device**, and select a script according to the sample type.
  - For cell, tissue, bacterial or yeast samples, select **JB\_007\_V1**.
  - For the human whole blood sample, tap the drop-down list of **Script**, select **JB\_007\_B\_V1**.
4. Place the pre-filled Plate A1 into Position A and Plate B1 into Position B, and install the plastic magnetic wand cover.
 

 **Tips** To avoid confusion, designate the left side of the device as Position A and the right side as Position B.
5. Tap **Run**. The device starts extraction according to the following table.
 

 **CAUTION** Please add reagents manually according to the prompts within 15 minutes. Otherwise, RNA may degrade, which causes experimental failure.

During the workflow, tap **Stop** to stop the workflow, **Pause** to pause the workflow, and tap **Pause/Reset** to suspend the workflow and raise the magnetic rods if required.

■ Extraction from tissue, cell, bacterial or yeast samples

When the workflow reaches step 6, open the compartment door, remove the deep-well plate, and manually add 700 µL of Wash Buffer I to wells in columns 3 and 9 within 30 seconds. The workflow automatically pauses at this point. After adding the reagent, return the deep-well plate to its corresponding position, close the compartment door, and the workflow will automatically resume.

Temperature setting:

Lysis Temp: OFF

Elution Temp: OFF

**Table 22 Automated extraction settings for tissue, cell, bacterial, or yeast samples**

Step No.	Col. No.	Name	Wait time (min:ss)	Mix time (min:ss)	Mag time (min:ss)	Volume (µL)	Mixing method	Collect method
1	1	Lysis	00:00	08:00	01:30	680	Fast	Strong
2	1	Lysis	00:00	00:00	00:10 x 5	680	Slow	Cycle
3	2	Wash I	00:00	01:00	01:30	700	Fast	Strong
4	2	Wash I	00:00	00:00	00:10 x 5	700	Slow	Cycle
5	3	Bind	05:00	10:00	00:00	80	Fast	Strong
6	3	Wash I	00:30	03:00	01:30	780	Fast	Strong
7	3	Wash I	00:00	00:00	00:10 x 5	780	Slow	Cycle
8	4	Wash II	00:00	01:00	01:30	700	Fast	Strong
9	4	Wash II	00:00	00:00	00:10 x 5	700	Slow	Cycle
10	6	Wash II	00:00	01:00	01:30	700	Fast	Strong
11	6	Wash II	00:00	00:00	00:10 x 5	700	Slow	Cycle
12	5	Elute	05:00	03:00	01:30	80	Fast	Strong
13	5	Elute	00:00	00:00	00:10 x 5	80	Slow	Cycle
14	2	Beads	00:00	00:30	00:00	700	Fast	Normal

■ Extraction from the human whole blood sample

When the workflow reaches step 9, open the compartment door immediately, and replace Plate A1 with Plate A2 and Plate B1 with Plate B2. Close the compartment door, and the workflow will automatically resume.


When the workflow reaches step 10, open the compartment door, remove Plate A2 and Plate B2, and manually add 700 µL of Wash Buffer II to wells in columns 1 and 7 of these two plates within 20 seconds. The workflow automatically pauses at this point. After adding the reagent, return Plate A2 and Plate B2 to Position A and Position B, respectively. Close the compartment door, and the workflow will automatically resume.




**Tips** To avoid confusion, designate the left side of the device as Position A and the right side as Position B. Place the pre-filled Plate A1 into Position A and Plate B1 into Position B, and replace Plate A1 with Plate A2 and Plate B1 with Plate B2 when the workflow reaches step 9,

**Table 23 Automated extraction settings for the human whole blood sample**

Plate name	Step No.	Col. No.	Name	Wait time (min:ss)	Mix time (min:ss)	Mag time (min:ss)	Volume (µL)	Mixing method	Collect method
A1/A2	1	1	Lysis	00:00	08:00	02:00	786	Fast	Strong
	2	1	Lysis	00:00	00:00	00:20 x 5	786	Slow	Cycle
	3	2	Wash I	00:00	01:00	02:00	700	Fast	Strong
	4	2	Wash I	00:00	00:00	00:20 x 5	700	Slow	Cycle
	5	3	Wash I	00:00	01:00	02:00	700	Fast	Strong
	6	3	Wash I	00:00	00:00	00:20 x 5	700	Slow	Cycle
	7	4	Wash II	00:00	01:00	01:00	700	Fast	Strong
	8	4	Wash II	00:00	00:00	00:10 x 3	700	Slow	Cycle
B1/B2	9	1	Bind	05:00	10:00	00:00	80	Fast	Strong
	10	1	Wash II	00:20	03:00	01:00	780	Fast	Strong
	11	1	Wash II	00:00	00:00	00:10 x 3	780	Slow	Cycle
B1/B2	12	2	Wash II	00:00	01:00	01:00	700	Fast	Strong
	13	2	Wash II	00:00	00:00	00:10 x 3	700	Slow	Cycle
	14	5	Elute	05:00	03:00	01:00	80	Fast	Strong
	15	5	Elute	00:00	00:00	00:10 x 5	80	Slow	Cycle
	16	2	Beads	00:00	00:30	00:00	700	Fast	Strong

 **CAUTION** After the workflow is completed, wait for the completion alert sound to finish and ensure that the robotic arm is no longer moving before opening the compartment door. Otherwise, the device may become jammed.

6. After the program ends, remove and transfer the plastic magnetic wand cover to the medical waste bag.
7. Take out the two plates, transfer the extracted nucleic acid from columns 5 and 11 into new 8-strip tubes, and securely cap the tubes. The extracted product can be used directly for subsequent experiments or stored at -80 °C until use.

 **CAUTION** After the workflow is completed, immediately take out the deep-well plates. Do not leave the product in the temperature control position for an extended period, as this may affect the product quality.

8. Dispose of the used deep-well plate in the designated waste area. Use lint-free paper moistened with 75% alcohol to wipe the interior surfaces of the compartment, and then close the compartment door. In the main interface, tap **UV Lamp**, and set the time to 30 minutes. Tap **Confirm** to turn on the UV lamp.

### 4.3.5 Extracting the nucleic acids automatically on MGISP-NEXRS


#### 4.3.5.1 Preparing consumables

According to the following table, prepare consumables for a workflow of automated extraction on MGISP-NEXRS and place them at room temperature until use:

**Table 24 Required automated consumables for extraction on MGISP-NEXRS**

Name	Brand	Cat. No.	Quantity
1000 µL Black Conductive Tips, Sterile, Filtered (Wide-Bore)	MGI	091-000020-00	Refer to <i>MGI Easy Total RNA Extraction Set Reagent Adding Volume and Consumable Consumption Calculation Table</i>
1000 µL Tips, Conductive, Filtered, Suspended	MGI	091-000223-00	
200 µL Tips, Transparent, Filtered, Suspended	MGI	091-000158-00	

Name	Brand	Cat. No.	Quantity
200 $\mu$ L Tips, Conductive, Filtered, Suspended	MGI	091-000156-00	Refer to <i>MGIEasy Total RNA Extraction Set Reagent Adding Volume and Consumable Consumption Calculation Table</i>
50 $\mu$ L Tips, Transparent, Filtered, Suspended	MGI	091-000157-00	
50 mL Single Well Reservoir	MGI	012-000780-00	
100 mL Single Well Reservoir	MGI	012-000779-00	
2.2 mL 96-well V Bottom deep well plate	MGI	091-000287-00	
Hard-shell Thin-wall 96-well Skirted PCR Plates, White Shell/Clear Well	MGI	091-000165-00	
96-well tips comb	MGI	1000025661	


 **Tips** *MGIEasy Total RNA Extraction Kit Reagent Adding Volume and Consumable Consumption Calculation Table* is included in the script package. If you are unable to obtain it, please contact technical support.

#### 4.3.5.2 Preparing samples

You can extract 8 to 32 samples on MGISP-NEXRS.

Perform the following steps:

1. Ensure that samples to be extracted have been pretreated.
2. For tissue, cell, bacterial or yeast samples, vortex the tube for 10 seconds. For the human whole blood sample, invert the tube more than 10 times.
3. After mixing the sample well, open the tube cap and place the sample in Tube Carrier A for later use.

 **Tips** This workflow is compatible with sample loaded with standard 6 mL blood collection tubes, with dimensions of 13 mm $\times$ 100 mm.

#### 4.3.5.3 Preparing reagents

Perform the following steps:

1. Remove all reagents from the set until use.




2. Take out 2 or 4 (if the number of samples is greater than 48) 100 mL Single Well Reservoir, one 2.2 mL 96-well V-bottom Deep-well Plate, one 2.0 mL cryovial, and four 50 mL Single Well Reservoir. Label the consumables according to the table below. Add an appropriate amount of reagents to corresponding consumables based on the calculation results from the *MGIEasy Total RNA Extraction Kit Reagent Adding Volume and Consumable Consumption Calculation Table*.




- Tips**
- *MGIEasy Total RNA Extraction Kit Reagent Adding Volume and Consumable Consumption Calculation Table* is included in the script package. If you are unable to obtain it, please contact technical support.
  - If the number of samples is not greater than 48, only 2 of the 100 mL single-well reservoirs are needed.
  - The volume graduation lines on the single-well reservoirs are not entirely accurate, so it is not recommended to use the single-well reservoirs as measuring containers.

**Table 25 Required reagent volume for each plate**

Reagent name	Consumable type	Consumable label	Adding volume/well
Buffer WB I	100 mL Single Well Reservoir	Buffer WB I-1	Refer to <i>MGIEasy Total RNA Extraction Set Reagent Adding Volume and Consumable Consumption Calculation Table</i>
Buffer WB I		Buffer WB I-2	
Buffer WB II		Buffer WB II-1	
Buffer WB II		Buffer WB II-2	
Add the following two reagents to the same well. For the well position, refer to <i>Figure 3 on Page 29</i> : <ul style="list-style-type: none"> <li>• Magnetic Beads T</li> <li>• Proteinase K</li> </ul>	2.2 mL 96-well V-bottom Deep-well Plate	Magnetic Beads T+Proteinase K	
DNase I	2.0 mL cryovial	DNase I	
Buffer RDD	50 mL Single Well Reservoir	Buffer RDD	
Buffer LY		Buffer LY	
Absolute ethanol		Absolute ethanol	
RNase Free Water		RNase Free Water	






**Figure 3 Adding well position for Magnetic Beads T and Proteinase K**

 **Tips** Magnetic beads T are prone to aggregation and should be thoroughly vortexed to mix well before adding to the wells.

3. Prepare an appropriate number of consumables according to the following table, and label them.

**Table 26 Consumable quantity and label reference table**


Consumable type	Quantity	Consumable label
1000 $\mu$ L Black Conductive Tips, Sterile, Filtered (Wide-Bore)  <b>Tips</b> This consumable is only applicable to human whole blood samples loaded with 6 mL blood collection tubes.	1	1000 $\mu$ L Tips Conductive Wide-Bore
200 $\mu$ L Tips, Transparent, Filtered, Suspended  <b>Tips</b> For the human whole blood sample, this consumable is only applicable to samples loaded with 1.5 mL centrifuge tubes.	1	200 $\mu$ L Tips
1000 $\mu$ L Tips, Conductive, Filtered, Suspended	1	1000 $\mu$ L Tips Conductive
200 $\mu$ L Tips, Conductive, Filtered, Suspended	1 or 2 (if the number of samples is greater than 72)	200 $\mu$ L Tips Conductive
50 $\mu$ L Tips, Transparent, Filtered, Suspended	1	50 $\mu$ L Tips

Consumable type	Quantity	Consumable label
Tube Adapter A, 2.0 mL, Black <b>💡 Tips</b> For the human whole blood sample, this consumable is only applicable to samples loaded with 1.5 mL centrifuge tubes.	n (the number of samples)	/
2.2 mL 96-well V-bottom deep-well plate	1	Wash Buffer I-1
2.2 mL 96-well V-bottom deep-well plate	1	Wash Buffer II-1
2.2 mL 96-well V-bottom deep-well plate	1	According to the sample type, label the plate as follows: <ul style="list-style-type: none"> <li>For cell, tissue, bacterial or yeast samples: Wash Buffer II-2</li> <li>For the human whole blood sample: Wash Buffer II-3</li> </ul>
2.2 mL 96-well V-bottom deep-well plate	1	DNase I+Buffer RDD
2.2 mL 96-well V-bottom deep-well plate	1	Combine according to the figure below and label it as "Mix+96-well tips comb":
96-well tips comb	1	 <p>96-well tips comb + Deep-well plate</p> <p>Mix+96-well tips comb</p>
2.2 mL 96-well V-bottom deep-well plate	1	RNA Wash
Hard-shell Thin-wall 96-well Skirted PCR Plates, White Shell/Clear Well	1	RNA

#### 4.3.5.4 Starting extraction

Perform the following steps:

1. Power on MGISP-NEXRS.
2. Double-tap the control software icon on the desktop to launch the software.

3. Select **user** from the drop-down list and input the password **123456** in the pop-up login interface and tap **Login** to enter the main interface.
4. Tap **Workflow** to enter the Workflow interface.
5. Tap **Initialize** on the top of the interface to start initializing.  
A prompt is displayed after completion.
6. Tap **Run** to enter the Run interface.
7. Tap **Browse**, and select a script according to the sample type.
  - For cell, tissue, bacterial or yeast samples, select **MGIEasy Total RNA Extraction Set. spproj**.
  - For the human whole blood sample, select **MGIEasy Total RNA Extraction Set (Blood-1). spproj**. if the sample is stored in 6 mL blood collection tubes, and select **MGIEasy Total RNA Extraction Set (Blood-2). spproj** if the sample is stored in 1.5 mL centrifuge tubes.
8. Tap  to open the deck preview picture. Place the sample, reagent and consumables according to the following figure and table.










Figure 4 Deck picture

Table 27 Sample, reagent, and consumable placement position

Consumable label	Position
Buffer WB I-1	Lane6, Single well reservoir 1
Buffer WB I-2	Lane6, Single well reservoir 2
Buffer WB II-1	Lane6, Single well reservoir 3
Buffer WB II-2	Lane6, Single well reservoir 4
Magnetic Beads T+Proteinase K	Lane13, Pos3
DNase I	Lane19, Pos2, 1A 
Buffer LY	Lane25, Pos5, Single well reservoir 1
Absolute ethanol	Lane25, Pos5, Single well reservoir 2
Buffer RDD	Lane25, Pos5, Single well reservoir 3
RNase Free Water	Lane25, Pos5, Single well reservoir 4
1000 $\mu$ L Tips Conductive Wide-Bore <b>Tips</b> This consumable is only applicable to human whole blood samples loaded with 6 mL blood collection tubes.	Lane7, Pos1
200 $\mu$ L Tips <b>Tips</b> For the human whole blood sample, this consumable is only applicable to samples loaded with 1.5 mL centrifuge tubes.	Lane7, Pos1
1000 $\mu$ L Tips Conductive	Lane7, Pos2
200 $\mu$ L Tips Conductive	Lane7, Pos3 或 Pos4
50 $\mu$ L Tips	Lane7, Pos5

Consumable label	Position
Tube Adapter A, 2.0 mL, Black  <b>Tips</b> For the human whole blood sample, this consumable is only applicable to samples loaded with 1.5 mL centrifuge tubes.	Lane1 to Lane4
Wash Buffer I	Lane19, Pos1
Wash Buffer II-1	Lane19, Pos3
Wash Buffer II-2  <b>Tips</b> This label is only applicable to cell, tissue, bacterial or yeast samples.	Lane19, Pos4
Wash Buffer II-3  <b>Tips</b> This label is only applicable to human whole blood samples.	Lane19, Pos4
DNase I+Buffer RDD	Lane25, Pos2
Mix+96-well tips comb	Lane25, Pos3
RNA Wash	Lane25, Pos4
RNA	Lane13, Pos4

9. After confirming that the samples, reagents, and consumables are correctly placed, close the door. Tap **Run** in the Run interface. Tap **OK** in the pop-up window.
10. Select the number of samples at the drop-down list in the pop-up window. The extraction workflow starts automatically.  
 (Optional) For the human whole blood sample, the sample type selection interface is displayed. Select **DNase I** at the drop-down list and tap **Continue**. The extraction workflow starts automatically.  
 **Tips** According to the number of samples, the run time ranges from 1 hour and 20 minutes for 8 samples to 2 hours and 30 minutes for 96 samples.
11. After the workflow is completed, follow the on-screen instructions to remove the RNA product plate from Lane 13, Pos 4. If the products are not to be used for subsequent experiments immediately, seal the plate and store it at -80 °C.
12. Tap **Continue** to end the workflow and return to the Run interface. Remove other deep-well plates, 96-well tips comb or waste bag, and dispose of them in the designated waste area.
13. Tap  to return to the Workflow interface. Tap **Clean** to clean the device.

---

## Chapter 5 Warnings and precautions

- This product is for research use only. Please read the instructions for use of the product carefully before use.
- Before experiment, be sure to be familiar with and master the operation methods and precautions of various devices to be used.
- Direct contact with skin and eyes should be avoided for all samples and reagents. Do not swallow. If accidental ingestion occurs, please get medical attention immediately. If skin exposure occurs, rinse with large amounts of water and get medical attention if irritation persists.
- All samples and wastes should be disposed of in accordance with relevant regulations.
- Do not use expired products.

---

## Appendix 1 Manufacturer information

Manufacturer	Wuhan MGI Tech Co., Ltd.
Address	Building 24, Stage 3.1, BioLake Accelerator, No.388, 2nd Gaoxin Road, East Lake High-Tech Development Zone, 430075, Wuhan, P.R. China
After-sale service address	Wuhan MGI Tech Co., Ltd.
E-mail	MGI-service@mgi-tech.com
Website	en.mgi-tech.com